

**510(K) SUMMARY****JUL 10 2008**

Submitter's name: Sonoma Orthopedic Products, Inc.

Address: 650 Larkfield Center, Suite C  
Santa Rosa, CA 95403

Phone Number: 707-526-1335 x-224

Fax Number: 707-526-2022

Contact Person: Amy Conuel  
Director Quality Assurance & Regulatory Affairs

Date Prepared: July 8, 2008

Trade Name: Ensplint™ Bone Screw

Common Name: Bone Screw

Classification Name: Smooth or threaded metallic bone fixation fastener

Predicate Device: Syntec-Taichung Non-Sterile Bone Plate and Screw Implants. K983495, cleared 12/16/1998

Device Description: HA Cortical Self-Tapping Bone Screws in diameters of 1.5mm to 5.0mm and lengths of 4mm to 90mm.

Intended Use: The Ensplint™ Bone Screws are intended to treat fractures of various small and long bones. The Ensplint™ Bone Screws can also be used with the EnsplintRx Distal Radius System.

Rationale for Substantial Equivalence / Comparison to Predicate: The Ensplint™ Bone Screw utilizes the same material and technology characteristics as the Syntec Scientific Corp. Bone Screws and conforms to the ASTM F543-2 Standard Specification for HA Bone Screws. The diameters and lengths are within the diameters and lengths of the predicate Syntec Scientific Corp bone screws. Therefore testing is not needed to demonstrate that the subject devices are substantially equivalent to other legally marketed bone screws.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Sonoma Orthopedic Products, LLC  
% Ms. Amy Conuel  
Director, Quality Assurance/Regulatory Affairs  
650 Larkfield Center, Suite C  
Santa Rosa, California 95403

JUL 10 2008

Re: K080778

Trade/Device Name: Ensplint™ Bone Screw  
Regulation Number: 21 CFR 888.3040  
Regulation Name: Smooth or threaded metallic bone fixation fastener  
Regulation Class: Class II  
Product Code: HWC  
Dated: June 24, 2008  
Received: June 25, 2008

Dear Ms. Conuel:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Amy Conuel

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at 240-276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at 240-276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## INDICATIONS FOR USE

510(k) Number (if known): K080778

Device Name: Ensplint™ Bone Screw

### Indications for Use:

The Ensplint™ Bone Screws are intended to treat fractures of various small and long bones. The Ensplint™ Bone Screws can also be used with the EnsplintRx Distal Radius System.

Prescription Use ☒ X  
(Part 21 CFR 801 Subpart D)

AND/OR


Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
OF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1

  
(Division Sign-Off)

**Division of General, Restorative,  
and Neurological Devices**

510(k) Number K080778