

510(k) SUMMARY**APR 23 2008****Orthocord #2-0 suture with or without needles****Submitter's Name and Address:**

DePuy Mitek
a Johnson & Johnson company
325 Paramount Drive
Raynham, MA 02767

Contact Person

Kristine Christo
Regulatory Affairs Project Lead
DePuy Mitek
a Johnson & Johnson company
325 Paramount Drive
Raynham, MA 02767
Telephone: 508 828-3359
Facsimile: 508-977-6955
e-mail: KChristo@Dpyus.jnj.com

Name of Medical Device

Classification Name: Non-absorbable poly(ethylene terephthalate) surgical suture; and Suture, Surgical, Absorbable, Polydioxanone
Common/Usual Name: Suture
Proprietary Name: Orthocord suture

Substantial Equivalence

Orthocord #2-0 suture is substantially equivalent to:

Orthocord #2 suture (K040004, K043298) manufactured by Mitek;
Mini QA+ / Minilok QA+ (K071622) manufactured by Mitek and
Fiberwire (K021434) manufactured by Arthrex

Device Classification

Sutures, classified by the FDA, are Class II Medical Devices.
PDS Suture carries an FDA product code NEW, and is classified as absorbable surgical suture, polydioxanone under 21 CFR 878.4840.
Polyethylene suture carries an FDA product code GAT, and is classified under 21 CFR 878.5000.

Device Description

ORTHOCORD suture is a synthetic, sterile, braided composite suture composed of dyed (D&C Blue #6 or D&C Violet #2) absorbable polydioxanone (PDS) and un-dyed non-absorbable polyethylene. The partially absorbable suture is coated with a copolymer of 90% caprolactone and 10% glycolide.

Indications for Use

ORTHOCORD suture is indicated for use in general soft tissue approximation and/or ligation, including use in orthopedic surgeries.

Safety and Performance

The determination of substantial equivalence for this device was based on a detailed device description, and conformance to consensus and voluntary standards. Non-clinical laboratory testing was performed demonstrating that the device conformed to the USP monograph for absorbable sutures.

Based on the indications for use, technological characteristics, and comparison to predicate devices, the ORTHOCORD suture has been shown to be substantially equivalent to predicate devices under the Federal Food, Drug and Cosmetic Act.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Depuy Mitek
% Ms. Kristine Christo
Regulatory Affairs Project Lead
325 Paramount Drive
Raynham, Massachusetts 02767

APR 23 2008

Re: K080918

Trade/Device Name: Orthocord #2-0 suture with or without needles
Regulation Number: 21 CFR 878.4840
Regulation Name: Absorbable polydioxanone surgical suture
Regulatory Class: II
Product Code: NEW, GAT
Dated: March 13, 2008
Received: April 1, 2008

Dear Ms. Christo:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Kristine Christo

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATIONS FOR USE

510(k) Number (if known): K080918

Device Names: Orthocord #2-0 suture with or without needles

Indications for Use:

Orthocord #2-0 suture is indicated for use in general soft tissue approximation and/or ligation, including use in orthopedic surgeries.

Prescription Use ✓
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Neil R. Ogle, Jr. FR MCM
(Division Sign-Off)

Division of General, Restorative
and Neurological Devices

Page 1 of 1

510(k) Premarket Notification (PMA) Special
ORTHOCORD #2-0 suture
with or without needles

510(k) Number K080918

Confidential