
JAN 23 2009

Summary

This summary of 510(k)-safety and effectiveness information is being submitted in accordance with the requirements of 21 CFR § 807.92

1 Company making the submission

Name	Gish BioMedical, Inc
Address	22942 Arroyo Vista Rancho Santa Margarita, CA 92688-2600
Telephone	949-635-6200 voice 949-635-6299 fax edw@gishbiomedical.com
Contact	Edward F. Waddell Director RA/QA

2 Device

Proprietary Name	Gish Tubing and Connectors with HA Coating
Common Name	Cardiopulmonary Bypass Tubing and Connectors
Classification Name	Catheter, Cannula or Tubing, Vascular Cardiopulmonary Bypass Adaptor, Stopcock, Manifold, or Fitting, Cardiopulmonary Bypass

3 Predicate Devices

Gish Tubing and Connectors Both manufactured by Gish Biomedical, Inc

4 Classifications Names & Citations

21 CFR 870.4210, Catheter, Cannula or Tubing, Vascular Cardiopulmonary Bypass, Class II, DWF, Cardiovascular

21 CFR 870.4290, Adaptor, Stopcock, Manifold, or Fitting, Cardiopulmonary Bypass, Class II, DTL, Cardiovascular

5 Description

The Gish Tubing and Connectors with HA Coating are used in surgical procedures to provide a conduit for extracorporeal blood flow when interconnecting components of the bypass circuit. Specifically, they are used in connecting oxygenators, reservoirs, filters, heat exchangers and other devices used in surgical procedures. The tubing is polyvinyl chloride (PVC). The connectors are polycarbonate of various configurations such as straight, "Y", Luer and reducer types.

The components of this system which have contact with the fluid path are sterile and nonpyrogenic

All materials of the Gish Tubing and Connectors are biocompatible and coated with a proprietary coating

The Gish Tubing and Connectors with HA Coating may be purchased separately or pre-connected with tubing and other components of an extracorporeal circuit

6 Indications for use

The Gish Tubing and Connectors with HA Coating are indicated for use in surgical procedures to provide a conduit for extracorporeal blood flow when interconnecting components of the bypass circuit. It is designed to operate at flow rates of one (1) to six (6) liters per minute for periods up to six (6) hours.

7 Contra-indications

For HA coated tubing and connectors, no contra-indications have been noted.

8 Comparison

The Gish Tubing and Connectors with HA Coating have the same device characteristics as the predicate devices.

9 Test Data

The Gish Tubing and Connectors with HA Coating has been subjected to extensive safety, performance, and validations prior to release. Final testing for the systems includes various performance tests designed to ensure that the device meets all of its functional requirements and performance specifications.

10 Literature Review

A review of literature pertaining to the safety and effectiveness has been conducted. Appropriate safeguards have been incorporated in the design of Gish Tubing and Connectors with HA Coating.

11 Conclusions

The conclusion drawn from these tests is that Gish Tubing and Connectors with HA Coating is equivalent in safety and efficacy to its predicated devices.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JAN 29 2009

Gish Biomedical, Inc
c/o Ms Janet Peets
Regulatory & Clinical Affairs Specialist
22942 Arroyo Vista
Rancho Santa Margarita, CA 92688

Re K081881
Gish Tubing and Connectors with HA Coating
Regulation Number 21 CFR 870.4210
Regulation Name Cardiopulmonary Bypass Cannula or Tubing
Regulatory Class Class II
Product Code DWF
Dated January 8, 2009
Received January 12, 2009

Dear Ms Peets

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to registration and listing (21 CFR Part 807), labeling (21 CFR Part 801), good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820), and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act), 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometrics' (OSB's) Division of Postmarket Surveillance at 240-276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at 240-276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number K 081881

Device Name: Gish Tubing and Connectors with HA Coating

Indications for use.

The Gish Tubing and Connectors with HA Coating are indicated for use in Surgical procedures to provide a conduit for extracorporeal blood flow when interconnecting components of the bypass circuit. It is designed to operate at flow rates of one (1.0) to six (6) liters per minute for periods up to six (6) hours.

Prescription Device

Federal Law (US) restricts this device to sale by or on the order of a physician

Prescription Use Yes OR Over-The-Counter Use No

PLEASE DO NOT WRITE BELOW THIS LINE- CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Diana R. Walker
(Division Sign-Off)
Division of Cardiovascular Devices

510(k) Number K081881