

OCT 08 2008

510(k) SUMMARY

BÂRRX Medical's HALO³⁶⁰ Energy Generator

1. Submitter's Name, Address, Telephone Number, Contact Person, and Date Prepared:

BÂRRX Medical Inc.
540 Oakmead Parkway
Sunnyvale, CA 94085

Phone: (408) 328-7302

Facsimile: (408) 328-7395

Contact Person: Viorica Filimon

Date Prepared: August 1, 2008

2. Name of device and Name/Address of Sponsor:

HALO³⁶⁰ Energy Generator

BÂRRX Medical Inc.
540 Oakmead Parkway
Sunnyvale, CA 94085

3. Common or Usual Name(s):

Electrosurgical Coagulation System

4. Classification Name:

Product code: GEI

CFR Section: 878.4400 Electrosurgical, cutting & coagulation & accessories

Device Class: II

Classification panel: General & Plastic Surgery

5. Predicate Devices

HALO³⁶⁰ Energy Generator model 1100C-115B software version 4.2 (K051168) manufactured by Stellartech Research;

Stellartech Coagulation System 2, model 1100C-115A (K050831) manufactured by Stellartech Research;

6. Intended Use / Indications for Use

The HALO³⁶⁰ Energy Generator intended use is for the coagulation of bleeding and non-bleeding sites in the gastrointestinal tract.

The HALO³⁶⁰ Energy Generator is indicated for use in the coagulation of bleeding and non-bleeding sites in the gastrointestinal tract including but not limited to, the esophagus. Indications include Esophageal Ulcers, Mallory-Weiss tears, Arteriovenous Malformations, Angiomata, Barrett's Esophagus, Dieulafoy Lesions, and Angiodysplasia.

7. Technological Characteristics

The HALO³⁶⁰ Coagulation System consists of the HALO³⁶⁰ Energy Generator model 1100C-115B with a disposable single-use HALO³⁶⁰ Coagulation Catheter, output cable, and an optional footswitch. The HALO³⁶⁰ Coagulation System performance and mode of operation is substantially equivalent to the already cleared HALO³⁶⁰ Coagulation System, and Stellartech Coagulation System 2.

HALO³⁶⁰⁺ Coagulation Catheter

There are no changes associated to the HALO³⁶⁰⁺ Coagulation Catheter.

HALO³⁶⁰ Energy Generator

The HALO³⁶⁰ Energy Generator model 1100C-115B is configured with an output cable (model CCC-001B), an optional footswitch (model FS-100B), and a power cord.

There are no changes to the hardware associated to HALO³⁶⁰ Energy Generator cleared by the 510(k) k050831 and k051168. There are minor changes implemented to the software version of the generator.

This submission addresses software changes for HALO³⁶⁰ Energy Generator cleared by the 510(k) k050831 and k051168. These minor changes implemented to the software version of the generator did not raise any new questions of safety or efficacy.

8. Substantial Equivalence

The HALO³⁶⁰ Energy Generator model 1100C-115B (software version V4.5) and the predicate devices HALO³⁶⁰ Energy Generator model 1100C-115B (software version 4.2), and Stellartech Energy Generator 2 have the same intended use, indications for use, technological characteristics, and principles of operation. The technological differences between the HALO³⁶⁰ Energy Generator model 1100C-115B software version 4.5 and its predicates are:

- Modification of the coagulation catheter operating pressure

The differences were evaluated on bench and did not raise questions regarding safety and efficacy. Thus the devices are equivalent.



OCT 08 2008

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Barrx Medical, Inc.
% Ms. Viorica Filimon
VP of Quality/Regulatory Affairs
540 Oakmead Parkway
Sunnyvale, California 94085

Re: K082202

Trade/Device Name: The HALO³⁶⁰ Energy Generator
Regulation Number: 21 CFR 878.4400
Regulation Name: Electrosurgical cutting and coagulation device and accessories
Regulatory Class: II
Product Code: GEI
Dated: September 5, 2008
Received: September 9, 2008

Dear Ms. Filimon:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K082202

g. 1 o + 1

Indications for Use Statement

510(k) Number (if known): _____

Device Name:

Indications for Use:

The HALO³⁶⁰ Energy Generator is indicated for use in the coagulation of bleeding and non-bleeding sites in the gastrointestinal tract including but not limited to, the esophagus. Indications include Esophageal Ulcers, Mallory-Weiss tears, Arteriovenous Malformations, Angiomata, Barrett's Esophagus, Dieulafoy Lesions, and Angiodysplasia.

Prescription Use X
(Part 21 C.F.R. 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 C.F.R. 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE -- CONTINUE ON ANOTHER
PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

[Signature]
(Division Sign-Off)
Division of General, Restorative,
and Neurological Devices

510(k) Number K082202

Page 1 of 1