

SEP - 2 2008

**2 Section 2**  
Summary of Safety and Effectiveness  
Prepared in accordance with 21 CFR Part 807.92(c).



GE Healthcare

General Electric Company  
P.O. Box 414, Milwaukee, WI 53201

**Section a):**

1. Submitter: GE Medical Systems Israel, Ultrasound LTD  
Haetgar Str. 4  
Tirat Carmel, ISRAEL 39120  
  
Contact Person: Israel Citron,  
Quality Assurance and Regulatory Affairs manager  
Telephone: 972-4-8519-555; Fax: 972-4-8519-500  
  
Date Prepared: July 31, 2008
2. Device Name: GE Vivid-i and Vivid q Diagnostic Ultrasound System  
Ultrasonic Pulsed Doppler Imaging System, 21 CFR 892.1550, 90-IYN  
Ultrasonic Pulsed Echo Imaging System, 21 CFR 892.1560, 90-IYO  
Diagnostic Ultrasonic Transducer, 21 CFR 892.1570, 90-ITX
3. Marketed Device: GE Vivid-i Ultrasound System, K061525 currently in commercial distribution.
4. Device Description: The GE Vivid-i and Vivid-q is compact and portable diagnostic ultrasound system with integrated keyboard, fold-up LCD type display and interchangeable electronic-array transducers. It has an overall size approximately 36 cm wide, 31.5 cm deep and 6 cm high in transport configuration and provides digital acquisition, processing and display capability. The user interface includes a computer keyboard, an intuitive layout of specialized controls, color GUI display and Doppler audio.
5. Indications for Use: The device is intended for use by a qualified physician for ultrasound evaluation of Fetal; Abdominal; Pediatric; Small Organ (breast, testes, thyroid); Neonatal Cephalic; Adult Cephalic; Cardiac (adult and pediatric); Peripheral Vascular; Musculo-skeletal Conventional and Superficial; Urology (including prostate); Transesophageal; Transrectal; Transvaginal; Intraoperative (abdominal, thoracic, and vascular), Intra-cardiac and Intra-luminal.
6. Comparison with Predicate Device: The modified GE Vivid-i is of a comparable type and substantially equivalent to the currently marketed GE Vivid-i. It is a compact and readily portable unit having the same design, construction, and materials; is comparable in key safety and effectiveness features. It has the same intended uses as the predicate device and additional software features are identical to that of other cleared GE Ultrasound systems.

**Section b):**

1. Non-clinical Tests: The device has been evaluated for acoustic output, biocompatibility, cleaning and disinfection effectiveness as well as thermal, electrical and mechanical safety, and has been found to conform with applicable medical device safety standards.
2. Clinical Tests: None required.
3. Conclusion: Intended uses and other key features are consistent with traditional clinical practice, FDA guidelines, and established methods of patient examination. The design and development process of the manufacturer conforms with 21 CFR 820, ISO 9001 and ISO13485 quality systems. The device conforms to applicable medical device safety standards and compliance is verified through independent evaluation with ongoing factory surveillance. Diagnostic ultrasound has accumulated a long history of safe and effective performance. Therefore, it is the opinion of GE Medical Systems that the GE Vivid-i and Vivid-q Diagnostic Ultrasound is substantially equivalent with respect to safety and effectiveness to devices currently cleared for market.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

GE Medical Systems, Ultrasound  
% Mr. Mark Job  
Responsible Third Party Official  
Regulatory Technology Services LLC  
1394 25<sup>th</sup> Street NW  
BUFFALO MN 55313

SEP - 2 2008

Re: K082374

Trade/Device Name: Vivid - i and Vivid - q Diagnostic Ultrasound System  
Regulation Number: 21 CFR 892.1550  
Regulation Name: Ultrasonic pulsed doppler imaging system  
Regulatory Class: II  
Product Code: IYN, IYO, and ITX  
Dated: August 15, 2008  
Received: August 18, 2008

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the Vivid - i and Vivid - q, as described in your premarket notification:

Transducer Model Number

Vivid i/q with M4S-RS	Vivid i/q with 3C-RS	Vivid i/q with 3S-RS	Vivid i/q with 9T-RS
Vivid i with 4C-RS	Vivid i/q with 8C-RS	Vivid i/q with 7S-RS	Vivid i/q with P2D
Vivid i with 12L-RS	Vivid i/q with 8L-RS	Vivid i/q with 10S-RS	Vivid i/q with P6D
Vivid i/q with 5S-RS	Vivid i/q with i12L-RS	Vivid i/q with 6T/6Tc-RS	Vivid i/q with 6S-RS
Vivid i/q with AcuNav™ 10F			

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any

Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This determination of substantial equivalence is granted on the condition that prior to shipping the first device, you submit a postclearance special report. This report should contain complete information, including acoustic output measurements based on production line devices, requested in Appendix G, (enclosed) of the Center's September 30, 1997 "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers." If the special report is incomplete or contains unacceptable values (e.g., acoustic output greater than approved levels), then the 510(k) clearance may not apply to the production units which as a result may be considered adulterated or misbranded.

The special report should reference the manufacturer's 510(k) number. It should be clearly and prominently marked "ADD-TO-FILE" and should be submitted in duplicate to:


Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
9200 Corporate Boulevard  
Rockville, Maryland 20850

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

If you have any questions regarding the content of this letter, please contact Ms. Lauren Hefner at (240) 276-3666.

Sincerely yours,

  
Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

**Diagnostic Ultrasound Indications for Use Form**  
**GE Vivid-i and Vivid-q Diagnostic Ultrasound System**

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse <sup>a</sup>	Other
Ophthalmic											
Fetal	P	P	P	P	P	P	P	P	P	P	
Abdominal <sup>[1]</sup>	P	P	P	P	P	P	P	P	P	P	
Pediatric	P	P	P	P	P	P	P	P	P	P	
Small Organ (specify) <sup>[2]</sup>	P	P	P		P	E	P	P	P	P	
Neonatal Cephalic	P	P	P	P	P	P	P	P	P	P	
Adult Cephalic	P	P	P	P	P	P	P	P	P	P	
Cardiac <sup>[3]</sup>	P	P	P	P	P	P	P	P	P	P	
Peripheral Vascular	P	P	P	P	P	E	P	P	P	P	
Musculo-skeletal Conventional	P	P	P		P	E	P	P	P	P	
Musculo-skeletal Superficial	P	P	P		P	E	P	P	P	P	
Other <sup>[4]</sup>	P	P	P		P	E	P	P	P	P	
<i>Exam Type, Means of Access</i>											
Transesophageal	P	P	P	P	P	P		P	P		
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify) <sup>[5]</sup>	P	P	P	N	P		P	P	P		
Intraoperative Neurological											
Intracardiac and Intraluminal	N	N	N	N	N	N		N	N		
Laparoscopic											

N = new indication; P = previously cleared by FDA; E = added under Appendix E

Notes: [1] Abdominal includes Renal.

[2] Small organ includes breast, testes, thyroid.

[3] Cardiac is Adult and Pediatric.

[4] Other use includes Urology.

[5] Intraoperative includes abdominal, thoracic (cardiac), and vascular (PV).


[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

[†] Coded Pulse includes Coded Octave Imaging (COI), and Coded Phase Inversion (CPI).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number R082374

## Diagnostic Ultrasound Indications for Use Form

### GE Vivid i/q with M4S-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse <sup>a</sup>	Other
Ophthalmic											
Fetal	N	N	N	N	N	N		N	N	N	
Abdominal <sup>[1]</sup>	N	N	N	N	N	N	N	N	N	N	
Pediatric	N	N	N	N	N	N	N	N	N	N	
Small Organ (specify)											
Neonatal Cephalic											
Adult Cephalic	N	N	N	N	N	N	N	N	N	N	
Cardiac <sup>[2]</sup>	N	N	N	N	N	N		N	N	N	
Peripheral Vascular											
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other (specify)											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication (transducer previously cleared with Vivid 7 system K051449); P = previously cleared by FDA; E = added under Appendix E

Notes: [1] Abdominal includes Renal

[2] Cardiac is Adult and Pediatric.


[\*] Combined modes are B/M, B/Color M, B/PWD or CWD, B/Color/PWD or CWD, B/Power/PWD.

[\*] Coded Pulse includes Coded Octave Imaging (COI), and Coded Phase Inversion (CPI).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number       K082374

**Diagnostic Ultrasound Indications for Use Form**  
**GE Vivid-i with 4C-RS Transducer**

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse <sup>1</sup>	Other
Ophthalmic											
Fetal	P	P	P		P	E	P	P	P	P	
Abdominal <sup>[1]</sup>	P	P	P		P	E	P	P	P	P	
Pediatric											
Small Organ (specify)											
Neonatal Cephalic											
Adult Cephalic											
Cardiac											
Peripheral Vascular											
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other <sup>[2]</sup>	P	P	P		P	E	P	P	P	P	
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K061525); E = added under Appendix E

Notes: [1] Abdominal includes Renal;

[2] Other use includes Urology.

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

[\*] Coded Pulse includes Coded Octave Imaging (COI), and Coded Phase Inversion (CPI).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number K082374

## Diagnostic Ultrasound Indications for Use Form

### GE Vivid-i with 12L-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse <sup>a</sup>	Other
Ophthalmic											
Fetal											
Abdominal											
Pediatric											
Small Organ <sup>[1]</sup>	P	P	P		P	E	P	P	P	P	
Neonatal Cephalic											
Adult Cephalic											
Cardiac											
Peripheral Vascular	P	P	P		P	E	P	P	P	P	
Musculo-skeletal Conventional	P	P	P		P	E	P	P	P	P	
Musculo-skeletal Superficial	P	P	P		P	E	P	P	P	P	
Other (specify)											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K061525); E = added under Appendix E

Notes: [1] Small organ includes breast, testes, thyroid.


[\*] Combined modes are B/M, B/PWD, B/Color/PWD, B/Power/PWD.

[<sup>a</sup>] Coded Pulse includes Coded Octave Imaging (COI), and Coded Phase Inversion (CPI).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number K082374

## Diagnostic Ultrasound Indications for Use Form

### GE Vivid-i/q with 5S-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse <sup>a</sup>	Other
Ophthalmic											
Fetal	P	P	P	P	P	P		P	P	P	
Abdominal											
Pediatric	P	P	P	P	P	P	P	P	P	P	
Small Organ (specify)											
Neonatal Cephalic											
Adult Cephalic											
Cardiac <sup>[1]</sup>	P	P	P	P	P	P		P	P	P	
Peripheral Vascular											
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K061525); E = added under Appendix E

Notes: [1] Cardiac is Adult and Pediatric.


[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

[\*] Coded Pulse includes Coded Octave Imaging (COI), and Coded Phase Inversion (CPI).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K082374

## Diagnostic Ultrasound Indications for Use Form

### GE Vivid-i/q with 3C-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse <sup>a</sup>	Other
Ophthalmic											
Fetal	P	P	P		P	E		P	P	P	
Abdominal <sup>[1]</sup>	P	P	P		P	E	P	P	P	P	
Pediatric	P	P	P		P	E	P	P	P	P	
Small Organ (specify)											
Neonatal Cephalic											
Adult Cephalic											
Cardiac											
Peripheral Vascular											
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other <sup>[2]</sup>	P	P	P		P	E	P	P	P	P	
Exam Type, Means of Access											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K033139); E = added under Appendix E

Notes: [1] Abdominal includes Renal.

[2] Other use includes Urology.

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

[†] Coded Pulse includes Coded Octave Imaging (COI), and Coded Phase Inversion (CPI).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

**K082374**

### Diagnostic Ultrasound Indications for Use Form

#### GE Vivid-i/q with 8C-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse*	Other
Ophthalmic											
Fetal											
Abdominal	P	P	P		P	E	P	P	P	P	
Pediatric	P	P	P		P	E	P	P	P	P	
Small Organ (specify) <sup>[1]</sup>	P	P	P		P	E	P	P	P	P	
Neonatal Cephalic	P	P	P		P	E	P	P	P	P	
Adult Cephalic											
Cardiac											
Peripheral Vascular	P	P	P		P	E	P	P	P	P	
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K033139); E = added under Appendix E

Notes: [1] Small organ includes breast, testes, thyroid.

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

[\*] Coded Pulse includes Coded Octave Imaging (COI), and Coded Phase Inversion (CPI).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*[Signature]*

K882374

### Diagnostic Ultrasound Indications for Use Form

#### GE Vivid i/q with 8L-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse*	Other
Ophthalmic											
Fetal											
Abdominal	P	P	P		P	E	P	P	P	P	
Pediatric	P	P	P		P	E	P	P	P	P	
Small Organ (specify) <sup>[1]</sup>	P	P	P		P	E	P	P	P	P	
Neonatal Cephalic											
Adult Cephalic											
Cardiac											
Peripheral Vascular	P	P	P		P	E	P	P	P	P	
Musculo-skeletal Conventional	P	P	P		P	E	P	P	P	P	
Musculo-skeletal Superficial	P	P	P		P	E	P	P	P	P	
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K033139); E = added under Appendix E

Notes: [1] Small organ includes breast, testes, thyroid.

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

[\*] Coded Pulse includes Coded Octave Imaging (COI), and Coded Phase Inversion (CPI).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


Prescription User (Per 21 CFR 801.109)

A-9

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

  
K082374

## Diagnostic Ultrasound Indications for Use Form

### GE Vivid i/q with i12L-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse	Other
Ophthalmic											
Fetal											
Abdominal	P	P	P		P	E	P	P	P		
Pediatric	P	P	P		P	E	P	P	P		
Small Organ (specify) <sup>[1]</sup>	P	P	P		P	E	P	P	P		
Neonatal Cephalic											
Adult Cephalic											
Cardiac <sup>[2]</sup>	P	P	P		P	E		P	P		
Peripheral Vascular	P	P	P		P	E	P	P	P		
Musculo-skeletal Conventional	P	P	P		P	E	P	P	P		
Musculo-skeletal Superficial	P	P	P		P	E	P	P	P		
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify) <sup>[3]</sup>	P	P	P		P	E	P	P	P		
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K033139); E = added under Appendix E

Notes: [1] Small organ includes breast, testes, thyroid.

[2] Cardiac is Adult and Pediatric.

[3] Intraoperative includes abdominal, thoracic (cardiac), and vascular (PV).

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*[Signature]*  
K082374

### Diagnostic Ultrasound Indications for Use Form

#### GE Vivid i/q with 3S-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse <sup>†</sup>	Other
Ophthalmic											
Fetal	P	P	P	P	P	P		P	P	P	
Abdominal <sup>[1]</sup>	P	P	P	P	P	P	P	P	P	P	
Pediatric	P	P	P	P	P	P	P	P	P	P	
Small Organ (specify)											
Neonatal Cephalic											
Adult Cephalic	P	P	P	P	P	P	P	P	P	P	
Cardiac <sup>[2]</sup>	P	P	P	P	P	P		P	P	P	
Peripheral Vascular											
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new Indication; P = previously cleared by FDA (K033139); E = added under Appendix E

Notes: [1] Abdominal includes Renal.

[2] Cardiac is Adult and Pediatric.


[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

[†] Coded Pulse includes Coded Octave Imaging (COI), and Coded Phase Inversion (CPI).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K082374

## Diagnostic Ultrasound Indications for Use Form

### GE Vivid i/q with 7S-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse	Other
Ophthalmic											
Fetal	E	E	E	E	E	E		E	E		
Abdominal	P	P	P	P	P	P	P	P	P		
Pediatric	P	P	P	P	P	P	P	P	P		
Small Organ (specify)											
Neonatal Cephalic	P	P	P	P	P	P	P	P	P		
Adult Cephalic											
Cardiac <sup>[1]</sup>	P	P	P	P	P	P	P	P	P		
Peripheral Vascular	E	E	E	E	E	E	E	E	E		
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K033139); E = added under Appendix E

Notes: [1] Cardiac is Adult and Pediatric.

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*[Signature]*

K082374

## Diagnostic Ultrasound Indications for Use Form

### GE Vivid i/q with 10S-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse	Other
Ophthalmic											
Fetal											
Abdominal	P	P	P	P	P	P	P	P	P		
Pediatric	P	P	P	P	P	P	P	P	P		
Small Organ (specify)											
Neonatal Cephalic	P	P	P	P	P	P	P	P	P		
Adult Cephalic											
Cardiac <sup>[1]</sup>	P	P	P	P	P	P		P	P		
Peripheral Vascular	E	E	E	E	E	E	E	E	E		
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K033139); E = added under Appendix E

Notes: [1] Cardiac is Adult and Pediatric

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*[Signature]*  
K082374

## Diagnostic Ultrasound Indications for Use Form

### GE Vivid i/q with 6T/6Tc-RS Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse	Other
Ophthalmic											
Fetal											
Abdominal											
Pediatric											
Small Organ (specify)											
Neonatal Cephalic											
Adult Cephalic											
Cardiac <sup>[1]</sup>	P	P	P	P	P	P		P	P		
Peripheral Vascular											
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal	P	P	P	P	P	P		P	P		
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K033139); E = added under Appendix E


Notes: [1] Cardiac is Adult and Pediatric.

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K082374

## Diagnostic Ultrasound Indications for Use Form

### GE Vivid i/q with 9T-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse	Other
Ophthalmic											
Fetal											
Abdominal											
Pediatric											
Small Organ (specify)											
Neonatal Cephalic											
Adult Cephalic											
Cardiac <sup>[1]</sup>	P	P	P	P	P	P		P	P		
Peripheral Vascular											
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal	P	P	P	P	P	P		P	P		
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K033139); E = added under Appendix E

Notes: [1] Cardiac is Adult and Pediatric.

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*[Signature]*  
K082374

## Diagnostic Ultrasound Indications for Use Form

### GE Vivid i/q with P2D Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse	Other
Ophthalmic											
Fetal / Obstetrics											
Abdominal											
Pediatric											
Small Organ (specify)											
Neonatal Cephalic											
Adult Cephalic											
Cardiac <sup>[1]</sup>			P	P							
Peripheral Vascular			P	P							
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K033139); E = added under Appendix E

Notes: [1] Cardiac is Adult and Pediatric.

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number



K082374

**Diagnostic Ultrasound Indications for Use Form**  
**GE Vivid i/q with P6D Transducer**

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse	Other
Ophthalmic											
Fetal / Obstetrics											
Abdominal											
Pediatric											
Small Organ (specify)											
Neonatal Cephalic											
Adult Cephalic											
Cardiac <sup>[1]</sup>			P	P							
Peripheral Vascular			P	P							
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication; P = previously cleared by FDA (K033139); E = added under Appendix E

Notes: [1] Cardiac is Adult and Pediatric.

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K082374

## Diagnostic Ultrasound Indications for Use Form

### GE Vivid i/q with 6S-RS Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse	Other
Ophthalmic											
Fetal	N	N	N	N	N	N	N	N	N		
Abdominal	N	N	N	N	N	N	N	N	N		
Pediatric	N	N	N	N	N	N	N	N	N		
Small Organ (specify)											
Neonatal Cephalic	N	N	N	N	N	N	N	N	N		
Adult Cephalic											
Cardiac <sup>[1]</sup>	N	N	N	N	N	N	N	N	N		
Peripheral Vascular											
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal											
Laparoscopic											

N = new indication (transducer previously cleared with LOGIQ-i/e and Vivid-e systems K072797); P = previously cleared by FDA; E = added under Appendix E

Notes: [1] Cardiac is Adult and Pediatric.

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number



K082374

### Diagnostic Ultrasound Indications for Use Form

#### GE Vivid i/q with AcuNav™ 10F Transducer

Intended Use: Ultrasound imaging, measurement and analysis of the human body as follows:

Clinical Application <i>Anatomy/Region of Interest</i>	Mode of Operation										
	B	M	PW Doppler	CW Doppler	Color Doppler	Color M Doppler	Power Doppler	Combined Modes	Harmonic Imaging	Coded Pulse	Other
Ophthalmic											
Fetal											
Abdominal											
Pediatric											
Small Organ (specify)											
Neonatal Cephalic											
Adult Cephalic											
Cardiac <sup>[1]</sup>	N	N	N	N	N	N		N	N		
Peripheral Vascular											
Musculo-skeletal Conventional											
Musculo-skeletal Superficial											
Other											
<i>Exam Type, Means of Access</i>											
Transesophageal											
Transrectal											
Transvaginal											
Transurethral											
Intraoperative (specify)											
Intraoperative Neurological											
Intracardiac and Intraluminal	N	N	N	N	N	N		N	N		
Laparoscopic											

N = new indication (transducer previously cleared as a stand alone medical device K033650); P = previously cleared by FDA; E = added under Appendix E

Notes: [1] Cardiac is Adult and Pediatric.

[\*] Combined modes are B/M, B/PWD, B/CWD, B/Color/PWD, B/Amplitude/PWD.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription User (Per 21 CFR 801.109)

(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*[Signature]*

K082374