

K082650

3.0 510(k) SummaryPage 1 of 1

Sponsor: Synthes (USA)
1301 Goshen Parkway
West Chester, PA 19380
(610) 719-5000

Contact: Jill R. Sherman
Synthes (USA)
1301 Goshen Parkway
West Chester, PA 19380
Phone: 610-719-6538 Fax: 484-356-9682

Device Name: Synthes Large External Fixation, MR Conditional

Classification: 21 CFR Part 888.3030; Single/multiple component metallic bone fixation appliances and accessories.

Predicate Devices: Synthes External Fixation Devices

Device Description: The Synthes Large External Fixation, MR Conditional consists of previously cleared clamps, rods, Schanz screws and Steinmann pins. Synthes will also offer Sterile Packaged Large External Fixation Kits, MR Conditional. These kits contain previously cleared external fixation devices and manual surgical instruments.

Intended Use: Synthes Large External Fixation is intended for use to provide treatment for long bone and pelvic fractures that require external fixation. Specifically, the components can be used for:

- Stabilization of soft tissues and fractures
- Polytrauma/multiple orthopedic trauma
- Vertically stable pelvic fractures, or as a treatment adjunct for vertically unstable pelvic fractures
- Arthrodeses and osteotomies with soft tissue problems; failures of total joints
- Neutralization of fractures stabilized with limited internal fixation
- Non-unions/septic non-unions
- Intra-operative reductions/stabilization tool to assist with indirect reduction
- Unilateral rectilinear bone segment transport or leg lengthening

Substantial Equivalence Information presented supports substantial equivalence.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Synthes (USA)
% Ms. Jill R. Sherman
1301 Goshen Parkway
West Chester, PA 19380

NOV 18 2008

Re: K082650

Trade/Device Name: Synthes Large External Fixation, MR Conditional
Regulation Number: 21 CFR 888.3030
Regulation Name: Single/multiple component metallic bone fixation appliances and accessories.
Regulatory Class: II
Product Code: KTT
Dated: September 11, 2008
Received: September 12, 2008

Dear Ms. Sherman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at toll-free number (800) 638-2041 or (240) 276-3150 or the Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

2.0

Indications for Use

510(k) Number (if known): K082650Device Name: Synthes Large External Fixation, MR Conditional

Indications for Use:

Synthes Large External Fixation is intended for use to provide treatment for long bone and pelvic fractures that require external fixation. Specifically, the components can be used for:


- Stabilization of soft tissues and fractures
- Polytrauma/multiple orthopedic trauma
- Vertically stable pelvic fractures, or as a treatment adjunct for vertically unstable pelvic fractures
- Arthrodeses and osteotomies with soft tissue problems; failures of total joints
- Neutralization of fractures stabilized with limited internal fixation
- Non-unions/septic non-unions
- Intra-operative reductions/stabilization tool to assist with indirect reduction
- Unilateral rectilinear bone segment transport or leg lengthening

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

AND/OR Over-The-Counter Use _____
(21 CFR 807 Subpart C)


(Division Sign-Off) Page 1 of 1

Division of General, Restorative,
and Neurological Devices

510(k) Number K082650