

1003254

DEC 19 2008

**510(k) Summary of Safety and Effectiveness**  
**Prepared in accordance with 21 CFR Part 807.92**

**Section a):**

1. Submitter: Aloka Co., Ltd., 10 Fairfield Boulevard, Wallingford, CT 06492  
  
Contact Person: Richard J. Cehovsky, RA/QA Coordinator,  
Tel: (203)269-5088 Ext. 346, Fax: 203-269-6075  
  
Date Prepared: 9/18/08
2. Device Name: Aloka Prosound 2 Diagnostic Ultrasound System  
Ultrasonic Pulsed Doppler Imaging System, 21 CFR 892.1550, 90 IYN  
Diagnostic Ultrasound Transducer, 21 CFR 892.1570, 90 ITX  
Ultrasonic Pulsed Echo Imaging System., 21 CFR 892.1560, 90 IYO
3. Marketed Device: Aloka SSD-500 Diagnostic Ultrasound System K900805, (90-IYN, ITX, IYO)  
( A device currently in commercial distribution)
4. Device Description: The Aloka Prosound 2 Diagnostic Ultrasound System is a light weight, full-digital portable imaging and analysis system. It consist of a high resolution LCD flat panel monitor that provides excellent image quality and processing . The user interface includes a computer type keyboard, specialized controls and a display.
5. Indications for Use: The device is intended for use by a qualified physician for ultrasound evaluation of Gyneological, Fetal, Peripheral Vascular, Cardiac, Neonatal Cephalic, Small Parts, Intra-operative, Transrectal and Abdominal applications. The device is not indicated for Ophthalmic applications.
6. Comparison w/ Predicate Device:  
The Aloka Prosound 2 is technically comparable and substantially equivalent to the current Aloka SSD-500-(K900805). It has the same technological characteristics, key safety and effectiveness features, and has the same intended uses and basic operating modes as the predicate device.

**Section b):**

1. Non-clinical Tests: The device and its transducers have been evaluated for acoustic output, biocompatibility, cleaning & disinfection effectiveness, electromagnetic compatibility, as well as electrical and mechanical safety, and have been found to conform with applicable medical device safety standards.
2. Clinical Tests: None Required.
3. Conclusion: Intended uses and other key features are consistent with traditional clinical practices, FDA guidelines and established methods of patient examination. The design, development and quality process of the manufacturer confirms with 21 CFR 820, ISO 9001:2000 and ISO 13485 quality systems. The device conforms to applicable medical device safety standards and compliance is verified through independent evaluation with ongoing factory surveillance. Diagnostic ultrasound has accumulated a long history of safe and effectiveness performance. Therefore, it is the opinion of Aloka Co., Ltd. that the Aloka Prosound 2 Diagnostic Ultrasound System and its transducers are substantially equivalent with respect to safety and effectiveness to its predicate and other currently cleared Aloka systems.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JUN 2 2009

Aloka Co., Ltd.  
% Mr. Tamas Borsai  
Division Manager, Medical Division  
TÜV Rheinland of North America  
12 Commerce Road  
NEWTOWN CT 06470

Re: K083254  
Trade/Device Name: Aloka Prosound 2  
Regulation Number: 21 CFR 892.1550  
Regulation Name: Ultrasonic pulsed doppler imaging system  
Regulatory Class: II  
Product Code: IYN, IYO, and ITX  
Dated: December 2, 2008  
Received: December 4, 2008

Dear Mr. Borsai:

This letter corrects our substantially equivalent letter of December 19, 2008.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the Aloka Prosound 2, as described in your premarket notification:

Transducer Model Number

UST-556I-7.5  
UST-556T-7.5  
UST-556TU-7.5  
UST-586-5  
UST587I-5

UST-660-7.5  
UST-934N-3.5  
UST-935N-5  
UST-944B-3.5  
UST-945B-5

UST-974-5  
UST-981-5  
UST-5512U-7.5  
UST-5551

UST-5711  
UST-5820-5  
UST-9111-5  
UST-9137

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

If you have any questions regarding the content of this letter, please contact Lauren Hefner at (240) 276-3666.

Sincerely yours,



Janine M. Morris  
Acting Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure(s)

**Section 1.3**  
**Indications for Use**

**Section 1.3.1**  
**510(k) Indications for Use Statement/Forms**

( Immediately follows this page)

## Indications for Use

**510(K) Number (if known):**

**Device Name:** Aloka Prosound 2

**Indications For Use:**

The device is intended for use by a qualified physician for ultrasound evaluation of Gynecological, Fetal, Peripheral Vascular, Cardiac, Neonatal Cephalic, Small Parts, Intra-operative, Transrectal and Abdominal applications.

The device is not indicated for Ophthalmic applications.

**Prescription Use** ✓  
(Part 21 CFR 801 Subpart D)


**AND/OR**

**Over-The Counter Use** \_\_\_\_\_  
(21 CFR 801 Subpart C)

**(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)**

---

**Concurrence of CDRH, Office of Device Evaluation (ODE)**

  
\_\_\_\_\_  
(Division Sign-Off)

Division of Reproductive, Abdominal and  
Radiological Devices

510(k) Number K083254

**Page 1 of 1**

## 1.3.1

**Diagnostic Ultrasound Indications for Use Form**  
**Prosound 2**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal		N	N						See Below	
Abdominal		N	N						See Below	
Intraoperative (specify)		N	N						See Below	
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		N	N						See Below	
Neonatal Cephalic		N	N						See Below	
Adult Cephalic										
Cardiac		N	N						See Below	
Transesophageal										
Transrectal		N	N						See Below	
Transvaginal		N	N						See Below	
Transurethral										
Intravascular										
Peripheral Vascular		N	N						See Below	
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other: Gynecological		N	N						See Below	

N= new indication; P= previously cleared by FDA; E= added under Appendix A

Additional Comments: Mixed mode operation includes- B/M

Applications: Small Parts-(breast, testes & thyroid), abdominal, gynecological-fetal, neonatal cephalic, transrectal, transvaginal, cardiac, intra-operative-(liver, pancreas, gall bladder, etc.), peripheral vascular.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
 Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number 6083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-556I-7.5**  
**(K840540)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal										
Abdominal										
Intraoperative (specify)		P	P						See Below	
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		P	P						See Below	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

Applications: Small Parts- (breasts, tests, thyroid, etc.), Intra-Operative- (Liver, pancreas, gall Bladder)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number K083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-556T-7.5**  
**(K840540)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal										
Abdominal										
Intraoperative (specify)		P	P						Sec Below	
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		P	P						Sec Below	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other:										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

Applications: Small Parts- (breasts, tests, thyroid, etc.), Intra-Operative- (Liver, pancreas, gall Bladder)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal and  
Radiological Devices

510(k) Number K083254



**Diagnostic Ultrasound Indications for Use Form**  
**UST-556TU-7.5**  
**(K870916)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal										
Abdominal										
Intraoperative (specify)		P	P						See Below	
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		P	P						See Below	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

Applications: Small Parts- (breasts, tests, thyroid, etc.), Intra-Operative- (Liver, pancreas, gall Bladder)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
 Prescription Use (Per 21 CFR 801.109)

  
 (Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number K083254

## Diagnostic Ultrasound Indications for Use Form

UST-586-5

(K861538)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		N	N						See Below	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E


Additional Comments: Mixed mode operation includes- B/M

Applications: Small Parts- (breasts, tests, thyroid, etc.)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)



(Division Sign-Off)

Division of Reproductive, Abdominal and Radiological Devices

510(k) Number 15083254

## Diagnostic Ultrasound Indications for Use Form

UST-587I-5

(K840540)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal										
Abdominal										
Intraoperative (specify)		P	P						See Below	
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		P	P						See Below	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

Applications: Small Parts- (breasts, tests, thyroid, etc.), Intra-Operative- (Liver, pancreas, gall Bladder)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal and  
Radiological Devices

510(k) Number

K083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-660-7.5**  
**(K870916)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal		P	P						See Below	
Transvaginal-										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										


N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

  
 (Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number K083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-934N-3.5**  
**(K900805)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		P	P						See Below	
Abdominal		P	P						See Below	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other: Gynecological		P	P						See Below	

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
 Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number K083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-935N-5**  
**(K900805)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:


Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal		P	P						See Below	
Abdominal		P	P						See Below	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other: Gynecological		P	P						See Below	

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
 Prescription Use (Per 21 CFR 801.109)

  
 (Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number 4083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-944B-3.5**  
**(K870916)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:


Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal										
Abdominal		P	P						See Below	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		P	P						See Below	
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other:										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

**(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)**

Concurrence of CDRH, Office of Device Evaluation (ODE)  
 Prescription Use (Per 21 CFR 801.109)



(Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number K083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-945B-5**  
**(K900805)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		P	P						See Below	
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal		P	P						See Below	
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other: Gynecological		P	P						See Below	

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

  
 (Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number 1083254



## Diagnostic Ultrasound Indications for Use Form

UST-974-5

(K910153)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic		P	P						See Below	
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other:										


N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)



(Division Sign-Off)  
 Division of Reproductive, Abdominal and  
 Radiological Devices  
 510(k) Number 14083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-981-5**  
**(K900805)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal		N	N						See Below	
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal		N	N						See Below	
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other: Gynecological		N	N						See Below	

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

**(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)**

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number

14083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-5512U-7.5**  
**(K861538)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		P	P						See Below	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		P	P						See Below	
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other:										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

Applications: Small Parts- (breasts, tests, thyroid, etc.)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
 Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number K083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-5551**  
**(K861538)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		N	N						See Below	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other:										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

Applications: Small Parts- (breasts, tests, thyroid, etc.)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
 Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number 14083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-5711**  
**(K861538)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		N	N						See Below	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		N	N						See Below	
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other:										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

Applications: Small Parts- (breasts, tests, thyroid, etc.)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
 Prescription Use (Per 21 CFR 801.109)

  
 (Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number 6083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-5820-5**  
**(K840540)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)		N	N						See Below	
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		N	N						See Below	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other:										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

Applications: Small Parts- (breasts, tests, thyroid, etc.)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
 Prescription Use (Per 21 CFR 801.109)

  
 (Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number K083254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-9111-5**  
**(K910153)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal										
Abdominal										
Intraoperative (specify)		N	N						See Below	
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic		N	N						See Below	
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other:										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

Applications: Intra-Operative- (Liver, pancreas, gall Bladder)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal and  
Radiological Devices

510(k) Number KC83254

**Diagnostic Ultrasound Indications for Use Form**  
**UST-9137**  
**(K900805)**

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Modes of operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Opthalmic										
Fetal		N	N						See Below	
Abdominal		N	N						See Below	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other: Gynecological		N	N						See Below	

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Mixed mode operation includes- B/M

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal and  
 Radiological Devices

510(k) Number 14083254



**1.3.2 New Indications for Use:** There are no new indications for use.

**1.3.3 Previously Cleared Indications for Use:**

The Prosound 2 diagnostic ultrasound system and its transducers are intended for use in diagnostic ultrasound examinations. These ultrasound applications and indications for use include: Gyneological, Fetal, Peripheral Vascular, Cardiac, Neonatal Cephalic, Small Parts, Intra-operative, Transrectal and Abdominal applications.

The Prosound 2 has the same indications for use as its predicate – Aloka SSD-500 (K910153) and other market cleared systems and transducers manufactured by Aloka.

Promotional information for the Prosound 2 is provided in Appendix D of this submission.