

**510(k) Summary**

[in accordance with SMDA of 1990, 21 CFR 807 92(c)]

JAN 22 2009

**Contact** Mr Hartmut Loch  
Regulatory Consultant  
Allez Spine, LLC  
2301 Dupont Drive, Suite 510  
Irvine CA 92612

**Trade name** Allez Spine *Laguna*® Polyaxial Pedicle Screw System

**Common name** Spinal Fixation System

**Classification name** Appliance, Fixation, Spinal Interlaminar - § 888 3050 (KWP)  
Appliance, Fixation, Spinal Intervertebral - § 888 3060 (KWQ)  
Orthosis, Spinal Pedicle Fixation - § 888 3070 (MNI)  
Orthosis, Spondylolisthesis Spinal Fixation - § 888 3070 (MNH)  
All Class II, Orthopedic Device-Panel 87

**Product Code** KWP, KWQ, MNI & MNH

**Device Description and Characteristics** The *Laguna*® Polyaxial Pedicle Screw System is intended to help provide correction, immobilization and stabilization of spinal segments as an adjunct to fusion of the thoracic, lumbar and/or sacral space

The LAGUNA® Polyaxial Pedicle Screw System consists of a variety of rods and screws, which can be rigidly locked into a variety of configurations, with each construct being tailor made for the individual case. Multi axial screws are supplied in winged and non-winged configurations, in a variety of different length, and in 5mm, 6mm, and 7mm diameter sizes. All sizes are able to receive 5.5mm connecting rods only.

The LAGUNA® Polyaxial Pedicle Screw System implant components are fabricated from medical grade titanium alloy per ASTM F136.

**Equivalence** The modified *Laguna*® Polyaxial Pedicle Screw System is substantially equivalent to the original *Laguna*® Pedicle Screw System, which is manufactured and marketed by Allez Spine, LLC.

**Indications** The LAGUNA® Polyaxial Pedicle Screw System is intended for posterior, non-cervical fixation for the following conditions: degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies), spondylolisthesis, trauma (i.e. fracture or dislocation), spinal stenosis, curvatures (i.e. scoliosis, kyphosis and/or lordosis), tumor pseudarthrosis, and/or failed previous fusion.

**Performance data** Biomechanical tests have been performed. The test results were equivalent to other similar implants and are sufficient for *in vivo* loading.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Allez Spine, LLC  
% Mr Hartmut Loch  
2301 Dupont Drive, Suite 510  
Irvine, CA 92612

JAN 22 2009

Re K083826  
Trade/Device Name Laguna Polyaxial Pedicle Screw System  
Regulation Number 21 CFR 888 3070  
Regulation Name Pedicle Screw Spinal System  
Regulatory Class II  
Product Code MNH, MNI, KWP, KWQ  
Dated December 16, 2008  
Received December 23, 2008

Dear Mr Loch

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to, registration and listing (21 CFR Part 807), labeling (21 CFR Part 801), good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820), and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act), 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at toll-free number (800) 638-2041 or (240) 276-3150 or the Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

Sincerely yours,



Mark N. Melkerson  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known) K083826

Device Name(s) LAGUNA® Polyaxial Pedicle Screw System

### Indications for Use

The LAGUNA® Polyaxial Pedicle Screw System is intended for posterior, non-cervical fixation for the following conditions degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies), spondylolisthesis; trauma (i.e. fracture or dislocation), spinal stenosis, curvatures (i.e. scoliosis, kyphosis and/or lordosis), tumor pseudarthrosis, and/or failed previous fusion

Prescription Use ✓ AND/OR Over-The-Counter-Use \_\_\_\_\_

(Part 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

*Neil R. Ogle for MRM*  
(Division Sign-Off)

Division of General, Restorative,  
and Neurological Devices

510(k) Number K083826

Page 1 of 1