

SECTION 5 - 510(k) Summary or 510(k) Statement

510(k) Summary

(As Required by Section 807.92 (c))

1. Submitter

Name:

Unicare Biomedical

Address:

25951 La Cuesta Avenue, Laguna Hills, CA 92653

Contact:

Stan Yang, 949-643-6707

Date:

January 7, 2009

2. Device Name

Trade Name:

CytoFlex® Resorb

Common Name:

Absorbable barrier membrane

Classification Name:

Barrier, synthetic, intraoral

Device Classification:

II

3. Predicate Devices

Gore Resolut (K973594; W.L. Gore & Associates)

Biomesh (K990363; Sam Yang Corporation)

Bioscaff Alvelac (K080308; Bioscaffold International)

4. Device Description

CytoFlex® Resorb barriers are made from polyglycolide, polylactide and poly(glycolide-co-lactide) copolymer. CytoFlex® Resorb is a resorbable barrier membrane and is supplied in a variety of shapes and sizes in sealed pouches. CytoFlex® Resorb membrane is tested, evaluated and found to be substantially equivalent to legally marketed predicate devices.

5. Indication

CytoFlex® Resorb membranes are intended for use as a space-making barrier in the treatment of periodontal defects and maxillofacial guided tissue regeneration procedures, including preservation and regeneration of alveolar bone height and volume, ridge and extraction site augmentation, sinus lifts, and treatment of associated cystic defects. It is also intended for use as a grafting material containment matrix.

DEPARTMENT OF HEALTH & HUMAN SERVICES



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

MAR 2 3 2009

Mr. Stan Yang Vice President Unicare Biomedical, Incorporated 22971Triton Way, Unit B Laguna Hills, California 92653

Re: K090083

Trade/Device Name: CytoFlex® Resorb Regulation Number: 21 CFR 872.3930 Regulation Name: Bone Grafting Material

Regulatory Class: II Product Code: LYC Dated: January 7, 2009 Received: January 12, 2009

Dear Mr. Yang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at 240-276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at 240-276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Ginette Y. Michaud, M.D.

Acting Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices

note 4. Michan mo

Office of Device Evaluation

Center for Devices and Radiological Health

SECTION 4 - Indications for Use Statement

510(k) Number (if known):

Device Name: CytoFlex® Resorb

Indications for Use

Indications for Use:			
CytoFlex® Resorb membranes are interestment of periodontal defects and mincluding preservation and regeneration extraction site augmentation, sinus lift also intended for use as a grafting mat	naxillofacial g on of alveolar ts, and treatm	guided tissue regeneration probone height and volume, rident of associated cystic defec	ocedures, ge and
	•		
·			
Prescription Use (Part 21 CFR 801 Subpart D)	AND/OR	Over-The-Counter Use(21 CFR 801 Subpart C)	 -
(PLEASE DO NOT WRITE BELOW THIS	S LINE-CONTI	NUE ON ANOTHER PAGE OF N	TEEDED)
	(Division Si Division of Infection Co	gn-Off) Anesthesiology, General Hospitz Ontrol, Dental Devices	1
Unicare Biomedical 510 (k) Notification	510(k) Num	Confi	dential
the first of the second second	* ;* ****,*		