

**510(k) SUMMARY**

DEC 29 2009

**V.A.C.® Therapy Systems**

<b>Date prepared</b>	December 23, 2009
<b>510(k) owner</b>	
<b>Name</b>	KCI USA, Inc. (Kinetic Concepts, Inc.)
<b>Address</b>	6203 Farinon Drive; San Antonio, Texas 78249
<b>Fax number</b>	210 255-6727
<b>Name of contact person</b>	Margaret Marsh
<b>Name of the device</b>	
<b>Trade or proprietary name</b>	V.A.C.® Therapy Systems
<b>Common or usual name</b>	Negative Pressure Wound Therapy Systems
<b>Classification name</b>	Negative Pressure Wound Therapy Powered Suction Pump (and components)
<b>Legally marketed device(s) to which equivalence is claimed</b>	V.A.C.® Therapy System Family of Products
<b>Device description</b>	The software-controlled therapy unit applies negative pressure to the wound bed. The open cells of the V.A.C.® Foam Dressing to which the therapy unit is connected enable distribution of the negative pressure across the surface of the wound, while the tubing transfers accumulated fluids to the canister. The software monitors and maintains target pressure and alarms as needed to help assure target pressure is maintained and constant therapy is delivered.

<b>Intended use of the device</b>	<p>The ActiV.A.C.<sup>®</sup>, InfoV.A.C.<sup>®</sup>, V.A.C. ATS<sup>®</sup>, and V.A.C. Freedom<sup>®</sup> Therapy Systems are integrated wound management systems for use in acute, extended and home care settings. They are intended to create an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudate and infectious material. They are indicated for patients with chronic, acute, traumatic, subacute and dehiscent wounds, partial-thickness burns, ulcers (such as diabetic, pressure or venous insufficiency), flaps and grafts.</p> <p>The V.A.C. GranuFoam Silver<sup>®</sup> Dressing is an effective barrier to bacterial penetration and may help reduce infection in the above wound types.</p> <p>The V.A.C. Instill<sup>®</sup> Therapy System is indicated for patients who would benefit from V.A.C.<sup>®</sup> Therapy coupled with controlled delivery of topical wound treatment solutions and suspensions over the wound bed. It is intended for patients with chronic acute, traumatic, subacute and dehiscent wounds, partial thickness burns, ulcers (such as diabetic, pressure, or venous insufficiency), flaps and grafts.</p>
<b>Differences in intended use from the predicate(s)</b>	<p>The only difference between the current and the predicate products is in the addition of venous insufficiency ulcers to the Indications for Use statement.</p>
<b>Summary of the technological characteristics of the device compared to the predicate device</b>	<p>The technological characteristics of the V.A.C.<sup>®</sup> Therapy System Family of Products is unchanged.</p>
<b>Summary of nonclinical tests</b>	<p>None required for the change in indication</p>
<b>Summary of clinical tests</b>	<p>Medical literature documents acceptable clinical experience using the V.A.C.<sup>®</sup> Therapy System to treat venous insufficiency ulcers.</p>



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-066-0609  
Silver Spring, MD 20993-0002

KCI USA  
% Ms. Margaret Marsh  
Regulatory Affairs Technical Director  
6203 Farinon Drive  
San Antonio, Texas 78249

DEC 29 2009

Re: K091585  
Trade/Device Name: V.A.C.<sup>®</sup> Therapy System; V.A.C. Instillamet (V.A.C.<sup>®</sup> Instill)  
Regulation Number: 21 CFR 878-4780  
Regulation Name: Powered suction pump  
Regulatory Class: Class II  
Product Code: OMP  
Dated: October 12, 2009  
Received: October 13, 2009

Dear Ms. Marsh:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical


Page 2 - Ms. Margaret Marsh

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/cdrh/mdr/> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

*for*  *ms mark*  
Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K091585  
pg 1 of 2

## INDICATIONS FOR USE

510(k) Number (if known): K091585

Device Name: V.A.C.® Therapy System

Indications for Use:

The ActiV.A.C.®, InfoV.A.C.®, V.A.C. ATS®, and V.A.C. Freedom® Therapy Systems are integrated wound management systems for use in acute, extended and home care settings. They are intended to create an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudate and infectious material. They are indicated for patients with chronic, acute, traumatic, subacute and dehiscent wounds, partial-thickness burns, ulcers (such as diabetic, pressure or venous insufficiency), flaps and grafts.

The V.A.C. GranuFoam Silver® Dressing is an effective barrier to bacterial penetration and may help reduce infection in the above wound types.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR Over-The-Counter Use             
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Posted November 13, 2003)

David Krane for MxM  
(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

Page   1   of   1  

510(k) Number   K091585

K091585  
pg 2 of 2

## INDICATIONS FOR USE

510(k) Number (if known):

K091585

Device Name: V.A.C. Instillamat (V.A.C.® Instill)

### Indications for Use:

The V.A.C.® Instill device is indicated for patients who would benefit from vacuum assisted drainage and controlled delivery of topical wound treatment solutions and suspensions over the wound bed.

The V.A.C.® Instill is intended for patients with chronic, acute, traumatic, subacute and dehiscent wounds, partial-thickness burns, ulcers (such as diabetic, pressure or venous insufficiency), flaps and grafts.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR Over-The-Counter Use             
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF  
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

David Krane for MXM  
(Division Sign-Off)

Division of Surgical, Orthopedic,  
and Restorative Devices

Page   1   of   1  

(Posted November 13, 2003)

510(k) Number K091585