

CryoPatch® SG Pulmonary Human Cardiac Patch 510(k) Submission  
CryoLife, Inc.

### 510(K) SUMMARY

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of 21 C.F.R. § 807.92. AUG 07 2009

**Submitter:** CryoLife, Inc.  
1655 Roberts Blvd., NW  
Kennesaw, GA 30144  
(770) 419-3355

Contact Person:	John D. Ferros Director, Regulatory Affairs		
Device Names:	Device Trade Name:	CryoPatch® SG Pulmonary Human Cardiac Patch	
	Common/Usual Name:	Allograft Patch	
	Classification Name:	Intracardiac patch or pledget	
Intended Use:	CryoPatch® SG Pulmonary Human Cardiac Patch is indicated for repair or reconstruction of the right ventricular outflow tract.		

**Predicate Devices:**

Device	Company	510 (k) Number(s), Clearance Date	Product Code
MatrACELL™ Pulmonary Artery Patch Allograft	LifeNet Health Virginia Beach, VA	K081438 – 10/17/2008	DXZ
CryoValve® SG Pulmonary Valve (and Conduit).	CryoLife, Inc. Kennesaw, GA	K033484 – 02/07/2008 K083106 – 02/06/2009	OHA

**Device Description:**

The CryoPatch® SG Pulmonary Human Cardiac Patch is derived from human pulmonary valve and artery tissue aseptically recovered from qualified donors. The patch is treated with an antimicrobial solution, and treated to remove the cells and cellular debris that have not already been removed during the post mortem period, harvesting, and the antimicrobial process. The patch is cryopreserved in a tissue culture medium, containing a cryoprotectant, within the innermost pouch of a three pouch packaging system. The packaging system not only withstands ultra cold temperatures, but also allows for aseptic introduction of the patch into the operating room. Supercooling by liquid nitrogen boost is begun prior to crystallization to minimize ice crystal damage to the patch matrix. Finally, the patch is transferred for long term storage at or below -135° C.

CryoPatch® SG is distributed in three anatomic configurations: pulmonary hemi-artery, pulmonary trunk, and pulmonary branch.

Implantation of the CryoPatch SG Pulmonary Human Cardiac Patch reduces the risk for induction of HLA class I and class II alloantibodies, based on Panel Reactive Antibody measured at up to one year, compared to standard-processed pulmonary cardiac tissues. Data have not been provided to evaluate the effect of reduced HLA class I and class II alloantibodies on the long-term durability, or long-term resistance to rejection by the patient, of the CryoPatch SG.

**Testing Supporting Substantial Equivalence:**

CryoPatch SG has undergone biomechanical, histological, DNA content, and decellularization evaluation. In addition, a large animal study and clinical data support the substantial equivalence of the device to previously marketed devices.

CONFIDENTIAL



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-0609  
Silver Spring, MD 20993-0002

AUG 07 2009

CryoLife, Inc.  
c/o Mr. John D. Ferros  
Director, Regulatory Affairs  
1655 Roberts Boulevard, NW  
Kennesaw, GA 30144

Re: K091626  
Trade/Device Name: CryoPatch® SG Pulmonary Human Cardiac Patch  
Regulation Number: 21 CFR 870.3470  
Regulation Name: Intracardiac Patch or Pledget  
Regulatory Class: Class: II (two)  
Product Code: DXZ  
Dated: June 1, 2009  
Received: June 3, 2009

Dear Mr. Ferros:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



 Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K091626

Device Name: CryoPatch® SG Pulmonary Human Cardiac Patch

Indications For Use: Indicated for repair or reconstruction of the right ventricular outflow tract.

Prescription Use ✓  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1

*Anna D. Kimes*  
(Division Sign-Off)  
Division of Cardiovascular Devices

510(k) Number K091626