

K091844

AUG 28 2009

4 510(k) Summary of Safety and Effectiveness

Manufacturer/Sponsor	Arthrex, Inc. 1370 Creekside Boulevard Naples, Florida 34108-1945
510(k) Contact	Ivette Galmez Regulatory Affairs Associate Telephone: 239/643.5553, ext. 1263 Fax: 239/598.5508 Email: igalmez@arthrex.com
Trade Name	Arthrex Bio-Composite SutureTak Anchors
Common Name	Suture Anchor
Product Code - Classification Name	MAI - Fastener, Fixation, Biodegradable, Soft Tissue HWC - Screw, Fixation, Bone
Predicate Device(s)	Arthrex Bio-Composite Suture Anchors: K071177
Device Description and Intended Use	<p>The Arthrex Bio-Composite SutureTak Anchors family is similar to the predicate devices in overall design. The difference lies in the introduction of smaller SutureTak models made with Bio-Composite material (PLLA or PLDLA combined with TCP).</p> <p>The Arthrex Bio-Composite SutureTak Anchors family is intended to be used for suture or tissue fixation in the foot, ankle, knee, hand, wrist, elbow, and shoulder. Please see indications for use form for specific indications.</p>
Substantial Equivalence Summary	The Arthrex Bio-Composite SutureTak Anchors family is substantially equivalent to the Arthrex Bio-Composite Suture Anchors predicate, in which the basic features and intended uses are the same. Any differences between the <i>Bio-Composite SutureTak Anchors</i> and the predicate K071177 are considered minor and do not raise questions concerning safety and effectiveness. Based on the information submitted, Arthrex, Inc. has determined that the new <i>Bio-Composite SutureTak Anchors</i> family is substantially equivalent to currently marketed predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration
10903 New Hampshire Avenue
Document Mail Center - WO66-G609
Silver Spring, MD 20993-0002

AUG 28 2009

Arthrex, Inc.
c/o Ms. Ivette Galmez
Regulatory Affairs Associate
1370 Creekside Blvd.
Naples, Florida 34108-1945

Re: K091844
Trade/Device Name: Arthrex Bio-Composite SutureTak Anchors
Regulation Number: 21 CFR 888.3030
Regulation Name: Smooth or threaded metallic bone fastener
Regulatory Class: Class II
Product Code: MAI, HWC
Dated: July 29, 2009
Received: July 30, 2009

Dear Ms. Galmez:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson
Director Division of Surgical, Orthopedic
and Restorative Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

3 Indications for Use Form

Indications for Use

510(k) Number:

K091844

Device Name:

Arthrex Bio-Composite Tak

The **Arthrex Bio-Composite Tak** is intended to be used for suture or tissue fixation in the foot, ankle, knee, hand, wrist, elbow, and shoulder. Specific indications are listed below and are size appropriate per patient needs:

- Elbow:** Biceps Tendon Reattachment, Ulnar or Radial Collateral Ligament Reconstruction.
- Shoulder:** Rotator Cuff Repairs, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair, Deltoid Repair, Capsular Shift or Capsulolabral Reconstruction.
- Hand/Wrist:** Scapholunate Ligament Reconstruction, Carpal Ligament Reconstruction, Repair/Reconstruction of Collateral Ligaments, Repair of Flexor and Extensor Tendons at the PIP, DIP and MCP joints for all digits, Digital Tendon Transfers.
- Foot/Ankle:** Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Metatarsal Ligament Repair, Hallux Valgus Reconstruction, Digital Tendon Transfers, Mid-foot reconstruction.
- Knee:** Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis.

Prescription Use ☒ AND/OR Over-The-Counter Use _____


(Per 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

PAGE 1 of 1


(Division Sign-Off) for MKM
Division of Surgical, Orthopedic,
and Restorative Devices

510(k) Number K091844