

BY DENTAL: Hi-Tech Implant/HT Physio Implant  
510(k): Original Submission Traditional

## SECTION 14: SUMMARY OF SAFETY AND EFFECTIVENESS

This 510(k) summary of safety and effectiveness information is being submitted in accordance with the requirements of SMD 1990 and CFR 807.92.

OCT 30 2009

### 14.1 SUBMITTER INFORMATION

- |    |                       |                                                                             |
|----|-----------------------|-----------------------------------------------------------------------------|
| a. | Company Name          | BY DENTAL SRL                                                               |
| b. | Company Address       | Via Vecchia Prov.le Lucchese<br>47/F<br>51030 Serravalle Pistoiese<br>Italy |
| c. | Company phone         | +39 0573 994 355                                                            |
|    | Company fax           | +39 0573 919 480                                                            |
| d. | Contact Person        | Daniele Poli<br>President                                                   |
| e. | Date Summary Prepared | June 22, 2009                                                               |

### 14.2 DEVICE IDENTIFICATION

- |    |                          |                                                      |
|----|--------------------------|------------------------------------------------------|
| a. | Trade/Proprietary Names: | HI TECH IMPLANT /HT<br>PHYSIO IMPLANT                |
| b. | Classification Name:     | Dental Handpieces and<br>Accessories 21 CFR 872.4200 |
| c. | Common Names:            | Surgical Micromotor for<br>implantology              |

### 14.3 IDENTIFICATION OF PREDICATE DEVICE

<u>Company</u>	<u>Device</u>	<u>510(k) No.</u>	<u>Date Cleared</u>
ATR	Implant System	K033597	25 Nov. 2003

### 14.4 DEVICE DESCRIPTION

The HI TECH IMPLANT and HT PHYSIO IMPLANT consist of a microprocessor controlled unit, foot pedal, electric micromotor, support rods and sterile, disposable irrigation tubes.

The unit also houses the peristaltic pump.

The HI TECH IMPLANT and HT PHYSIO IMPLANT provide electronic control of velocity and torque.

The HI TECH IMPLANT and HT PHYSIO IMPLANT can be programmed and retain programs into memory.

The HI TECH IMPLANT and HT PHYSIO IMPLANT are fully operational from foot pedal.

#### **14.5 SUBSTANTIAL EQUIVALENCE**

The HI TECH IMPLANT and HT PHYSIO IMPLANT are substantially equivalent to the ATR Implant System Surgical Micromotor in commercial distribution by ATR.

The fundamental technical characteristics of the HI TECH IMPLANT and HT PHYSIO IMPLANT are similar to those of the predicate device and are listed on the comparison chart provided in this 510(k) submission.

HI TECH IMPLANT and HT PHYSIO IMPLANT and the predicate device have adjustable speed, torque and reduction rates, and are programmable. The micromotor handpiece of HI TECH IMPLANT and HT PHYSIO IMPLANT and the one of the predicate device are autoclavable.

#### **14.6 INTENDED USE**

HI TECH IMPLANT and HT PHYSIO IMPLANT are intended for the preparation of intra-oral bone for implantology procedures.

#### **14.7 TECHNICAL CHARACTERISTICS**

HI TECH IMPLANT and HT PHYSIO IMPLANT were designed and developed to provide a microprocessor controlled surgical system with similar performances compared to predicate device.

- HI TECH IMPLANT and HT PHYSIO IMPLANT
- is equivalent in functions to the predicate device,
- has adjustable speeds and torque values (related to the reduction rate of the handpiece selected) that are fully customisable by the end user,
- has automatic motor shutdown system which provides to switch off the motor whenever set torque is reached,
- has a sterilizable micromotor according the recommended protocols,
- can be operated by foot control.

Any of these feature is found in the predicate device also.

#### **14.8 PERFORMANCE DATA**

HI TECH IMPLANT and HT PHYSIO IMPLANT were tested in accordance with the technical requirements of IEC 60601-1 and IEC 60601-1-2.

All evaluation of the HI TECH IMPLANT and HT PHYSIO IMPLANT were performed by Notified Laboratory, and all the results comply to standard listed below.

The conclusions drawn from the performance test are that HI TECH IMPLANT and HT PHYSIO IMPLANT comply with:

- IEC 55011 EMC Conducted RF emissions
- IEC 60529 Degrees of protection provided by enclosures (IP Code)
- IEC 60601-1 Medical Electrical Equipment Part 2,"General Safety Norms"
- IEC 60601-1-2 Medical Electrical Equipment Part 1 EMC
- IEC 61000-4-2 EMC Electrostatic Discharge Immunity
- IEC 61000-4-3 Radiated RF immunity
- IEC 61000-4-4 Fast Transient Immunity
- IEC 61000-4-5 Pulse immunity
- IEC 61000-4-6 Conducted RF immunity
- IEC 61000-4-11 Supply Voltage Hole Immunity
- IEC 61558-1 Transformer safety
- IEC 61558-2-6 Transformer safety-Particular prescriptions
- ISO 11498:1997 Dental handpieces – Dental low-voltage electrical motors

and are effective and safe to use.

Declaration of conformity to a standard in Chapter 12.1 of this submission.

#### **14.9 510(k) CHECKLIST**

This notification contains all information required by 21 CFR 807.87.

A completed copy of the Premarket Notification 510(k) Reviewer's Checklist is provided in this submission.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-O66-0609  
Silver Spring, MD 20993-0002

Dr. Daniele Poli  
President  
By Dental S.R.L.  
Via Vecchia Prov. le Lucchese 49/FG  
Serravalle Pistoiese, Pistoia  
ITALY 51030

**OCT 30 2009**

Re: K091913  
Trade/Device Name: Hi-Tech Implant/HT Physio Implant  
Regulation Number: 21 CFR 872.4200  
Regulation Name: Dental Handpiece and Accessories  
8Regulatory Class: I  
Product Code: EKX  
Dated: September 21, 2009  
Received: September 29, 2009

Dear Dr. Poli:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please contact the CDRH/Office of Surveillance and Biometrics/Division of Postmarket Surveillance at 240-276-3464. For more information regarding the reporting of adverse events, please go to <http://www.fda.gov/cdrh/mdr/>.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "for" followed by a stylized signature.

Susan Runner, D.D.S., M.A.

Acting Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

K091913

BY DENTAL: Hi-Tech Implant/HT Physio Implant  
510(k): Original Submission Traditional

INDICATIONS FOR USE

510(k) Number:

To be Assigned by FDA

Device Name:

Hi-Tech Implant/HT Physio Implant

Indications for Use:

Hi-Tech Implant/HT Physio Implant are intended to prepare Intraoral bone for implantology procedures.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒

OR Over-The-Counter Use ☐

(Per 21 CFR 801.109)

Ken Muly for MSR  
(Division Sign-Off)

Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K091913