

## 5.0 510(k) Summary

### 1. Sponsor

SpineFrontier, Inc.  
500 Cummings Center  
Suite 3500  
Beverly, MA 01915

NOV 19 2009

**Primary Contact:** John Sullivan  
**Telephone:** 1- 978-232-3990

**Date Prepared:** November 12, 2009

### 2. Device Name and Classification:

**Proprietary Names:** **KRD1™ Pedicle Screw System**

**Common/Usual Name:** Spine Fixation System

**Classification Name:** Orthosis, Spondylolisthesis Spinal Fixation  
Orthosis, Spinal Pedicle Fixation  
(21 CFR 888.3070(b)(1)), Class II

**Product Code:** MNH, MNI

### 3. Predicate Devices

K950099 – Synergy Posterior Spinal System  
K992739 – Synthes USS Click-X Variable Axis System  
K052123 – Synthes Pangea System

### 4. Device Description

The **KRD1™ Pedicle Screw System** consists of longitudinal rods, polyaxial screws, transverse connectors, and instrumentation to facilitate installation of this system.

The **KRD1™ Pedicle Screw** is offered in multiple screw diameters, with multiple lengths to accommodate individual patient needs. **KRD1™ Pedicle Screw** assemblies consist of the screw, a tulip, and a washer that secures the tulip to the screw and also serves as the saddle for the longitudinal rod. The

tulip and screw designs allow the screw to have a polyaxial rotation relative to the tulip. A locking cap is placed on top of the pedicle screw assembly to secure the position of the implant and to retain the longitudinal rod.

The **KRD1™ Pedicle Screw** assembly serves as the central fixation device to which various rods and cross connectors are secured to provide desired fixation.

Longitudinal rods are provided in two configurations, straight and lordotic, have a fixed diameter (5.5mm), and vary by length.

Cross connector assemblies are provided in multiple configurations, varying by length. Cross connectors are used to provide additional fixation support.

The **KRD1™ Pedicle Screw System** components are fabricated from medical grade titanium alloy conforming to ASTM F-136 specifications.

## 5. Intended Use

The **KRD1™ Pedicle Screw System** is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute instabilities or deformities of the thoracic, lumbar, and sacral spine: degenerative spondylolisthesis with objective evidence of neurological impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor, and failed previous fusion (Pseudoarthrosis).

In addition, the **KRD1™ Pedicle Screw System** is indicated for the treatment of severe spondylolisthesis (Grade 3 and 4) of the L5-S1 vertebra in skeletally mature patients receiving fusion by autogenous bone graft.

## 6. Technological Characteristics & Substantial Equivalence Determination

The SpineFrontier **KRD1™ Pedicle Screw System** was shown to be substantially equivalent to predicate devices through comparison of indications for use, function, operating principles, and materials.

The **KRD1™ Pedicle Screw System** has been shown to be substantially equivalent to predicate devices in terms of performance (mechanical testing). Clinical data was not required for this device.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-O66-0609  
Silver Spring, MD 20993-0002

**NOV 19 2009**

SpineFrontier, Inc.  
% Mr. John Sullivan  
Director, QA and Regulatory Compliance  
500 Cummings Center, Suite 3500  
Beverly, Massachusetts 01915

Re: K092420  
Trade/Device Name: KRD1™ Pedicle Screw System  
Regulation Number: 21 CFR 888.3070  
Regulation Name: Pedicle screw spinal system  
Regulatory Class: Class II  
Product Code: MNH, MNI  
Dated: October 20, 2009  
Received: October 20, 2009

Dear Mr. Sullivan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

Page 2 – Mr. John Sullivan


device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucml115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/cdrh/mdr/> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



 Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

#### 4.0 Indications for Use Statement

510(k) Number (if Known): K092420

##### Indications For Use:


The **KRD1™ Pedicle Screw System** is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute instabilities or deformities of the thoracic, lumbar, and sacral spine: degenerative spondylolisthesis with objective evidence of neurological impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor, and failed previous fusion (Pseudoarthrosis).

In addition, the **KRD1™ Pedicle Screw System** is indicated for the treatment of severe spondylolisthesis (Grade 3 and 4) of the L5-S1 vertebra in skeletally mature patients receiving fusion by autogenous bone.

Prescription Use: X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use: \_\_\_\_\_  
(Part 21 CFR 807 Subpart C)

  
(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number K092420