

K092507  
162

5. 510(K) SUMMARY AS REQUIRED BY 21 CFR 807.92

NOV 12 2009

Submitter:

Ash Access Technology, Inc.  
3601 Sagamore Pkwy N, Suite B  
Lafayette, IN 47904  
Tel. 765.742.4813  
Fax 765.742.4823

Contact:

Roland Winger  
Ash Access Technology, Inc.  
3601 Sagamore Pkwy N, Suite B  
Lafayette, IN 47904  
Tel. 765.742.4813

Date Summary was Prepared: August 20, 2009

Device Proprietary Name: Centros® Chronic Hemodialysis Catheter Set  
Classification Name: Catheter, Hemodialysis, Implanted  
Device Product Code: MSD

Predicate Device(s): The Centros® Chronic Hemodialysis Catheter Set is substantially equivalent to the following devices:

Device Name	Manufacturer	510(k) No.
Centros® Chronic Hemodialysis Catheter Set	Ash Access Technology, Inc.	K070572, cleared 6/26/2007
AngioDynamics Duramax™ Catheter Kits	AngioDynamics®	K080400, cleared 5/13/08
AngioDynamics® EvenMore® Chronic Hemodialysis Catheter	AngioDynamics®	K040402, cleared 4/19/2004

**Device Description:**

The Centros® Chronic Hemodialysis Catheter Set is a dual lumen, 15FR catheter available in multiple lengths, in straight and pre-curved configurations. The catheter lumens are D-shaped and made from radiopaque Carbothane. The distal end design is a fixed length split-tip, without side-holes. The distal venous lumen extends past the arterial lumen, and includes a guidewire slit for insertion by the optional guidewire placement technique. The proximal device contains a fixed polyester cuff, an integrated hub, suture wing, and extension set with color coded occlusion clamps and luer connectors (red and blue for the arterial and venous lumens respectively). The lumen priming volumes are printed on the clamps. The procedure kit includes the necessary accessories to correctly insert the catheter.

000016

**Statement of Intended Use:**

The Centros® Chronic Hemodialysis Catheter Set is indicated for use in attaining long-term vascular access for hemodialysis and apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. This catheter is indicated for > 30 days long-term placement.

**Discussion of Nonclinical Tests:**

The safety and performance of the Centros® Chronic Hemodialysis Catheter Set have been substantiated through extensive non-clinical testing, including tensile strength, joint strength, leakage, flow rate, kit component compatibility.

Results of testing show that the Centros® Chronic Hemodialysis Catheter Set can reliably perform as a conventional hemodialysis catheter for obtaining blood access for hemodialysis and apheresis. No new questions of safety or effectiveness have been raised.

**Substantial Equivalence:**

The Centros® Chronic Hemodialysis Catheter Set product information, technological comparison to predicate products, and test results demonstrate that the Centros® Chronic Hemodialysis Catheter Set is safe and performs as intended. The Centros® Chronic Hemodialysis Catheter Set is substantially equivalent to the currently marketed predicate devices with respect to intended use, materials, technological characteristics, performance, insertion method, anatomical location, kit components, and labeling.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-066-0609  
Silver Spring, MD 20993-0002

Mr. Roland K. Winger  
VP, Clinical and Product Development  
Ash Access Technology, Inc.  
3601 Sagamore Pkwy. N., Suite B  
LAFAYETTE IN 47904

NOV 12 2009

Re: K092597

Trade/Device Name: Centros™ Chronic Hemodialysis Catheter Set  
Regulation Number: 21 CFR §876.5540  
Regulation Name: Blood access device accessories  
Regulatory Class: III  
Product Code: MSD  
Dated: August 21, 2009  
Received: August 24, 2009

Dear Mr. Winger:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. However, you are responsible to determine that the medical devices you use as components in the kit have either been determined as substantially equivalent under the premarket notification process (Section 510(k) of the act), or were legally on the market prior to May 28, 1976, the enactment date of the Medical Device Amendments. *Please note:* If you purchase your device components in bulk (i.e., unfinished) and further process (e.g., sterilize) you must submit a new 510(k) before including these components in your kit. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, and labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

  
for Janine M. Morris, Director (Acting)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

K092597  
12/1

4. INDICATIONS FOR USE STATEMENT

**Indications for Use**

510(k) Number (if known): K092597

Device Name: Centros™ Chronic Hemodialysis Catheter Set

Indications for Use: The Centros™ Chronic Hemodialysis Catheter Set is indicated for use in attaining long-term vascular access for hemodialysis and apheresis.

It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. This catheter is indicated for > 30 days long-term placement.


Prescription Use ✓  
(Part 21 CFR 801 Subpart  
D)

AND/OR

Over-The-Counter Use  
(21 CFR 801 Subpart  
C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER  
PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K092597

000015