

1092630

**510(k) SUMMARY**

**OCT 22 2009**

**Submitter:**

Name: ZENCERA, INC.  
Address: Room 1108, ACE Techno Tower 5<sup>th</sup> Building, 197-22, Guro3-Dong, Guro-Gu, Seoul, Republic of Korea, 152-766

**US Agent:**

Name: Haeryang Joo  
Address: 4525 Wilshire Blvd., Los Angeles, CA, 90010, U.S.A.  
Phone No.: +1-323-935-0606  
Fax No.: +1-323-935-7779  
E-Mail: sylee68@gmail.com

**Official Correspondent:**

Name: San-ho Hyun  
Address: Room 1108, ACE Techno Tower 5<sup>th</sup> Building, 197-22, Guro3-Dong, Guro-Gu, Seoul, Republic of Korea, 152-766  
Phone No.: +82-2-2109-0955  
Fax No.: +82-2-2109-0957  
E-Mail: sanho0517@hanmail.net

**Device Identification**

Proprietary Name: ZIROX  
Common/Usual Name: Porcelain Powder  
Classification Name: Porcelain Powder for Clinical Use  
Product Code: EIH  
Review Panel: Dental  
Regulation Number: 872.6660

**Substantially Equivalent Predicate Legally Marketed Devices**

The subject device is deemed to be substantially equivalent to those following devices manufactured and currently available in commercial distribution.

Device Name	Dentsply-Cercon Base	3M-LAVA Zirconia	Sagemax Bioceramics- Sagemax Z-Blank
510(k) Number	K013230	K011394	K062695
Decision Date	10/25/2001	06/29/2001	10/20/2006
Decision	Substantially Equivalent	Substantially Equivalent	Substantially Equivalent
Product Code	EIH	EIH	EIH
Regulation Number	872.6660	872.6660	872.6660

#### **Device Description**

ZIROX is a pre-formed machineable dental blank composed of zirconium oxide. ZIROX is available in partially-sintered. ZIROX is available in different shapes, and dimensions. ZIROX has two types that are 51x23x16 and 51x23x18. The only difference between 51x23x16 and 51x23x18 is the size(height).

ZIROX is a pre-formed ceramic dental blank intended for CAD/CAM fabrication of zirconia frameworks for all-ceramic dental restorations. ZIROX is designed for manufacturing ceramic dental restorations such as single crowns or bridgeworks. The blank is machined by the customers/dental laboratories on their milling centers or similar equipment using CAD/CAM techniques for design.

#### **Indications for Use**

The ZIROX is intended for CAD/CAM fabrication of all-ceramic dental restorations. The ZIROX is used for the manufacturing of inlays, onlays, veneers, crowns and bridges.

#### **Technological Characteristics and Substantial Equivalence**

ZIROX and predicate devices are identical in intended use and material. Therewith, ZIROX and predicate devices are biocompatible and have similar biomechanical strength and properties.

Based on the discussion above, ZENCERA, INC. believes that ZIROX is substantially equivalent in comparison with predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-066-0609  
Silver Spring, MD 20993-0002

Mr. Sanho Hyun  
Assistant Manager  
Zencera, Incorporated  
Ace Techno Tower 5<sup>th</sup>, Room 1108  
197-22 Guro-3Dong, Guro-Gu  
Seoul 152766  
REPUBLIC OF KOREA

OCT 22 2009

Re: K092630  
Trade/Device Name: ZIROX  
Regulation Number: 21 CFR 872.6660  
Regulation Name: Porcelain Powder for Clinical Use  
Regulatory Class: II  
Product Code: EIH  
Dated: August 13, 2009  
Received: August 27, 2009

Dear Mr. Hyun:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please contact the CDRH/Office of Surveillance and Biometrics/Division of Postmarket Surveillance at 240-276-3464. For more information regarding the reporting of adverse events, please go to <http://www.fda.gov/cdrh/mdr/>.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Susan Runner", followed by the word "for" in a smaller, cursive script.

Susan Runner, D.D.S., M.A.

Acting Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known): K092630

Device Name: ZIROX

Indications for Use:

The ZIROX is intended for CAD/CAM fabrication of all-ceramic dental restorations. The ZIROX is used for the manufacturing of inlays, onlays, veneers, crowns and bridges.

Prescription Use ✓  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

*Rei Muly for KSR*  
(Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K092630