

K092741

510(k) Summary
as required by 807.92

NOV 10 2009

1. Company Identification

EIZO NANA CORPORATION
153 Shimokashiwano-cho, Hakusan, Ishikawa-ken, 924-8566, Japan
Tel: +81-76-274-2468
Fax: +81-76-274-2484

2. Official Correspondent

Hiroaki Hashimoto (Mr.)
Manager of Engineering Management Section

3. Date of Submission

September 4, 2009

4. Device Trade name

Color LCD Monitor, RadiForce MX241W

5. Common/Usual Name:

Image display system, medical image workstation, image monitor/display, and others

6. Classification Number:

Medical displays classified in Class II per 21 CFR 892.2050.

7. Predicate Device

Manufacturer : EIZO NANA CORPORATION
Device Name : Color LCD Monitor
Model Name : FlexScan MX300W
510(k) No. : K073340

8. Description of Device

RadiForce MX241W is 61cm Color LCD display for medical image viewing digital images viewing.

9. Intended Use

RadiForce MX241W is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners. The device is not specified for digital mammography system.

10. Substantial Equivalence to Predicate Device

RadiForce MX241W is substantially equivalent to FlexScan MX300W (K073340). The panel size is changed to 61 cm (24.1") from 76cm (29.8"). MX241W improved the brightness. The brightness improved to 320 cd/m2 from 300 cd/m2. DisplayPort is added as new input terminal.

Comparison table of the principal characteristics of 2 devices is shown in the Attachment 1.
Appendix 1: Comparison Table with Predicate Device



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

NOV 10 2009

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

Mr. Hiroaki Hashimoto
Manager
EIZO NANA O Corporation, Engineering Management Section
153 Shimokashiwano-cho,
Hakusan, Ishikawa-ken 924-8566
JAPAN

Re: K092741

Trade/Device Name: Color LCD Monitor, RadiForce MX241W
Regulation Number: 21 CFR 892.2050
Regulation Name: Picture archiving and communications system
Regulatory Class: II
Product Code: LLZ
Dated: October 21, 2009
Received: October 22, 2009

Dear Mr. Hashimoto:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical


Page 2 --

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Janine M. Morris", is written over the typed name.

Janine M. Morris
Acting Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

K092741

510(k) Number:

K092741

Device Name : Color LCD Monitor, RadiForce MX241W

Indications for Use:

RadiForce MX241W is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners. RadiForce MX241W does not support the display of mammography images for diagnosis.


Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

Division of Reproductive, Abdominal,
and Radiological Devices

510(k) Number

K092741