

**510(k) SUMMARY OF SAFETY AND EFFECTIVENESS**

As required by section 807.92(c)

NOV - 4 2009

Submitter	SPINEART International Center Cointrin 20 route de pré-bois CP1813 1215 GENEVA 15 SWITZERLAND
Contacts	Franck PENNESI Director of Industry & Quality Phone : +41 22 799 40 25 Fax : +41 22 799 40 26 Mail : <a href="mailto:fpennesi@spineart.ch">fpennesi@spineart.ch</a> Regulatory contact : Dr Isabelle DRUBAIX (Idée Consulting) <a href="mailto:idrubaix@nordnet.fr">idrubaix@nordnet.fr</a>
Trade Name	ROMEO posterior osteosynthesis system
SPECIAL 510k	Modification to ELLIPSE posterior osteosynthesis system
Classification Name	Pedicle screw spinal system
Class	II
Product Code	MNI orthosis, spinal pedicle fixation
Subsequent product codes	MNH orthosis, spondylolisthesis spinal fixation KWP Spinal interlaminar fixation orthosis
CFR section	888.3070
Device panel	Orthopedic
Legally marketed predicate devices	ELLIPSE posterior osteosynthesis system (K081165) manufactured by SPINEART
Description	ROMEO posterior osteosynthesis system includes pedicular screws, spondylolisthesis screws longitudinal rods and transverse connector rods, connector and nut. All components of ROMEO posterior osteosynthesis system are made of TA6V4ELI conforming to ISO 5832.3 and ASTM F 136. ROMEO components are supplied either sterile or not sterile and with a complete set of surgical instruments.

Intended Use	ROMEO posterior osteosynthesis system is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of the thoracic, lumbar, and sacral spine: severe spondylolisthesis (grades 3 and 4) of the L5-S1 vertebra; degenerative spondylolisthesis with objective evidence of neurologic impairment; fracture; dislocation; scoliosis; kyphosis; spinal tumor; and failed previous fusion (pseudarthrosis).
Performance data	ROMEO posterior osteosynthesis system conforms to special control established for Pedicle screw spinal system and to « Spinal System 510(k)s - Guidance for Industry and FDA Staff Document » issued on: May 3, 2004.
Substantial equivalence	ROMEO posterior osteosynthesis system is substantially equivalent to its predicate device in terms of intended use, material, design, mechanical properties and function.

Preparation date, September 30, 2009



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-066-0609  
Silver Spring, MD 20993-0002

Spineart  
% Mr. Franck Pennesi  
Director of Industry and Quality  
International Center Cointrin  
20, route de Pre-Bois  
CP 1813  
1215 Geneva 15  
Switzerland

NOV - 4 2009

Re: K093170

Trade/Device Name: ROMEO posterior osteosynthesis system  
Regulation Number: 21 CFR 888.3070  
Regulation Name: Pedicle screw spinal system  
Regulatory Class: Class II  
Product Code: MNI, MNH, KWP  
Dated: October 5, 2009  
Received: October 7, 2009

Dear Mr. Pennesi:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act

or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/cdrh/mdr/> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Mark N. Melkerson

Director

Division of Surgical, Orthopedic  
and Restorative Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

## INDICATIONS FOR USE

**510(k) Number (if known):** K093170

**Device Name:** ROMEO posterior osteosynthesis system

**Indications for Use:**

ROMEO posterior osteosynthesis system is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of the thoracic, lumbar, and sacral spine: severe spondylolisthesis (grades 3 and 4) of the L5-S1 vertebra; degenerative spondylolisthesis with objective evidence of neurologic impairment; fracture; dislocation; scoliosis; kyphosis; spinal tumor; and failed previous fusion (pseudarthrosis).

Prescription Use <input checked="" type="checkbox"/>	Over-The-Counter Use
(Part 21 CFR 801 Subpart D)	(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number K093170