

510(k) Summary
LDR Spine SpineTune™ TL Spinal System

MAR 31 2010

1. Owner's Name & Address

LDR Spine USA

4030 West Braker Lane, Suite 360
Austin, TX 78759

Phone: (512) 344-3333

Fax: (512) 344-3350

2. Contact Person

Noah Bartsch, MS, RAC

Manager, Clinical, Regulatory & Quality Affairs

LDR Spine USA

4030 West Braker Lane, Suite 360

Austin, TX 78759

Phone: (512) 344-3319

Fax: (512) 344-3350

Email: noahbartsch@ldrspine.com

3. Date 510(k) Summary Prepared: February 26, 2010

- 4. Trade Name:** LDR Spine SpineTune™ TL Spinal System
Common Name: Spinal Fixation System (MNH, MNI, KWP)
Classification: MNH 888.3070 – Orthosis, Spondilolisthesis Spinal Fixation
MNI 888.3070 – Orthosis, Spinal Pedicle Fixation
KWP 888.3050 – Orthosis, Spinal Interlaminar Fixation

5. Legally Marketed Equivalent Predicate Device:

Perfix Spinal System (K091725)

6. Device Description

The Spinetune TL Spinal System is a top-loading posterior spinal fixation system consisting of pedicle screws, rods, set screws, transverse connectors and lateral, axial and domino type connectors.

7. Intended Use of the device

The SpineTune™ TL Spinal System is a posterior, noncervical pedicle fixation system indicated to provide immobilization and stabilization of spinal segments in skeletally-mature patients as an adjunct to fusion by autogenous bone graft in the

treatment of the following acute and chronic instabilities or deformities of thoracic, lumbar and sacral spine:

- Spondylolisthesis (Grade 3 and 4)
- Degenerative spondylolisthesis with objective evidence of neurological impairment
- Trauma (i.e., fracture or dislocation)
- Spinal stenosis
- Deformities or curvatures (i.e., scoliosis, kyphosis, and/or lordosis)
- Tumor
- Pseudoarthrosis
- Failed previous fusion

8. Non-Clinical Performance Data

Non Clinical testing was not required.

The outcomes of verification and validation activities indicate that the propose SpineTune™ TL Spinal System is substantially equivalent to the predicate device.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room W-066-0609
Silver Spring, MD 20993-0002

MAR 31 2010

LDR Spine USA
% Mr. Noah Bartsch, MS, RAC
Manager, Clinical, Regulatory and
Quality Affairs
4030 West Braker Lane, Suite 360
Austin, Texas 78759

Re: K100575

Trade/Device Name: SpineTune™ TL Spinal System
Regulation Number: 21 CFR 888.3070
Regulation Name: Pedicle screw spinal system
Regulatory Class: Class II
Product Code: MNI, MNH, KWP
Dated: February 26, 2010
Received: March 01, 2010

Dear Mr. Bartsch:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must

Page 2 - Mr. Noah Bartsch, MS, RAC

comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson".

Mark N. Melkerson

Director

Division of Surgical, Orthopedic

And Restorative Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

INDICATIONS FOR USE

510(k) Number (if known):

Device Name: LDR Spine SpineTune™ TL Spinal System

Indications for Use:

The SpineTune™ TL Spinal System is a posterior, noncervical pedicle fixation system indicated to provide immobilization and stabilization of spinal segments in skeletally-mature patients as an adjunct to fusion by autogenous bone graft in the treatment of the following acute and chronic instabilities or deformities of thoracic, lumbar and sacral spine:

- Spondylolisthesis (Grade 3 and 4)
- Degenerative spondylolisthesis with objective evidence of neurological impairment
- Trauma (i.e., fracture or dislocation)
- Spinal stenosis
- Deformities or curvatures (i.e., scoliosis, kyphosis, and/or lordosis)
- Tumor
- Pseudoarthrosis
- Failed previous fusion


Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of Surgical, Orthopedic,
and Restorative Devices

510(k) Number: K100575