



K 101756

---

**SYBRON DENTAL SPECIALTIES**

**Section III - 510(k) Summary of Safety and Effectiveness**

**Submitter:**

Sybron Dental Specialties, Inc.  
1717 W. Collins Avenue  
Orange, California 92867  
(714) 516-7602 - Phone  
(714) 516-7472 - Facsimile  
Wendy Garman - Contact Person

OCT 26 2010

Date Summary Prepared: October 2010

**Device Name:**

- Trade Name – TSC
- Common Name – Gingival Retraction/Hemostatic Paste
- Classification Name – Unclassified

**Device for Which Substantial Equivalence is Claimed:**

- *Expa-syl*, Produits Dentaires Pierre Rolland

**Device Description:**

TSC is a paste containing a hemostatic agent which is intended to be used for the temporary retraction and hemostasis of the gingival margin during dental procedures such as, but not limited to, dental impressions, seating of temporary and permanent restorations, restorations of cavities and placement of a rubber dam. TSC is a tissue management solution that allows clinicians to quickly and easily obtain sulcular expansion in clinical situations prior to an impression. Additionally, TSC will help stop bleeding and prevent the flow of crevicular fluid upon removal, further assuring accurate and complete impressions.

**Intended Use of the Device:**

TSC is a paste containing a hemostatic agent which is intended to be used for the temporary retraction and hemostasis of the gingival margin during dental procedures

such as, but not limited to, dental impressions, seating of temporary and permanent restorations, restorations of cavities and placement of a rubber dam.

---

**Substantial Equivalence:**

TSC is substantially equivalent to another legally marketed device in the United States. TSC functions in a manner similar to *Expa-syl*, which is currently marketed by Kerr Corporation. TSC is different in formulation from the predicate device in that it contains a glass filler which has been used in numerous other Kerr products. The application and mode of operation is the same for both products.

**Nonclinical Test Data:**

TSC was tested according to the following biocompatibility studies per ISO 10993-5 and 10993-10: in vitro cytotoxicity, sensitization and oral irritation. The results of the studies are listed in the table below.

Biocompatibility Study	Result
In Vitro Cytotoxicity	Non-cytotoxic
Sensitization	Elicited no skin reaction / weak allergenic potential
Oral Irritation	Non-irritant

This 510(k) submission also includes data from bench testing used to evaluate the performance characteristics of TSC compared to the predicate device, *Expa-syl*. The characteristics evaluated include viscosity and rinse time.

**Clinical Test Data:**

Clinical testing has not been conducted on this product.

**Conclusion:**

Based upon the biocompatibility test and bench testing, the clinical performance of TSC is substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

Kerr Corporation  
C/O Ms. Wendy Garman  
Director, Regulatory Affairs  
Sybron Dental Specialties  
1717 West Collins Avenue  
Orange, California 92867

OCT 26 2010

Re: K101756  
Trade/Device Name: TSC  
Regulation Number: Unclassified  
Regulation Name: None  
Regulatory Class: Unclassified  
Product Code: MVL  
Dated: October 13, 2010  
Received: October 14, 2010

Dear Ms. Garman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to

<http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Anthony Watson".

Anthony Watson, B.S., M.S., M.B.A.  
Director

Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K101756

OCT 26 2010

Device Name: TSC

### Indications For Use:

TSC is a paste containing a hemostatic agent which is intended to be used for the temporary retraction and hemostasis of the gingival margin during dental procedures such as, but not limited to, dental impressions, seating of temporary and permanent restorations, restorations of cavities and placement of a rubber dam.

Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Susan Bussell  
(Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

Page 1 of 1

510(k) Number: K101756