

EXHIBIT #3

510(k) Summary

DEC 20 2010

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR §807.92.

Date: March 8, 2010

Updated: November 11, 2010

- 1. Applicant:**
Sunmax Vietnam Co. Ltd
Km 8, Pham Van Dong Road,
Hai Thanhdong Kinh District Hai Phong,
Hai Phong, Tahn Pho, 18671, Vietnam
- 2. Manufacturer:**
Sunmax Vietnam Co. Ltd
Km 8, Pham Van Dong Road,
Hai Thanhdong Kinh District Hai Phong,
Hai Phong, Tahn Pho, 18671, Vietnam
- 3. Submitter:**
Mr. Jigar Shah
Official Correspondent for
Sunmax Vietnam Co. Ltd
- 4. Address:**
mdi Consultants, Inc.
55 Northern Blvd., Suite 200
Great Neck, New York 11021
Tel: 516-482-9001
Fax: 516-482-0186
Jigar@mdiconsultants.com
- 5. Trade/proprietary Name:**
Sunmax Vietnam Co. Ltd Powder free Nitrile Patient Examination Gloves.
- 6. Common Names:**
POWDER-FREE Patient Examination Glove
- 7. Classification name:**
Patient Examination Glove
- 8. Classification number:**
21 CFR 880.6250

9. Device Description:

Sunmax Vietnam Co. Ltd Powder free Nitrile Examination Glove is a class I device having product code 80LZA. It is a disposable device that meets all requirements of ASTM D 631900a-05.

10. Intended Use:

A powder free patient examination glove is a disposable device intended for medical purposes that is worn on the examiners hand or finger to prevent contamination between patient and examiner.

11. Substantial Equivalence Discussion:

A powder free patient examination glove is substantially equivalent to the predicate devices.

| Characteristic and parameters | Sunmax Vietnam Co., LTD Powder Free Nitrile Examination Glove (New Device) | Sunmax Enterprise Shanghai Co., LTD Powder free Blue Nitrile Patient Examination Glove tested with chemotherapy drugs. (K090336) | ULTRAWIN SDN BHD Non-Sterile Powder Free Nitrile Examination Gloves (K 090828) | PT. MAHAKA RYA INTI BUANA Powder Free Black Nitrile Examination gloves. (K090464) | Substantially Equivalence Comparison |
|--------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|
| Product Code | LZA | LZA/LZC | LZA | LZA | |
| Intended Use | A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner. | A patient examination gloves (Tested for Use with Chemotherapy Drugs) is a disposable device intended for medical purpose that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner. | A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner. | A powder-free patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner. | SE |
| Width (size medium) | 89 | 89 | 93-98 | 97.7 | SE |
| Overall length | 240 | 240 | 240 | 240.9 | SE |
| Palm thickness | 0.12 | 0.15 | Min 0.08 | 0.109 | Minor Difference |

| | | | | | |
|-----------------------------------|------------------------------------------------------------------------------------------------------------------------|-----------------------------------|-----------|---------|-------------------|
| Finger thickness | 0.12 | 0.16 | Min 0.08 | 0.148 | Minor Difference |
| Tensile strength pre aging min | 22 | 16.4 | 15 - 21 | 18.8 | Minor Difference |
| Tensile strength after aging min | 23.6 | 18.2 | 14-22 | 21.3 | Minor Difference |
| Ultimate elongation pre aging min | 500 | 510 | 550 - 630 | 679.4 | Minor Difference |
| Ultimate elongation after aging | 500 | 520 | 520 - 610 | 767.4 | Minor Difference |
| Meets Biocompatibility | yes | yes | Yes | Yes | SE |
| Duration of bio-compatibility | Limited | Limited | Limited | Limited | SE |
| Skin irritation test | Passes | Passes | Passes | Passes | SE |
| Dermal sensitization | Passes | Passes | Passes | Passes | SE |
| Residual powder test | Passes | Passes | Passes | Passes | SE |
| Labeling | Guidance document "Medical Glove Guidance Manual" has been thoroughly followed with respect to Labeling of the device. | Specialty Medical Gloves Labeling | NA | NA | Minor Differences |

12. Summary of Testing:

| Test | | Results |
|-------------|---------------------------------------|------------------------------------------------------|
| a. | Dermal Sensitization Test | Passes |
| b. | Primary Skin irritation | Passes |
| c. | Permeation testing per ASTM D 6978-05 | Passes |
| d. | Iodine Test | Passes |
| e. | Tensile strength | Gloves meets the requirements of ASTM D63 19-00a. |
| f. | Barrier strength | Gloves meets the requirements of ASTM D63 19-00a. |

The standards used by Sunmax Vietnam Co. Ltd to determine substantial equivalence are based on ASTM D 631900a-2005. All testing meets requirements for physical specifications and dimensions conducted on gloves, Inspection level S-2, AQL 4.0, pinholes at AQL 2.5

There are no special labeling claims and we do not claim our gloves to be hypoallergenic.

13. Conclusion:

Powder free Nitrile Patient Examination Glove performance was equivalent to any other conventional method evaluated. Our evaluation concluded that our device raises no new issues of Safety and effectiveness.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room -WO66-G609
Silver Spring, MD 20993-0002

Sunmax Vietnam Company, Limited
C/O Mr. Jigar Shah
Official Correspondent
MDI Consulting, Incorporated
55 Northern Boulevard, Suite 200
Great Neck, New York 11021

DEC 20 2010

Re: K101870
Trade/Device Name: Powder Free Nitrile Patient Examination Gloves
Regulation Number: 21 CFR 880.6250
Regulation Name: Patient Examination Glove
Regulatory Class: I
Product Code: LZA
Dated: November 17, 2010
Received: November 19, 2010

Dear Mr. Shah:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Anthony D. Watson, B.S., M.S., M.B.A.
Director

Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number: K 101870

Applicant: Sunmax Vietnam Co. Ltd

DEC 20 2010

Device Name: Powder Free Nitrile Patient Examination Gloves.


Indications for Use:

A powder free patient examination glove is a disposable device intended for medical purposes that is worn on the examiners hand or finger to prevent contamination between patient and examiner.

Prescription Use _____ AND/OR Over-The-Counter Use X
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

Division of Anesthesiology, General Hospital
Infection Control, Dental Devices

510(k) Number: K101870