

510(k) Summary

AUG 26 2011

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirement of SMDA and 21 CFR 807.92.

1.0 submitter's information

Name: Andon Medical Co., Ltd.
Address: No.04-23-3 AIRPORT INDUSTRIAL PARK, TIANJIN
Phone number: 86-22-8761 2426
Fax number: 86-22-6052 6162
Contact: Yi Liu
Date of Application: 9/10/2010

2.0 Device information

Trade name: Andon Health Care Management System Software
Common name: Data management software
Classification name: Data management software

3.0 Classification

Production code: NBW- Blood Glucose Monitoring System.
Regulation number: 862.1345
Classification: II
Panel: Clinical Chemistry

4.0 Predict device information

Manufacturer: Taidoc Technology Corporation
Device: Clever Chek Health Care System Software
510(k) number: k070941

5.0 Device description

Andon Health Care Management System Software is an optional software accessory for use with AG-608 Single Blood Glucose Monitoring System and AG-608 Multi Blood Glucose Monitoring System. When used with AG-608 Single Blood Glucose Monitoring System and AG-608 Multi Blood Glucose Monitoring System, Andon Health Care Management System Software transfers data from the device's memory into a computer for enhanced data management.

6.0 Intended use

Andon Health Care Management System Software is an optional software accessory for use with the Andon blood glucose meters with data management capacities. When used with one of these meters, Andon Health Care Management System Software transfers data from the device's memory into a computer for enhanced data management. Andon Health Care Management System Software is intended for use in home and clinical settings via the internet to assist people with diabetes and their healthcare professionals in uploading, storing, analyzing, and communicating historical blood glucose test results and other biological statistics to support diabetes management. Andon Health Care System Software is not intended to provide treatment decisions or to be used as a substitute for professional healthcare judgment. All patient medical diagnoses and treatment are to be performed under the supervision and oversight of an appropriate healthcare professional.

The AG-608 MULTI Blood Glucose Monitoring System is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips. The AG-608 MULTI Blood Glucose Monitoring System is intended to be used by healthcare professionals for multiple patients in a professional healthcare setting as an aid in monitoring the effectiveness of diabetes control.

The AG-608 MULTI Blood Glucose Monitoring System should not be used for the diagnosis of or screening for diabetes, nor for neonatal use.

The AGS-1000 MULTI Blood Glucose Test Strips are for use with the AG-608 MULTI Blood Glucose Meter to quantitatively measure glucose (sugar) in fresh capillary whole blood sample.

This system should only be used with single-use, auto-disabling lancing devices.

The AG-608 Single Blood Glucose Monitoring System is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips. The AG-608 Single Blood Glucose Monitoring System is intended to be used by a single person and should not be shared.

The AG-608 Single Blood Glucose Monitoring System is intended for self testing outside the body (in vitro diagnostic use) by people with diabetes at home as an aid to monitor the effectiveness of diabetes control. The AG-608 Single Blood Glucose Monitoring System should not be used for the diagnosis of or screening for diabetes, nor for neonatal use.

The AGS-1000 Single Blood Glucose Test Strips are for use with the AG-608 Single Blood Glucose Meter to quantitatively measure glucose (sugar) in fresh capillary whole blood sample.

7.0 Summary comparing technological characteristics with predicate device

item	Andon	Clever
Indications for use	The ANDON Health Care Management System Software is intended for use in home and clinical settings as an aid for people with users and their health care professionals to review, analyze and evaluate the historical test results to support health management effectively.	The CLEVER CHEK Health Care System Software is intended for use in home and clinical settings as an aid for people with users and their health care professionals to review, analyze and evaluate the historical test results to support health management effectively.
Installation method	Exe file	Exe file
Package Contents	N/A	N/A
Capable of deleting results	Delete all results in meter	Delete all results in meter
Language capabilities	English, Spanish	English, Spanish, Traditional Chinese
Customizable schedule	N/A	N/A
Types of graphs etc.	Coordinates graph	Coordinates graph
Auto-detect COM port	Yes	Yes
System components	PC, USB cable, meter	PC, USB cable, meter
Software platform	Microsoft	Windows XP, Windows Vista, Windows 7
Hardware requirements	CPU: optimal at 1,2 GHz+ Main memory: optimal at	Personal computer with 400 megahertz (MHz) or higher

	256 MB+ RAM Disk space: optimal 200 MB+ free space - at least 100 MB Graphic resolution starting from 1024 x 768, CD-ROM drive, USB interface	processor clock speed recommended; the software does not run on Apple computers. _ Random access memory (RAM) of 64 megabytes (MB) or more Available hard disk space of 30 MB for running the program. Monitor with 1024 x 768 or higher resolution. Keyboard and mouse.
Technology	Visual Basic	Visual Basic
Performance specifications, including any testing	Read memories in meter. Delete all memories in meter Set time to meter Draw table and graph Print Set the personal information	Read memories in meter. Delete all memories in meter Set time to meter Draw table and graph Print Set the personal information

8.0 Performance summary

Testing of Andon Health Care Management System Software included system test and unit test.

9.0 Comparison to the predict device and the conclusion

Andon Health Care Management System Software is very similar with the predicted device

Clever Chek Health Care System Software, However, they are intended to use together with the different meters, their port are also different.

However, the test in this submission provides demonstration that these small differences do not raise any new questions of safety and effectiveness.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Andon Medical Co., Ltd
c/o Mr. Yi Liu
1st Floor Building C
No.3 Jin Ping Street
Nankai District
Tianjin, China 300190

Food and Drug Administration
10903 New Hampshire Avenue
Silver Spring, MD 20993

AUG 26 2011

Re: k102678
Trade name: Andon Health Care Management System Software
AG-608 Single Blood Glucose Monitoring System
AG-608 Multi Blood Glucose Monitoring System

Regulation Number: 21 CFR 862.1345
Regulation Name: Glucose Test System
Regulatory Class: Class II
Product Code: NBW, JQP
Dated: August 2, 2010
Received: August 2, 2010

Dear Mr. Liu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'CH', followed by a long horizontal line extending to the right.

Courtney Harper, Ph.D.
Director
Division of Chemistry and Toxicology
Office of *In Vitro* Diagnostic Device
Evaluation and Safety
Center for Devices and Radiological Health

Enclosure

Indication for Use

510(k) Number (if known): k102678

Device Name: Andon Health Care Management System Software

Indication For Use:

Andon Health Care Management System Software is an optional software accessory for use with the Andon blood glucose meters with data management capacities. When used with one of these meters, Andon Health Care Management System Software transfers data from the device's memory into a computer for enhanced data management.

Andon Health Care Management System Software is intended for use in home and clinical settings via the internet to assist people with diabetes and their healthcare professionals in uploading, storing, analyzing, and communicating historical blood glucose test results and other biological statistics to support diabetes management. Andon Health Care System Software is not intended to provide treatment decisions or to be used as a substitute for professional healthcare judgment. All patient medical diagnoses and treatment are to be performed under the supervision and oversight of an appropriate healthcare professional.

Prescription Use Yes
(21 CFR Part 801 Subpart D)

And/Or

Over the Counter Use Yes
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety (OIVD)



Division Sign-Off
Office of In Vitro Diagnostic Device
Evaluation and Safety

510(k) k102678

Indication for Use

510(k) Number (if known): K102678

Device Name: AG-608 MULTI Blood Glucose Monitoring System

Indication For Use:

The AG-608 MULTI Blood Glucose Monitoring System is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips. The AG-608 MULTI Blood Glucose Monitoring System is intended to be used by healthcare professionals for multiple patients in a professional healthcare setting as an aid in monitoring the effectiveness of diabetes control.

The AG-608 MULTI Blood Glucose Monitoring System should not be used for the diagnosis of or screening for diabetes, nor for neonatal use.

The AGS-1000 MULTI Blood Glucose Test Strips are for use with the AG-608 MULTI Blood Glucose Meter to quantitatively measure glucose (sugar) in fresh capillary whole blood sample.

This system should only be used with single-use, auto-disabling lancing devices.

Prescription Use ☒
(21 CFR Part 801 Subpart D)

And/Or

Over the Counter Use ☒
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety (OIVD)



Division Sign-Off
Office of In Vitro Diagnostic Device
Evaluation and Safety

510(k) K102678

Indication for Use

510(k) Number (if known): K102678

Device Name: AG-608 Single Blood Glucose Monitoring System

Indication For Use:

The AG-608 Single Blood Glucose Monitoring System is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips. The AG-608 Single Blood Glucose Monitoring System is intended to be used by a single person and should not be shared.

The AG-608 Single Blood Glucose Monitoring System is intended for self testing outside the body (in vitro diagnostic use) by people with diabetes at home as an aid to monitor the effectiveness of diabetes control. The AG-608 Single Blood Glucose Monitoring System should not be used for the diagnosis of or screening for diabetes, nor for neonatal use.

The AGS-1000 Single Blood Glucose Test Strips are for use with the AG-608 Single Blood Glucose Meter to quantitatively measure glucose (sugar) in fresh capillary whole blood sample.

Prescription Use ☒
(21 CFR Part 801 Subpart D)

And/Or

Over the Counter Use ☒
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety (OIVD)



Division Sign-Off
Office of In Vitro Diagnostic Device
Evaluation and Safety

510(k) K102678