

---

**APPENDIX A: 510(k) SUMMARY**

**Sponsor/Submitter:** Arstasis, Inc.  
1021 Howard Avenue, Suite C  
San Carlos, CA 94070

**Contact Person:** Debra Cogan  
Director, Regulatory & Clinical Affairs  
Phone: (650) 508-1549 x273  
Fax: (650) 594-4326

**Date of Submission:** September 20, 2010

**Device Trade Name:** Arstasis<sup>one</sup> Access System

**Common Name:** Catheter Introducer

**Device Classification:** Class II

**Regulation Number:** 21 CFR 870.1340

**Classification Name:** Catheter Introducer

**Product Code:** DYB

**Predicate Device:** Arstasis<sup>one</sup> Access System (K100615)

**Device Description:** Arstasis<sup>one</sup> is a device that is comprised of a sheath, anchor mechanism, shaft and handle with control features.

**Indications for Use:** Arstasis<sup>one</sup> Access System is intended to provide access for the percutaneous introduction of devices into the peripheral vasculature and to promote hemostasis at the arteriotomy site as an adjunct to manual compression. The Arstasis<sup>one</sup> Access System is indicated for use in patients undergoing diagnostic femoral artery catheterization procedures using 5F or 6F introducer sheaths.

**Technological Characteristics** Arstasis<sup>one</sup> is designed to create a shallow access path through the arterial wall for the guidewire to enter the vessel lumen.

**Summary of Substantial Equivalence:** Selective bench testing was performed on the subject device as follows: functionality, torque loading, flexibility and tensile strength of the core wire assembly and sheath.

Prior bench testing included deployment forces, flexibility, tensile, compression, and torque loading, were conducted on the predicate device. Additional prior testing included biocompatibility testing pursuant to ISO-10993-1, Biological Evaluation of Medical Devices Part 1: Evaluation and Testing (1995), preliminary animal studies (non-GLP) and cadaver assessments, as well as clinical

OCT 7 2010

investigations.<sup>1</sup> Multiple clinical evaluations were conducted. The short term safety and clinical performance of the device were established. The long term safety, as well as the ability to access and re-access, was retrospectively studied in a smaller cohort of patients.

In summary, the cumulative data provided herein demonstrates that the Arstasis<sup>one</sup> Access System is substantially equivalent to its predicate in providing access to the arterial lumen and facilitating the introduction and placement of devices into the peripheral vasculature and achievement of hemostasis.

---

<sup>1</sup> The preliminary Animal Studies and Cadaver Assessments were conducted using prototypes of a similar design and configuration.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

OCT 7 2010

Arstasis, Inc.  
c/o Ms. Su-Mien Chong  
Acting Vice President Research and Development  
1021 Howard Avenue, Suite C  
San Carlos, CA 94070

Re: K102728  
Trade/Device Name: Arstasis Access System  
Regulation Number: 21 CFR §870.1340  
Regulation Name: Catheter, Introducer  
Regulatory Class: Class II (Two)  
Product Code: DYB  
Dated: September 20, 2010  
Received: September 22, 2010

Dear Ms. Su-Mien Chong:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



 Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

---

**APPENDIX B: INDICATIONS FOR USE STATEMENT**510(k) Number (if known): K102728Trade Name: Arstasis<sup>one</sup> Access System

Common Name: Catheter Introducer

Indications For Use: The Arstasis<sup>one</sup> Access System is intended to provide access for the percutaneous introduction of devices into the peripheral vasculature and to promote hemostasis at the arteriotomy site as an adjunct to manual compression. The Arstasis<sup>one</sup> Access System is indicated for use in patients undergoing diagnostic femoral artery catheterization procedures using 5F or 6F introducer sheaths.

Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

*Donna R. K. [Signature]*  
(Division Sign-Off)  
Division of Cardiovascular Devices

Page 1 of 1

(Posted November 13, 2003)

510(k) Number K102728