

## 2. 510(k) SUMMARY

DEC 20 2010

**Submitted by:** Merete Medical GmbH  
Alt Lankwitz 102  
12247 Berlin, Germany

**FDA Registration Number:** 3002949614

**Contact Person:** Emmanuel Anapliotis  
Merete Medical, Inc.  
49 Purchase Street  
Rye, New York 10580  
Phone: 914 967 1532

**Proprietary Name:** TwistCut™ EndoSorb™ Screws

**Device Classification:** Smooth or threaded metallic bone fixation  
Fastener (888.3040)

**Product Code:** JDW

**Subsequent Product Code:** HWC

**Proposed Regulatory Class:** Class II

### Legally Marketed Devices To Which Substantial Equivalence is Claimed:

ProToe™ EndoSorb™ Small Hammer Toe Pin, Merete (K100414)  
Reunite Screws, Biomet (K992301)

### Intended Use:

The TwistCut™ EndoSorb™ Bone Screws are indicated for use in the presence of appropriate immobilization in the following procedures:

1. for metacarpal and phalangeal fusion and fracture
2. or repair of hallux valgus deformity (bunion).

### Device Description:

The TwistCut™ EndoSorb™ Bone Screws are made out of the absorbable EndoSorb™ material. EndoSorb™ is a polyester derivative of L-Lactic and glycolic acids. Poly(L-lactide-co-glycolide) material (PLGA) degrades and absorbs in vivo by hydrolysis into L-lactic and glycolic acids, which are then metabolized by the body. The screw head provides a hexagonal twist off portion to enable screw head to be flush with the bone.

**Substantial Equivalence:**

The TwistCut™ EndoSorb™ Bone Screws similar to legally marketed predicate device listed above in that it shares similar indications for use, is manufactured from similar materials and incorporate similar technological characteristics. Any differences have been found to have no obvious effect on the performance, function, or intended use of the prosthesis.

The following bench tests were performed:

- *In vitro* degradation test (determination of decreasing strength and inherent viscosity)
- Determination of inherent viscosity, glass transition temperature and crystallinity
- Brakeage strength and insertions torque
- Long-term stability of package and subject device
- Stability test under extreme transport conditions
- Temperature limit for storage
- A geometrical and a material comparison of the subject and the predicate devices were provided. This information served as a basis for a determination of substantial equivalence.

The predicates for the subject device are ProToe™ EndoSorb™ Small Hammer Toe Pin, Merete (K100414); Reunite Screws, Biomet (K992301).

**Software Documentation:**

No software is needed for the use of this device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

Merete Medical GmbH  
% Mr. Emmanuel Anapliotis  
President & CEO  
Alt Lankwitz 102  
12247 Berlin, Germany

DEC 20 2010

Re: K102777

Trade/Device Name: TwistCut EndoSorb Bone Screws  
Regulation Number: 21 CFR 888.3040  
Regulation Name: Smooth or threaded metallic bone fixation fastener  
Regulatory Class: II  
Product Codes: JDW, HWC  
Dated: September 21, 2010  
Received: September 24, 2010

Dear Mr. Anapliotis:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucml115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Mark N. Melkerson', with a stylized flourish at the end.

Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

Submission Date: 08 / 24 / 2010

1. INDICATIONS FOR USE STATEMENT

Indications for Use

FDA CDRH DMC  
SEP 24 2010  
Received

510(k) Number (if known): K102777

DEC 20 2010

Device Name: TwistCut™ EndoSorb™ Bone Screw

Indications for Use:

The TwistCut™ EndoSorb™ Bone Screws are indicated for use in the presence of appropriate immobilization in the following procedures:

1. for metacarpal and phalangeal fusion and fracture
2. or repair of hallux valgus deformity (bunion).

Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE) for M. Melkerson

(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices