

K103759

APR - 8 2011

5. 510(k) Summary

K103759
APR - 8 2011

510(k) Summary

Astra Tech Inc.
Atlantis™ Abutment in Zirconia for Osseospeed TX Profile Implant

ADMINISTRATIVE INFORMATION

510K Summary preparation date: December 20, 2010

Manufacturer Name: Astra Tech Inc.
590 Lincoln Street
Waltham, Massachusetts 02541
Telephone: 781-810-6462
Fax: 781-810-6719

Official Contact: Franklin Uyleman

~~Representative/Consultant:~~ ~~Betsy A. Brown~~
B.A. Brown and Associates Inc.
Telephone: 847-560-4406
Fax: 847-677-0177

DEVICE NAME AND CLASSIFICATION

Trade/Proprietary Name: Atlantis™ Abutment in Zirconia for Osseospeed TX
Profile Implant

Common Name: Endosseous dental implant abutment
21 CFR 872.3630

Product Code: NHA

Classification Panel: Dental Products Panel

Reviewing Branch: Dental Devices Branch

INTENDED USE

The Atlantis Abutment is intended for use with an endosseous implant to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be cemented or screw retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.

This device is compatible with the following manufacturers' implant systems:
The Atlantis Abutment in zirconia is compatible with the Astra Tech Osseospeed TX Profile 4.5, 5.0 and 5.0 S Implants.

INTENDED USE (continued)

Please note: This device may be used in an early load situation, but is dependent on the specific implant system and protocol used by the dental professional.

Highly angulated abutments on small implants are intended for the anterior region of the mouth only.

DEVICE DESCRIPTION

The devices covered in this submission are abutments which are placed into the dental implant to provide support for a prosthetic restoration. The subject abutments are indicated for cemented or screw retained restorations. The **Atlantis™ Abutment in Zirconia for Osseospeed TX Profile Implant** for the Astra Tech Osseospeed Profile 4.5, 5.0 and 5.0 S Implants is made of biocompatible material, yttria-stabilized tetragonal for the zirconia polycrystals (Y-TZP) (meets

~~ISO Standards 6972 & 13356). Zirconia may have variation in shade. The abutment screw is~~
made of Titanium grade Ti-6Al-4V ELI (meets ASTM Standard F-136). The zirconia abutments are placed over the implant shoulder and are mounted into the implant with a titanium screw.

EQUIVALENCE TO MARKETED DEVICE

Astra Tech Inc. demonstrated that, for purposes of the FDA's regulations of medical devices, the **Atlantis™ Abutment in Zirconia for Osseospeed TX Profile Implant** is substantially equivalent in indication and design principles to Astra Tech's Osseospeed TX Profile Implant cleared under K#080156 and The Atlantis Abutment in Zirconia for Astra Implants cleared under K#071946 which has been determined by FDA to be substantially equivalent to preamendment devices.

Table 1: Substantial Equivalence Summary

Technological Characteristics	Atlantis™ Abutment in Zirconia for Osseospeed TX Profile Implant	Astra Tech Osseospeed Profile Implant TX System
Material	-Biocompatible ceramic material	-comparable compatible titanium grade Ti-6A-4V ELI material
Performance characteristics	Allows the prosthesis to be cemented or screw retained to abutment. While the abutment screw is intended to secure the abutment to the endosseous implant.	Allows the prosthesis to be cemented or screw retained to abutment. While the abutment screw is intended to secure the abutment to the endosseous implant.
Intended Use	Intended for use with an endosseous implant to support a prosthetic device in a partially or completely edentulous patient. Intended for use to support single or multiple tooth prosthesis, in mandible or maxilla.	Intended for use with an endosseous implant to support a prosthetic device in a partially or completely edentulous patient. Intended for use to support single or multiple tooth prosthesis, in mandible or maxilla.

Summary of Non-clinical Testing

Static and fatigue compression testing was conducted on “worst case scenario” implant assemblies using Atlantis angled zirconia abutments with the Astra Tech Osseospeed TX Profile implant. Test results demonstrated that the Atlantis Abutments are compatible with the Astra Tech Osseospeed TX Profile implants and the implant system supported appropriate static and fatigue test loads demonstrating that the implant system performs as intended.

Conclusion for Substantial Equivalence:

The Atlantis™ Abutment in Zirconia for Osseospeed TX Profile Implant is substantially equivalent to Astra Tech’s Osseospeed TX Profile System (K080156) and the Atlantis Abutment in Zirconia for Astra Implant (K071946) predicate devices based on noted similarities in indication, manufacturing material, generated design principle and performance characteristics data.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room -WO66-G609
Silver Spring, MD 20993-0002

Astra Tech, Incorporated
C/O Ms. Betsy Brown
Regulatory Consultant
B.A. Brown & Associates
8944 Tamaroa Terrace
Skokie, Illinois 60076

APR - 8 2011

Re: K103759

Trade/Device Name: Atlantis™ Abutment in Zirconia for Osseospeed TX
Profile Implant

Regulation Number: 21 CFR 872.3630

Regulation Name: Endosseous Dental Implant Abutment

Regulatory Class: II

Product Code: NHA

Dated: March 21, 2011

Received: March 23, 2011

Dear Ms. Brown:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to

<http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "AW" followed by a flourish, and the word "for" written in a cursive script to the right.

Anthony Watson, B.S., M.S., M.B.A.
Director

Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Premarket Notification

Section 4: Page - 3

Indications for Use

510(k) Number (if known) K103759

Device Name: Atlantis™ Abutment in Zirconia for Osseospeed TX Profile Implant

Indication for Use:

~~The Atlantis Abutment is intended for use with an endosseous implant to support a prosthetic~~
device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be cemented or screw retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.

This device is compatible with the following manufacturers' implant systems:

The Atlantis Abutment in Zirconia is compatible with the Astra Tech Osseospeed TX Profile 4.5, 5.0 and 5.0 S Implants.

Please note: This device may be used in an early load situation, but is dependent on the specific implant system and protocol used by the dental professional.


Highly angled abutments on small diameter implants are intended for the anterior region of the mouth only.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR Over-the-Counter Use
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of Anesthesiology, General Hospital
Infection Control and Dental Devices
510(k) Number: K103759