



JUL 22 2011

510(k) SUMMARY

A summary of 510(k) safety and effectiveness information in accordance with the requirements of 21 CFR 807.92

Submitter Information	
Name	Biomet Manufacturing Corp.
Address	56 East Bell Drive Warsaw, IN 46581-0857
Phone number	(574) 267-6639
Fax number	(574) 371-1027
Establishment Registration Number	1825034
Name of contact person	Patricia Sandborn Beres Senior Regulatory Specialist
Date prepared	January 10, 2011
Name of device	
Trade or proprietary name	Maestro® Wrist Plating System
Common or usual name	<ul style="list-style-type: none"> • plate, fixation, bone • screw, fixation, bone
Classification name	<ul style="list-style-type: none"> • Single/multiple component metallic bone fixation appliances and accessories • Smooth or threaded metallic bone fixation fastener
Classification panel	Orthopedics
Regulation	<ul style="list-style-type: none"> • 21 CFR 888.3030 • 21 CFR 888.3040
Product Code(s)	<ul style="list-style-type: none"> • HRS • HWC
Legally marketed device(s) to which equivalence is claimed	K093761 - OptiLock® VL Distal Plating System K040908 - EBI® Distal Radius Plating System
Reason for 510(k) submission	New device
Device description	The Maestro® Wrist Plating System is comprised of anatomic plates in four styles: Volar, Dorsal, Radial Lateral and Ulnar. Both locking and non-locking screws in multiple lengths as well as pegs are designed for use with the plates. Plate sizing and contouring was developed through the use of IntelliFIT, a Biomet technology which uses contour analysis to map patterns in complex bone on cadaveric specimens to determine plate countouring. (Note, the software was used to determine a set of pre-defined plate sizes and is not used to create individual, patient matched plates.)

Mailing Address:
P.O. Box 587
Warsaw, IN 46581-0587
Toll Free: 800.348.9500
Office: 574.267.6639
Main Fax: 574.267.8137
www.biomet.com

Shipping Address:
56 East Bell Drive
Warsaw, IN 46582

K110271(2/2)

Intended use of the device	Bone fixation		
Indications for use	The Maestro® Wrist Plating System is indicated for fractures, fracture dislocations, osteotomies and non-unions of the distal radius and ulna.		
Summary of the technological characteristics of the device compared to the predicate			
Characteristic	New Device	Predicate Device*	
Plate Design	Volar, Dorsal, Lateral, Ulnar	K093761, K040908	
Plate Material	Stainless Steel ASTM F-138, F-139	K093761, K040908	
Plate Lengths	Length: 44-170mm	K093761, K040908	
Screw Design	Locking, Non-Locking, Peg	K093761, K040908	
Screw Material	Stainless Steel ASTM F-138, F-139	K093761, K040908	
Screw Dimensions	Diameter: 2.7 and 3.5mm Length: 10-30mm	K093761, K040908	
PERFORMANCE DATA			
SUMMARY OF NON-CLINICAL TESTS CONDUCTED FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE			
Performance Test Summary-New Device			
Characteristic	Standard/Test/FDA Guidance		Results Summary
Plate Strength	Engineering Analysis		Meet or exceed predicate
Comparative Performance Information Summary			
Characteristic	Requirement	New Device	Predicate Device*
Plate Strength	Meet or exceed predicate	Meet	K093761
SUMMARY OF CLINICAL TESTS CONDUCTED FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE AND/OR OF CLINICAL INFORMATION			
Clinical Performance Data/Information: None			
MAGNETIC RESONANCE (MR) ENVIRONMENT			
Biomet® has performed non-clinical Magnetic Resonance Imaging (MRI) studies on Plating Systems manufactured of 316L Stainless Steel per ASTM F-138. These Plating Systems are determined to be MR Conditional in accordance to ASTM F-2503-08 Standard Practice for Marking Devices and Other Items for Safety in the Magnetic Resonance Environment. MR Conditional refers to an item that has been demonstrated to pose no known hazards in a specified MR environment with specified conditions of use.			
CONCLUSIONS DRAWN FROM NON-CLINICAL AND CLINICAL DATA			
No mechanical or clinical testing was necessary for a determination of substantial equivalence. The results of engineering analysis indicated the devices performed within the intended use, did not raise any new safety and efficacy issues and were found to be substantially equivalent to the predicate devices.			

*Any statement made in conjunction with this submission regarding and/or a determination of substantial equivalence to any other product is intended only to relate to whether the product can be lawfully marketed without pre-market approval or reclassification and is not intended to be interpreted as an admission or any other type of evidence in patent infringement litigation. [Establishment Registration and Premarket Notification Procedures, Final Regulation, Preamble, August 23, 1977, 42 FR 42520 (Docket No. 76N-0355)]



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room -WO66-G609
Silver Spring, MD 20993-0002

Biomet Manufacturing Corporation
% Ms. Patricia Sandborn Beres
Senior Regulatory Specialist
56 East Bell Drive
Warsaw, Indiana 46581

JUL 22 2011

Re: K110271

Trade/Device Name: Maestro® Wrist Plating System

Regulation Number: 21 CFR 888.3030

Regulation Name: Single/multiple component metallic bone fixation appliance and accessories

Regulatory Class: Class II

Product Code: HRS, HWC

Dated: June 22, 2011

Received: June 24, 2011

Dear Ms. Beres:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Mark N. Melkerson *per Ch. Pr*
Director
Division of Surgical, Orthopedic,
and Restorative Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): _____

Device Name: Maestro® Wrist Plating System

Indications For Use:

The Maestro® Wrist Plating System is indicated for fractures, fracture dislocations, osteotomies and non-unions of the distal radius and ulna.


Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use NO
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of Surgical, Orthopedic,
and Restorative Devices

510(k) Number K110271