



K110862

510(k) Summary

APR 26 2011

Submitted By: Chris Stukel
Hollister Incorporated
2000 Hollister Drive
Libertyville, IL 60018
847-680-1000

Date Summary Prepared: March 24, 2011

Device Name: Classification Name- Urological catheter and accessories

Common/Usual Name- Catheter, Urethral

Proprietary Name- VaPro™ Plus Intermittent Catheter

Predicate Device: The VaPro Plus intermittent catheter is substantially equivalent to the following products:

| Product | 510(k) |
|---|---------|
| VaPro™ Intermittent Catheter | K090960 |
| Incare Advance/Advance Plus Intermittent Catheter | K013483 |

Device Description: The VaPro Plus intermittent catheter is a hydrophilic coated single use catheter to be used as a means of managing urinary incontinence by draining urine from the bladder. The catheter comes in a protective sleeve and is offered with a protective introducer tip as a way to shield the sterile catheter from bacteria in the distal urethra during insertion. The packaging contains a vapor strip that hydrates the catheter coating which then lubricates the catheter. The outer packaging was designed to facilitate access for those with limited dexterity. The catheter is connected to a collection bag for use when drainage into a suitable receptacle is not feasible or practical.

Intended Use: The VaPro Plus intermittent catheter is a flexible tubular device that is inserted through the urethra by male, female and pediatric patients who need to drain urine from the bladder.

Technological
Characteristics:

The table below summarizes the technological characteristics of the device as compared to the predicate devices.

| Characteristics | VaPro Plus Intermittent Catheter | VaPro Intermittent Catheter (K090960) | Incare Advance/ Advance Plus (K013483) |
|-------------------------|---|--|--|
| Intended Use | The VaPro Plus intermittent catheter is a flexible tubular device that is inserted through the urethra by male, female and pediatric patients who need to drain urine from the bladder. | The VaPro intermittent catheter is a flexible tubular device that is inserted through the urethra by male, female and pediatric patients who need to drain urine from the bladder. | Indicated for use by male, female and pediatric patients for the purpose of bladder drainage. The Catheter is a flexible tubular device that is inserted through the urethra and used to pass fluids to or from the urinary tract. |
| Condition of Use | Single Use | Single Use | Single Use |
| Prelubricated | Yes-by water vapor hydration | Yes-by water vapor hydration | Yes-by hydrogel |
| Ready to use | Yes | Yes | Yes |
| End Design | Catheter funnel attached to collection bag | Funnel | Advance – Funnel Advance Plus - Catheter funnel inserted into collection bag |
| Sterile | Yes | Yes | Yes |
| No touch design | Yes-contains sleeve | Yes-contains sleeve | Yes-contains sleeve |
| Lubricant | PVP Based (polyvinylpyrrolidone) Coating | PVP Based (polyvinylpyrrolidone) Coating | Hydrogel |
| Protective Tip | Yes | Yes | Yes |
| Collection Bag | Yes | No | No- Advance Yes- Advance Plus |

Performance Testing
Conclusions:

Biocompatibility testing was performed based on the United States Food and Drug Administration Office of Device Evaluation General program Memorandum #G95-1 and ISO 10993 biocompatibility testing standards. Results indicate compliance to the standard.

Product evaluation also supports device functionality.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration
10903 New Hampshire Avenue
Document Mail Center - WO66-G609
Silver Spring, MD 20993-0002

Ms. Christine Stukel
Sr. Global Regulatory Affairs Analyst
Hollister, Inc.
2000 Hollister Drive
LIBERTYVILLE IL 60048-3781

APR 26 2011

Re: K110862
Trade/Device Name: VaPro™ Plus Intermittent Catheter
Regulation Number: 21 CFR §876.5130
Regulation Name: Urological catheter and accessories
Regulatory Class: II
Product Code: GBM
Dated: March 24, 2011
Received: March 29, 2011

Dear Ms. Stukel:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related

adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in dark ink, appearing to read "Herbert P. Lerner, M.D.", with a stylized flourish at the end.

Herbert P. Lerner, M.D., Director (Acting)
Division of Reproductive, Gastro-Renal
and Urological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

K110862

510(k) Number (if known): K110862

Device Name: VaPro™ Plus Intermittent Catheter

Indications for Use:

The VaPro Plus intermittent catheter is a flexible tubular device that is inserted through the urethra by male, female and pediatric patients who need to drain urine from the bladder.

Prescription Use x
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

Division of Reproductive, Gastro-Renal, and
Urological Devices

510(k) Number K110862