

Pollogen™

Page 1 of 3

NOV 16 2011

## 510(K) Summary

**apollo™ system**

**510(k) Number K** 111026

**Applicant's Name:** Pollogen Ltd.  
6 Kaufman St.  
Gibor House, P.O.B. 50320  
Tel Aviv  
ISRAEL 68012  
Tel. (972)3-510-4110  
Fax (972)3-510-4112

**Contact Person:** Yoram Levy, Qsite  
31 Haavoda St.  
Binyamina, Israel 30500  
Tel (972)4-638-8837; Fax (972)4-638-0510  
[Yoram@qsite-med.com](mailto:Yoram@qsite-med.com)

**Trade Name:** *apollo™*

### **Summary**

**Preparation Date:** April 14, 2011

**Classification:** **Name:** Electrosurgical, cutting & coagulation device  
& accessories  
**Product Code:** GEI  
**Regulation No:** 21 CFR 878.4400  
**Class:** II  
**Panel:** General and Plastic Surgery

### **Device Description:**

The *apollo™* system delivers RF current into the skin to generate heat through electrical impedance in the dermis and subcutaneous layers. The system consists of:

- Main Unit (includes the Controller)
- Control Panel (User Interface)
- RF Generator
- Three Treatment Applicators
- Foot Switch
- Patient-Controlled Manual Switch

Pollogen™

K111026

Page 2 of 3

The three applicators differ in size and configuration and are indicated for treatment of various size facial areas. The operator can adjust treatment parameters such as the power level and treatment time from the user interface on the Main Unit. The Applicator is applied with low pressure and a rubbing/massaging technique (linear, circular, etc., depending on the area). The applicator should be moved continuously on the skin. No active cooling of the electrodes or the skin is required.

The RF power module provides RF energy to the selected applicator at a frequency of 1 MHz and a maximum output RMS power of 50 watts for applicators No. 1 and 2 and a maximum RMS power of 15 watts for applicator No. 3.

**Intended Use Statement:**

The *apollo*™ is a non-invasive device intended for use in Dermatologic and general surgical procedures for the non-invasive treatment of mild to moderate facial wrinkles and rhytides.

**Predicate Devices:** Substantial equivalence to the following predicate devices is claimed:

Device Name	510k No	Date of Clearance
Venus Concept Ltd. Venus Freeze	K100586	Nov 29, 2010
EndyMed Imagine TC Skin Treatment System	K083461	Jul 24, 2009

**Performance Standards**

*apollo*™ complies with the European Medical Directive 93/42/EEC concerning medical devices (Annex II) and with the following standards:

- *EN 60601-1* Medical Electrical Equipment-Part 1: General Requirements for Safety. Collateral Standard: Safety Requirements for Medical Electrical Systems.
- *IEC 60601-1-2* Medical Electrical Equipment Part 1-2: Collateral Standard: Electromagnetic Compatibility – Requirements and Tests.
- *ANSI AAMI 60601-2-2* - safety of high frequency surgical equipment.

A detailed description appears in Section 14.

**Summary of Clinical performance data**

Pollogen™

K111026

Page 3 of ③

The safety and efficacy of radiofrequency devices emitting energy with a frequency of 1 MHz and power of up to 50 W is well established in scientific research and clinical studies. Multiple studies with similar systems have shown safety in dermatologic therapy and the devices were cleared by the FDA for therapy of wrinkles and rhytides.

Due to the comprehensive animal and clinical study performed in scientific research and published in the literature, and since the power and frequency of the *apollo*™ are well within the previously cleared values, Pollogen believes that animal and clinical studies are not required to determine the safety and efficacy of the device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

NOV 16 2011

Pollogen Ltd.  
% Qsite  
Mr. Yoram Levy  
31 Haavoda Street  
Binyamina, Israel 30500

Re: K111026  
Trade/Device Name: apollo<sup>TM</sup>  
Regulation Number: 21 CFR 878.4400  
Regulation Name: Electrosurgical cutting and coagulation device and accessories  
Regulatory Class: Class II  
Product Code: GEI  
Dated: November 01, 2011  
Received: November 04, 2011

Dear Mr. Levy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act

or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.


If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Mark N. Melkerson', is written over the typed name and title. The signature is fluid and cursive.

Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

Pollogen™

K 111026

### INDICATIONS FOR USE STATEMENT

510(k) Number (if known):

Device Name:

*apollo*™

Indications for Use:

The *apollo*™ is a non-invasive device intended for use in Dermatologic and General Surgical procedures for the non-invasive treatment of mild to moderate facial wrinkles and rhytides.

Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE -CONTINUE ON ANOTHER PAGE IF  
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-off)

Division of General and Plastic Surgery Devices

510(k) Number

Neil R. Ogden for mxm

(Division Sign-Off)

Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number

K 111026