

k 111454

**510(K) Summary**

JUL 5 2012

Submitter: Cynosure, Inc.  
5 Carlisle Road  
Westford, MA 01886

Contact: George Cho  
Sr. VP of Medical Technology and Regulatory Affairs

Date Summary Prepared: May 24, 2011

Device Trade Name: Illuminage Diode Laser

Common Name: Medical Laser System

Classification Name: Instrument, Surgical, Powered, Laser  
79-GEX  
21 CFR 878.4810

Equivalent Device: Cynosure Affirm DO Diode Laser

Device Description: Illuminage is a diode laser composed of a handpiece, base unit and a charger unit.

Intended Use: The Illuminage Diode Laser is indicated for use in the treatment of periorbital and perioral wrinkles.

Comparison: The Illuminage Diode laser has the same indications for use, the same principle of operation, the same wavelength and the same technology as the predicate devices.

Nonclinical Performance Data: none

Clinical Performance Data: none

Conclusion: The Illuminage Diode laser is a safe and effective device for the intended uses.

Additional Information: none



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room –WO66-G609  
Silver Spring, MD 20993-0002

Cynosure, Incorporated  
% Mr. George Cho  
Senior Vice President of Medical Technology and Regulatory Affairs  
5 Carlisle Road  
Westford, Massachusetts 01886

JUL 5 2012

Re: K111454  
Trade/Device Name: Illuminage Diode Laser  
Regulation Number: 21 CFR 878.4810  
Regulation Name: Laser surgical instrument for use in general and plastic surgery and in dermatology  
Regulatory Class: II  
Product Code: OHS, GEX  
Dated: June 29, 2012  
Received: July 02, 2012

Dear Mr. Cho:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

  
Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K 111454

Device Name: Illuminage Laser

Indications For Use: The Illuminage laser device is indicated for use in the treatment of periorbital and perioral wrinkles.

Prescription Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D)

OR

Over-The-Counter Use X  
(Part 21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Neil R. P. D. for mxm  
(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number K 111454