

K112040

DEC 15 2011

 SPECIAL 510(k): Arthrex BioComposite TransFix

# 1 510(k) Summary of Safety and Effectiveness

<b>Date Summary Prepared</b>	December 14, 2011
<b>Purpose of Submission</b>	To obtain clearance of new <i>Arthrex BioComposite TransFix</i> devices.
<b>Manufacturer/Distributor /Sponsor</b>	Arthrex, Inc. 1370 Creekside Boulevard Naples, FL 34108-1945 USA
<b>510(k) Contact</b>	Courtney Smith Manager, Regulatory Affairs Arthrex, Inc. 1370 Creekside Boulevard Naples, FL 34108-1945 USA Telephone: 239/643.5553, ext. 1270 Fax: 239/598.5508 Email: <a href="mailto:courtney.smith@arthrex.com">courtney.smith@arthrex.com</a>
<b>Trade Name</b>	<i>BioComposite TransFix</i>
<b>Common Name</b>	Pin, Fixation, Bone
<b>Product Code - Classification Name</b>	HWC – Screw, fixation, bone MAI – Fastener, fixation, biodegradable 21 CFR 888.3030: Single/multiple component metallic bone fixation appliances and accessories.
<b>Predicate Devices</b>	K062466 ( <i>Arthrex Interference Screw</i> ): Arthrex Bio-TransFix
<b>Device Description and Intended Use</b>	<p>The <i>Arthrex BioComposite TransFix</i> is very similar to the predicate devices in diameter, and length. The difference lies in the biodegradable material used to manufacture the new model and the addition of a longer device. The implants will be manufactured from PLDLA combined with <i>biphasic</i> calcium phosphate.</p> <p>The Arthrex BioComposite TransFix is intended to be used for fixation of tissue including ligament or tendon to bone and bone tendon bone. Fixation with the BioComposite TransFix is appropriate for surgeries of the knee, shoulder, elbow, ankle, foot, and hand/wrist where the sizes offered are patient appropriate. Specifically:</p> <p><i>Shoulder:</i> Rotator Cuff Repairs, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair, Deltoid Repair, Capsular Shift or Capsulolabral Reconstruction</p> <p><i>Foot/Ankle:</i> Lateral Stabilization, Medial Stabilization, Achilles</p>

K112040

	<p>Tendon Repair, Hallux Valgus Reconstruction, Midfoot Reconstruction, Metatarsal Ligament Repair, Flexor Hallucis Longus for Achilles Tendon reconstruction, tendon transfers in the foot and ankle</p> <p><i>Knee:</i> Anterior Cruciate Ligament Repair, Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis</p> <p><i>Elbow:</i> Biceps Tendon Reattachment, Ulnar or Radial Collateral Ligament Reconstruction</p> <p><i>Hand/Wrist:</i> Scapholunate Ligament Reconstruction, Ulnar Collateral Ligament Reconstruction, Radial Collateral Ligament Reconstruction, Carpometacarpal joint arthroplasty (basal thumb joint arthroplasty), Carpal Ligament Reconstructions and repairs, tendon transfer in the hand/ wrist</p>
<b>Substantial Equivalence Summary</b>	<p>The <b>Arthrex BioComposite TransFix</b> is substantially equivalent to the <b>Arthrex Bio-TransFix</b> predicates, in which the basic features and intended uses are the same. Any differences between the <b>BioComposite TransFix</b> and the predicates are considered minor and do not raise questions concerning safety and effectiveness.</p> <p>The proposed devices are composed of PLDLA/<i>biphasic</i> calcium phosphate and are ISO 10993-1 tested to confirm the material's biocompatibility.</p> <p>The submitted mechanical testing data demonstrated that the degradation shear strength, through 16 weeks, of the proposed devices meets or exceeds the shear strength of the predicate device.</p> <p>Based on the indication for use, technological characteristics, and the comparison to the predicate devices, Arthrex, Inc. has determined that the <b>Arthrex BioComposite TransFix</b> are substantially equivalent to currently marketed predicate devices.</p>



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

Arthrex, Incorporated  
% Ms. Courtney Smith  
Manager, Regulatory Affairs  
1370 Creekside Boulevard  
Naples, Florida 34108-1945

DEC 15 2011

Re: K112040

Trade/Device Name: Arthrex Biocomposite Transfix  
Regulation Number: 21 CFR 888.3030  
Regulation Name: Single/multiple component metallic bone fixation appliances and accessories  
Regulatory Class: Class II  
Product Code: MAI, HWC  
Dated: November 16, 2011  
Received: November 17, 2011

Dear Ms. Smith:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.


Sincerely yours,

A handwritten signature in black ink, appearing to read 'Mark N. Melkerson', is written over the typed name.

Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic,  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K112040

Arthrex  SPECIAL 510(k): Arthrex BioComposite TransFix

## Indications for Use Form

### Indications for Use

510(k) Number (if known): K112040

Device Name: Arthrex BioComposite TransFix

#### Indications For Use:

The Arthrex BioComposite TransFix is intended to be used for fixation of tissue including ligament or tendon to bone and bone tendon bone. Fixation with the BioComposite TransFix is appropriate for surgeries of the knee, shoulder, elbow, ankle, foot, and hand/wrist where the sizes offered are patient appropriate. Specifically;

*Shoulder:* Rotator Cuff Repairs, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair, Deltoid Repair, Capsular Shift or Capsulolabral Reconstruction

*Foot/Ankle:* Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus Reconstruction, Midfoot Reconstruction, Metatarsal Ligament Repair, Flexor Hallucis Longus for Achilles Tendon reconstruction, tendon transfers in the foot and ankle

*Knee:* Anterior Cruciate Ligament Repair, Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis

*Elbow:* Biceps Tendon Reattachment, Ulnar or Radial Collateral Ligament Reconstruction

*Hand/Wrist:* Scapholunate Ligament Reconstruction, Ulnar Collateral Ligament Reconstruction, Radial Collateral Ligament Reconstruction, Carpometacarpal joint arthroplasty (basal thumb joint arthroplasty), Carpal Ligament Reconstructions and repairs, tendon transfer in the hand/ wrist

Prescription Use ☒ AND/OR Over-The-Counter Use ☐

(Per 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)

Division of Surgical, Orthopedic,  
and Restorative Devices

PAGE 1 of 1

510(k) Number K112040