

JAN 30 2012

## PREMARKET NOTIFICATION [510(k)] Summary

This Summary of Safety and Effectiveness is prepared in accordance with 21 CFR Part 807.92(c).

### 1. Company Name:

Chison Medical Imaging Co., Ltd.

No.8, Xiang Nan Road, Shuo Fang, New District, Wuxi, China 214142

Contact: Ms. Ruoli Mo

Tel: +86-510-85311707, 85310593

Fax: +86-510-85310726

U.S. Agent:

Leiker Regulatory & Quality Consulting

7263 Cronin Circle

Dublin, CA 94568

Contact: Bob Leiker

Tel: (925) 556-1302 Fax: (866) 718-3819

### 2. Device Name: SONO TOUCH Series (Portable) Diagnostic Ultrasound System

**Common/Usual Name:** Diagnostic Ultrasound System with Accessories

**Classification:** Regulatory Class: II

Review Category: Tier II

Classification Name	21 CFR Section	Product Code
Ultrasonic pulsed echo imaging system	892.1560	90-IYO
Diagnostic ultrasonic transducer	892.1570	90-ITX

### 2. Marketed Device:

A6 Portable Ultrasonic Diagnostic System

### 3. Device Description:

The SONO TOUCH Series ultrasound system is an integrated preprogrammed ultrasound imaging system, capable of producing high detail resolution intended for clinical diagnostic imaging applications.

The CHISON ultrasound system is configured as a portable model (SONO TOUCH Series). These systems are designed with the latest technology, using the same quality procedure as ultrasound systems, which have been available in the market for years.

This CHISON ultrasound system is a general purpose, software controlled, diagnostic ultrasound system. Its basic function is to acquire ultrasound echo data and display the image in B-Mode (including Tissue Harmonic Imaging), M-Mode, or a combination of these modes.

The SONO TOUCH Series Models, have been designed to meet the following product safety standards: NEMA UD 2, IEC 60601-1, IEC 60601-1-2, IEC 60601-2-37, IEC 10993-1.

#### **4. Indications for Use:**

The system is a general-purpose ultrasonic imaging instrument intended for use by a qualified physician for evaluation of Abdomen, Cardiac, Small Organ (Thyroid, parathyroid, parotid, submaxillary gland, testes and breast.), Peripheral Vascular, Transvaginal, Transrectal, Musculo-skeletal (Conventional and Superficial), Pediatric, Fetal, OB/Gyn and Urology.

#### **Comparison to Predicate Device:**

The SONO TOUCH Series Models is of comparable type and substantially equivalent to the A6 Portable Ultrasonic Diagnostic System. All systems transmit ultrasonic energy into patients, then perform post processing of received echoes to generate on-screen display of anatomic structures and fluid flow within the body, and have the same intended uses and basic operating modes as the predicate device. All systems allow for specialized measurements of structures and flow, and calculations.

#### **5. Conclusion:**

The SONO TOUCH Series Models is substantially equivalent in safety and effectiveness to the predicate systems. The systems are intended for diagnostic ultrasound imaging and fluid flow analysis. The systems have the same gray-scale. The systems have acoustic output levels below the applicable FDA limits. The systems are designed to applicable electrical and physical safety standards.

**End of 510(k) Summary.**



Food and Drug Administration  
10903 New Hampshire Avenue  
Silver Spring, MD 20993

Chison Medical Imaging CO., Ltd.  
% Mr. Bob Leiker  
U.S. Agent  
Leiker Regulatory & Quality Consulting  
7263 Cronin Circle  
DUBLIN CA 94568

MAR - 8 2012

Re: K112539  
Trade/Device Name: SONOTOUCH Series (Portable) Diagnostic Ultrasound System  
Regulation Number: 21 CFR 892.1560  
Regulation Name: Ultrasonic pulsed echo imaging system  
Regulatory Class: II  
Product Code: IYO, ITX  
Dated: January 14, 2012  
Received: January 18, 2012

Dear Mr. Leiker:

This letter corrects our substantially equivalent letter of February 1, 2012.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the SONOTOUCH Series (Portable) Diagnostic Ultrasound System, as described in your premarket notification:

Transducer Model Number

C3, 3.5 MHz Convex Array  
MC3, 3.0MHz Micro-convex Array  
V6, 6.0MHz Micro-convex Array

L7M, 7.5MHz Linear Array  
L7S, 7.5MHz Linear Array  
R7, 7.5MHz Linear Array

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

If you have any questions regarding the content of this letter, please contact Shahram Vaezy, Ph.D. at (301) 796-6242.

Sincerely Yours,

A handwritten signature in black ink, appearing to read "Janine M. Morris".

Janine M. Morris  
Acting Director  
Division of Radiological Devices  
Office of In Vitro Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

## Diagnostic Ultrasound Indications For Use

### 1.3 Indications for Use

The device is a general-purpose ultrasonic imaging instrument intended for use by a qualified physician for evaluation of Abdomen, Cardiac, Small Organ (Thyroid, parathyroid, parotid, submaxillary gland, testes and breast.), Peripheral Vascular, Transvaginal, Transrectal, Musculo-skeletal (Conventional and Superficial), Pediatric, Fetal, OB/Gyn and Urology.

Prescription Use ☒  
(Part 21.CFR 801 Subpart D)

AND/OR

Over-The-Counter Use  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Mary S. Patel  
(Division Sign-Off)

Division of Radiological Devices  
Office of In Vitro Diagnostic Device Evaluation and Safety

510K

K112539

## Diagnostic Ultrasound Indications For Use

**System:** SONO TOUCH Series Diagnostic Ultrasound Systems  
Diagnostic Ultrasound Pulsed Echo System

**Intended Use:** Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Other* Combined
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal	N	N					Note 1
	Abdominal	N	N					Note 1
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric	N	N					Note 1
	Small Organ <sup>(M)</sup> (Specify)	N	N					Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal	N	N					Note 1
	Trans-vaginal	N	N					Note 1
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)	N	N					Note 1
	Musculo-skeletal (Superficial)	N	N					Note 1
	Intravascular							
	Other (Urology)	N	N					Note 1
	Other (Ob/GYN)	N	N					Note 1
Cardiac	Cardiac Adult							
	Cardiac Pediatric	N	N					Note 1
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
Peripheral Vessel	Other (Specify)							
	Peripheral vessel	N	N					Note 1
	Other (Specify)							

N = new indication; P = previously cleared by FDA;

E = added under this appendix

Note 1: B/M

Comments:

Small Organ: Thyroid, parathyroid, parotid, submaxillary gland, testes and breast

Additional Comments:

Prescription Use ☒  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use ☐  
(21 CFR 801, Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

*Mary S. [Signature]*  
(Division Sign-Off)  
Division of Radiological Devices  
Office of In Vitro Diagnostic Device Evaluation and Safety

510K

*K112539*

System: SONO TOUCH Series Ultrasound Systems

Transducer: C3, 3.5 MHz Convex Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Other* Combined
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal	N	N					Note 1
	Abdominal	N	N					Note 1
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric							
	Small Organ <sup>(1)</sup> (Specify)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)							
	Musculo-skeletal (Superficial)							
	Intravascular							
	Other (Urology)	N	N					Note 1
	Other (Ob/GYN)	N	N					Note 1
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
	Other (Specify)							
Peripheral Vessel	Peripheral vessel							
	Other (Specify)							

N = new indication; P = previously cleared by FDA;

E = added under this appendix

Note 1: B/M

Comments:

Small Organ: Thyroid, parathyroid, parotid, submaxillary gland, testes and breast

Additional Comments:

Prescription Use ☒

AND/OR

Over-The-Counter Use ☐

(Part 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRII, Office of Device Evaluation (ODE)

*Mary S Patel*

(Division Sign-Off)

Radiological Devices

Office of Device Evaluation and Safety

K112539

System: SONO TOUCH Series Ultrasound Systems  
 Transducer: MC3, 3.0MHz Micro-convex Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Other* Combined
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal							
	Abdominal	N	N					Note 1
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric	N	N					Note 1
	Small Organ <sup>(1)</sup> (Specify)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)							
	Musculo-skeletal (Superficial)							
	Intravascular							
	Other (Urology)							
	Other (Ob/GYN)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric	N	N					Note 1
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
	Other (Specify)							
Peripheral Vessel	Peripheral vessel							
	Other (Specify)							

N = new indication; P = previously cleared by FDA;

E = added under this appendix

Note 1: B/M

Comments:

Small Organ: Thyroid, parathyroid, parotid, submaxillary gland, testes and breast

Additional Comments:

Prescription Use √

AND/OR

Over-The-Counter Use

(Part 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Mary S. Patel  
 (Division Sign-Off)

Division of Radiological Devices

Office of In Vitro Diagnostic Device Evaluation and Safety

Section 1.3

6118539

Indications For Use

Page 4 of 8

System: SONO TOUCH Series Ultrasound Systems

Transducer: V6, 6.0MHz Micro-convex Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Other* Combined
Ophthalmic Fetal Imaging & Other	Ophthalmic							
	Fetal							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric							
	Small Organ <sup>(1)</sup> (Specify)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal	N	N					Note 1
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)							
	Musculo-skeletal (Superficial)							
	Intravascular							
Cardiac	Other (Urology)	N	N					Note 1
	Other (Ob/GYN)	N	N					Note 1
	Cardiac Adult							
	Cardiac Pediatric							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
Peripheral Vessel	Intra-cardiac							
	Other (Specify)							
	Peripheral vessel							
	Other (Specify)							

N = new indication; P = previously cleared by FDA;

E = added under this appendix

Note 1: B/M

Comments:

Small Organ: Thyroid, parathyroid, parotid, submaxillary gland, testes and breast

Additional Comments:

Prescription Use ☒

AND/OR

Over-The-Counter Use

(Part 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

*Mary S. Patel*  
 (Division Sign-Off)  
 Division of Radiological Devices  
 Office of In Vitro Diagnostic Device Evaluation and Safety

Section 1.3

Indications For Use

Page 5 of 8

System: SONO TOUCH Series Ultrasound Systems  
 Transducer: L7M, 7.5MHz Linear Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Other* Combined
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric	N	N					Note 1
	Small Organ <sup>(1)</sup> (Specify)	N	N					Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)	N	N					Note 1
	Musculo-skeletal (Superficial)	N	N					Note 1
	Intravascular							
	Other (Urology)							
	Other (Ob/GYN)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
Peripheral Vessel	Other (Specify)							
	Peripheral vessel	N	N					Note 1
	Other (Specify)							

N = new indication; P = previously cleared by FDA;

E = added under this appendix

Note 1: B/M

Comments:

Small Organ: Thyroid, parathyroid, parotid, submaxillary gland, testes and breast

Additional Comments:

Prescription Use ☒  
 (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use ☐  
 (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRL, Office of Device Evaluation (ODE)

*Mary S. Patel*  
 (Division Sign-Off)

Division of Radiological Devices  
 Office of In Vitro Diagnostic Device Evaluation and Safety

510K  
 Section 1.3

6112539

Indications For Use

System: SONO TOUCH Series Ultrasound Systems  
 Transducer: L7S, 7.5MHz Linear Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD Note 3	Color Doppler	Power (Amplitude) Doppler	Other* Combined
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric	N	N					Note 1
	Small Organ <sup>(1)</sup> (Specify)	N	N					Note 1
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)	N	N					Note 1
	Musculo-skeletal (Superficial)	N	N					Note 1
	Intravascular							
	Other (Urology)							
	Other (Ob/GYN)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
Peripheral Vessel	Other (Specify)							
	Peripheral vessel	N	N					Note 1
	Other (Specify)							

N = new indication; P = previously cleared by FDA;

E = added under this appendix

Note 1: B/M

Comments:

Small Organ: Thyroid, parathyroid, parotid, submaxillary gland, testes and breast

Additional Comments:

Prescription Use ☒  
 (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use ☐  
 (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

*Mary Spatel*  
 (Division Sign-Off)  
 Division of Radiological Devices  
 Office of In Vitro Diagnostic Device Evaluation and Safety  
 Section 1.3 6112539

Indications For Use

System: SONO TOUCH Series Ultrasound Systems  
 Transducer: R7, 7.5MHz Linear Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD Note 3	Color Doppler	Power (Amplitude) Doppler	Other* Combined
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric							
	Small Organ <sup>(1)</sup> (Specify)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal	N	N					Note 1
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)							
	Musculo-skeletal (Superficial)							
	Intravascular							
	Other (Urology)	N	N					Note 1
	Other (Ob/GYN)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
	Other (Specify)							
Peripheral Vessel	Peripheral vessel							
	Other (Specify)							

N = new indication; P = previously cleared by FDA;

E = added under this appendix

Note 1: B/M

Comments:

Small Organ: Thyroid, parathyroid, parotid, submaxillary gland, testes and breast

Additional Comments:

Prescription Use ☒   
 (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use ☐   
 (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

*Mary S. Postel*  
 (Division Sign-Off)

Division of Radiological Devices  
 Office of In Vitro Diagnostic Device Evaluation and Safety

510K

*K112539*