



510(k) Notification:  
MEDISS Reprocessed  
Arthroscopic Shavers and Burs

Section 5: 510(K) Summary K113028 DEC 20 2011

<b>Submitter/ Owner</b>	MEDISS 2747 SW 6th St. Redmond, OR 97756	
<b>Contact Name</b>	Joyce Elkins Director of Quality Assurance and Regulatory Affairs P: 541-923-3310 F: 541-923-3375 E: JElkins@MEDISS.com	
<b>Date Prepared</b>	October 7, 2011	
<b>Device Names</b>	Proprietary Name: MEDISS Reprocessed Arthroscopic Shavers and Burs Common Name: Arthroscope and Accessories	
<b>Classification</b>	Arthroscope, Class II, 21 CFR 888.1100, product code HRX	
<b>Predicate Devices</b>	K940075	Arthrex Shaver Blade Set
	K012536	SterilMed Reprocessed Powered Arthroscopic Accessories
	K012667	Reprocessed Arthroscopic Shavers
<b>Device Description</b>	Arthroscopic devices reprocessed by MEDISS include burs and blades at the end of a long rod that rotates within a hollow stainless steel housing. The housing has an opening on one side of the distal end, allowing the cutting tip to resect tissue while protecting adjacent material with the housing on the opposite side of the bur or blade. This system attaches to a motorized hand piece that drives the internal bur or blade inside the outer housing and provides suction to pull the cut tissue away from the surgical site. Devices are provided sterile.	
<b>Intended Use</b>	MEDISS Reprocessed Shavers and Burs are powered arthroscopic accessories intended for use in orthopedic joint surgery.	
<b>Technological Characteristics</b>	The technological characteristics of the subject devices are substantially equivalent to the predicate devices listed in this submission. The subject devices have the same functionality and very similar indications as the predicate devices.	



510(k) Notification:  
MEDISSS Reprocessed  
Arthroscopic Shavers and Burs

---

<b>Performance Testing</b>	The functional characteristics of the subject devices have been evaluated and found to be equivalent to the predicate devices. The patient contacting materials are either identical to the materials of the predicate devices, or have been fully evaluated for biocompatibility and functionality.
<b>Conclusion</b>	Based on comparison of the indications for use, technological characteristics, and performance data to the predicate devices, MEDISSS Reprocessed Arthroscopic Shavers and Burs have been shown to be substantially equivalent to the predicate devices.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

Medisiss  
% Ms. Joyce Elkins  
Director of Quality Assurance and  
Regulatory Affairs  
2747 SW 6<sup>th</sup> Street  
Richmond, Oregon 97756

DEC 20 2011

Re: K113028  
Trade/Device Name: Medisiss Reprocessed Arthroscopic Shavers and Burs  
Regulation Number: 21 CFR 888.1100  
Regulation Name: Arthroscope  
Regulatory Class: Class II  
Product Code: HRX  
Dated: November 22, 2011  
Received: December 02, 2011

Dear Ms. Elkins:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Mark N. Melkerson

Director

Division of Surgical, Orthopedic  
and Restorative Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure



510(k) Notification:  
MEDISSS Reprocessed  
Arthroscopic Shavers and Burs  
All Information on This Page is Confidential

---

#### Section 4: Indications for Use

---

510(k) Number: ~~TBD~~ K113028

Device Name: MEDISSS Reprocessed Arthroscopic Shavers and Burs

Indications For Use:

MEDISSS Reprocessed Shavers and Burs are powered arthroscopic accessories intended for use in orthopedic joint surgery.

Prescription Use   X    
(Part 21 CFR 801.109)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Neil R. Ogden for nkm  
(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number K113028