

Submitter Name Pie Medical Imaging BV  
Submitter Address Becanusstraat 13 D; 6216 BX Maastricht, The Netherlands  
Contact Person: Florie Daniels, Product Registration Coordinator  
Phone Number +31 43 32 81 328  
Fax Number +31 43 32 81 329  
Email Address Florie.Daniels@pie.nl  
Preparation Date 22 December 2011

JAN 30 2012

Trade Name CAAS A-Valve  
Common Name Cardiovascular Angiography Analysis System  
CAAS A-Valve  
Regulation Class Class II (21 CFR, part 892.2050, LLZ)

Predicate Devices CAAS QxA 3D, cleared under K100292  
IC-Pro, cleared under K110256

Device Description The CAAS A-Valve is a stand-alone software application intended to run on a PC with a Windows operating system. The CAAS A-Valve is designed for objective, accurate and reproducible assessment of the aortic root geometry from a set of angiographic X-ray images from different projections.

On each of the 2D images the aortic root is segmented. The segmented 2D aortic root contours in each image are used to generate a 3D reconstruction of the aortic root. A number of analysis results can be calculated:

1. Optimal C-arm projection to optimize visualization during treatment;
2. Dimensions of the aortic root.

Results are corrected for out-of-plane magnification and foreshortening errors.

Intended Use CAAS A-valve has been developed to support the interventionalist during or in preparation of treatment of the aortic root. CAAS A-Valve quantifies dimensions of the aortic root and assists in C-arm positioning based on a set of angiographic images. The software is used by or under supervision of a cardiologist or radiologist.

Indications for use:

CAAS A-Valve has been developed to support the interventionalist during or in preparation of treatment of the aortic root. Based on a set of angiographic X-ray images an analysis is performed:

- To assist in C-arm projection selection to optimize visualization during treatment;
- To calculate dimensions of the aortic root corrected for out-of-plane magnification and foreshortening errors.

Performance Data CAAS A-Valve is developed, tested, validated and produced under the same Quality Assurance system applicable to the development and production of products currently marketed by Pie Medical Imaging.

Substantial Equivalence The intended use and technological characteristics of CAAS A-Valve are substantial equivalent to a combination of the intended use and technological characteristics of the predicate devices.

## Conclusion

The testing reported in this 510(k) establishes that CAAS A-Valve is substantial equivalent to a combination of predicate devices and is safe and effective for its intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room – WO66-G609  
Silver Spring, MD 20993-0002

Ms. Florie Daniels  
Product Registration Coordinator  
Pie Medical Imaging BV  
Becanusstraat 13D  
6216 BX MAASTRICHT LIMBURG  
THE NETHERLANDS

JAN 30 2012

Re: K113076  
Trade/Device Name: CAAS A-Valve  
Regulation Number: 21 CFR 892.2050  
Regulation Name: Picture archiving and communications system  
Regulatory Class: II  
Product Code: LLZ  
Dated: December 22, 2011  
Received: December 28, 2011

Dear Ms. Daniels:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of

medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely Yours,

A handwritten signature in black ink, appearing to read "Mary S. Pastel", is written over a horizontal line.

Mary S. Pastel, Sc.D.  
Director  
Division of Radiological Devices  
Office of In Vitro Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

## Indications for Use Form

510(k) Number (if known): K113076

Device Name: CAAS A-Valve

Indications for Use:

CAAS A-Valve has been developed to support the interventionalist during or in preparation of treatment of the aortic root. Based on a set of angiographic X-ray images an analysis is performed:

- To assist in C-arm projection selection to optimize visualization during treatment;
- To calculate dimensions of the aortic root corrected for out-of-plane magnification and foreshortening errors.

Prescription Use   ✓    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use             
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Mary Spital

Division Sign-Off  
Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k) K113076