

## Section 6

510(k) Summary

---

**6. 510(k) Summary**

This 510(k) summary information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

**APPLICANT:** Arges Imaging, Inc.  
**DATE PREPARED:** September 9, 2011  
**CONTACT PERSON:** Rebecca K Pine  
129 N Hill Ave  
Pasadena, CA 91106  
Phone: 760.809.5178  
Fax: 760.290.3216  
**TRADE NAME:** Arges Imaging Contrast Spray  
**CLASSIFICATION NAME:** Impression Material  
**DEVICE CLASSIFICATION:** Class II  
**REGULATION NUMBER:** 872.3660  
**PRODUCT CODE:** ELW  
**PREDICATE DEVICES:** CEREC Opti Spray (K080882)

**Substantially Equivalent To:**

The Arges Imaging Contrast Spray is substantially equivalent in intended use, principal of operation and technological characteristics to the CEREC Opti Spray.

**Description of the Device Subject to Premarket Notification:**

The Arges Imaging Contrast Spray is a pigment suspension in ethanol with a fluorinated hydrocarbon propellant. The Arges Imaging Contrast Spray is a coating medium which is applied to teeth and surrounding areas. A 3D impression of the coated dentition is taken using a commercially available CAD/CAM system. Upon completion of the impression scan, the patient's mouth is rinsed and suctioned to remove the Arges Imaging Contrast Spray.

**Indication for Use:**

Arges Imaging Contrast is indicated as a coating medium for optical impressions. It aids in intraoral topographical recordings of prepared and unprepared teeth and their surroundings.

**Technological Characteristics:**

The Arges Imaging Contrast Spray has the same technological characteristics and is similar in overall design, materials and configuration compared to the predicate device. It is substantially equivalent to the predicate in terms of:

- Indications for Use
- Basic design/configuration
- Where used
- Target population
- Delivery method

**Testing and Performance Data:**

The biocompatibility of the Arges Imaging Contrast Spray materials was evaluated following the principles of ISO 10993-1.

Bench testing, including dimensional verification and functional testing was completed to support substantial equivalence of the Arges Imaging Contrast Spray. In all cases the Arges Imaging Contrast Spray performed as intended with acceptable results.

Conclusion: The Arges Imaging Contrast Spray was found to be biocompatible for its intended use and substantially equivalent in performance characteristics when compared to the predicate devices.

**Basis for Determination of Substantial Equivalence:**

Upon reviewing and comparing intended use, design, materials, principle of operation and overall technological characteristics, the Arges Imaging Contrast Spray is determined by Arges Imaging, Inc. to be substantially equivalent to existing legally marketed devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

Arges Imaging, Inc.  
C/O Mr. Mark Job  
Regulatory Technology Services, LLC  
1394 25<sup>th</sup> Street, NW  
Buffalo, MN 55313

NOV 14 2011

Re: K113263  
Trade/Device Name: Arges Imaging Contrast Spray  
Regulation Number: 21 CFR 872.3660  
Regulation Name: Impression Material  
Regulatory Class: II  
Product Code: ELW  
Dated: November 3, 2011  
Received: November 4, 2011

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 – Mr. Job

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

5. *Indications for Use Statement*

**INDICATIONS FOR USE STATEMENT**

510(k) Number (if known): K113263

Device Name: **Arges Imaging Contrast Spray**

Indications for Use:

Arges Imaging Contrast is indicated as a coating medium for optical impressions. It aids in intraoral topographical recordings of prepared and unprepared teeth and their surroundings.

OR

Prescription Use   X    
(Per 21 CFR 801 Subpart D)


Over-The-Counter Use         
(Per 21 CFR 801 Subpart C)

**(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)**

---

**Concurrence of CDRH, Office of Device Evaluation (ODE)**

Page    of   

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K113263