

K 11 3 5 71

JAN 27 2012



**510(k) Summary**

Applicant/Sponsor: Medacta International SA  
Strada Regina  
6874 Castel San Pietro (CH)  
Switzerland  
Phone (+41) 91 696 60 60  
Fax (+41) 91 696 60 66

Contact Person: Mr. Adam Gross  
Director of Regulatory and Quality  
Medacta USA  
4725 Calle Quetzal, Unit B  
Camarillo, CA, 93012  
Phone: (805)437-7085  
Fax: (805)437-7553  
Email: [AGross@medacta.us.com](mailto:AGross@medacta.us.com)

Date Prepared: December 1, 2011

**DEVICE INFORMATION**

Trade/Proprietary Name: GMK Resurfacing Patella Size 4  
Common Name: Resurfacing Patella  
Classification Name: Prosthesis, Knee, Patellofemorotibial, Semi-constrained,  
Cemented, Polymer/Metal/Polymer

21 CFR 888.3560  
Class II  
Device Product Codes: JWH

Predicate Devices: K090988 GMK Total Knee System (Medacta Intl)  
K951987 Genesis II (Smith & Nephew)  
K102437 GMK Total Knee System- Revision (Medacta Intl)  
K103170 GMK Revision SC Liners (Medacta Intl)  
K081023 Evolis Total Knee System (Medacta Intl)

### Product Description

This modification to the original Medacta GMK® (Global Medacta Knee) Total Knee System is a line extension to include the GMK Resurfacing Patella Size 4.

### Indications for Use

The GMK knee prosthesis is designed for cemented use in total knee arthroplasty, if there is evidence of sufficient sound bone to seat and support the components.

This knee replacement system is indicated in the following cases:

- Severely painful and/or disabled joint as a result of arthritis, traumatic arthritis, rheumatoid arthritis or polyarthritis.
- Avascular necrosis of femoral condyle.
- Post traumatic loss of joint configuration.
- Primary implantation failure.

Tibial augments are to be attached to the tibial baseplate with both the fixing cylinders and bone cement.

In case a semi-constrained liner is used, an extension stem must be implanted both on the tibial and on the femoral components.

### Comparison to Predicate Devices

The GMK Resurfacing Patella Size 4 has the same indications for use, material, design, and performance characteristics as the previously cleared Medacta GMK Resurfacing Patellas (Sizes 1,2,3). The GMK Resurfacing Patella Size 4 is similar to the Smith & Nephew Genesis II Resurfacing patella in terms of indications for use, material, and size.

Performance Testing

A review of the mechanical data indicates that the GMK Resurfacing Patella Size 4 are equivalent to devices currently cleared for use and are capable of withstanding expected in vivo loading without failure.

The modification to the device system to include the addition of the GMK Resurfacing Patella Size 4 was evaluated by risk analysis to identify any new risks associated with the change. Based on the risk analysis, design verification was conducted to written protocols with pre-defined acceptance criteria. The protocols and pre-defined acceptance criteria were based on the standards, FDA guidance, and comparison to the predicate device system. The testing was conducted on the worst case component size and option/design based on engineering analysis. The GMK Resurfacing Patella Size 4 was compared to the GMK Resurfacing Patella Size 1 (worst case) in terms of volume for cement insertion, contact surface area, and risk of dislocation. The GMK Resurfacing Patella Size 4 has a larger volume for cement insertion which reduces the risk of loosening of cement fixation. The GMK Resurfacing Patella Size 4 also has a larger contact surface area which reduces the contact stress between the patella and the femoral component. The GMK Resurfacing Patella Size 4 also has an adequate level of constraint to prevent dislocation because of the larger thickness and contact area as compared to the worst case.

Conclusion:

Based on indications for use, materials, design, and risk verification of volume for cement insertion, contact surface area, and risk of dislocation as compared to the worst case, the GMK Resurfacing Patella Size 4 can be considered as substantially equivalent to the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

JAN 27 2012

Medacta International  
% Medacta USA  
Mr. Adam Gross  
4725 Calle Quetzal, Unit B  
Camarillo, California, 93012

Re: K113571

Trade/Device Name: GMK Resurfacing Patella Size 4

Regulation Number: 21 CFR 888.3560

Regulation Name: Knee joint, patellofemorotibial, polymer/metal/polymer, semi-constrained cemented prosthesis.

Regulatory Class: Class II

Product Code: JWH

Dated: December 29, 2011

Received: December 30, 2011

Dear Mr. Gross:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set


forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "E. Melkerson".

 Mark N. Melkerson

Director

Division of Surgical, Orthopedic  
and Restorative Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K 113571

Device Name: GMK Resurfacing Patella, Size 4

### Indications for Use:

The GMK knee prosthesis is designed for cemented use in total knee arthroplasty, if there is evidence of sufficient sound bone to seat and support the components.

This knee replacement system is indicated in the following cases:

- Severely painful and/or disabled joint as a result of arthritis, traumatic arthritis, rheumatoid arthritis or polyarthritis.
- Avascular necrosis of femoral condyle.
- Post traumatic loss of joint configuration.
- Primary implantation failure.

Tibial augments are to be attached to the tibial baseplate with both the fixing cylinders and bone cement.

In case a semi-constrained liner is used, an extension stem must be implanted both on the tibial and on the femoral components.

Prescription Use   x    
(21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use             
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

*[Signature]*  
(Division Sign-Off)

Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number K113571

GMK Resurfacing Patella Size 4 510(k)  
November 30, 2011

Section 4 - Page 2 of 2