

K113723

FEB 15 2012



840 Memorial Drive
Cambridge, MA 02139
tel: (617) 995-5400
fax: (617) 995-5401

Section IX
Summary of Safety and Effectiveness
(Prepared on December 16, 2011)

Pursuant to §513(i)(3)(A) of the Food, Drug, and Cosmetic Act, Tepha, Inc. is submitting the following summary of information respecting safety and effectiveness:

Trade Name: TephaFLEX® Mesh

Sponsor: Tepha, Inc.
99 Hayden Avenue, Suite 360
Lexington, MA 02421

Contact Person: Mary P. LeGraw, V.P., Regulatory Affairs

Device Classification Name: CFR §878.3300
Surgical Mesh

Classification: According to Section 13 of the Federal Food, Drug and Cosmetic Act, the device classification is Class II, Performance Standards.

Predicate Devices: TephaFLEX® Surgical Mesh, TephaFLEX® Composite Mesh, SupraMesh

Device Description: The TephaFLEX mesh is a resorbable mesh prepared from poly-4-hydroxybutyrate (P4HB). The mesh is prepared from non-dyed or dyed (D&C Violet #2) monofilament P4HB fiber that is knitted into a surgical mesh. It is provided in single sheets of varying widths and lengths and may be cut to the shape or size desired for a specific application.

Indications for Use: TephaFLEX mesh is intended to reinforce soft tissue where weakness exists in patients undergoing plastic and reconstructive surgery, or for use in procedures involving soft tissue repair, such as for the repair of hernia or other fascial defects that require the addition of a reinforcing or bridging material to obtain the desired surgical result.

Safety and Performance: Mechanical testing, *in vivo* animal testing, and biocompatibility testing, were performed based on recommendations identified in the FDA surgical mesh guidance document: The Guidance for the Preparation of a Pre-market Notification Application for a Surgical Mesh. Specifically, comparative burst strength, suture pull-out strength, tensile strength and tear resistance strength was characterized. *In vivo* strength retention was characterized via a subcutaneous implantation study. The mechanical and *in vivo* data collected determined the mesh to be substantially equivalent to the predicate devices.

Conclusion: Based on the indications for use, technological characteristics, and safety and performance testing, the TephaFLEX mesh has been shown to be substantially equivalent to predicate devices under the Federal Food, Drug and Cosmetic Act.

188



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room – WO66-G609
Silver Spring, MD 20993-0002

Tepha, Inc.
% Ms. Mary P. LeGraw
V.P., Regulatory Affairs
99 Hayden Avenue, Suite 360
Lexington, Massachusetts 02421

FEB 15 2012

Re: K113723
Trade/Device Name: TephaFLEX® Mesh
Regulation Number: 21 CFR 878.3300
Regulation Name: Surgical mesh
Regulatory Class: II
Product Code: OOD
Dated: January 17, 2012
Received: January 18, 2012

Dear Ms. LeGraw:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21

CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



for

Mark N. Melkerson

Director

Division of Surgical, Orthopedic
and Restorative Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): Unknown

Device Name: TephaFLEX® Mesh

Indications for Use:

TephaFLEX mesh is intended to reinforce soft tissue where weakness exists in patients undergoing plastic and reconstructive surgery, or for use in procedures involving soft tissue repair, such as for the repair of hernia or other fascial defects that require the addition of a reinforcing or bridging material to obtain the desired surgical result.

Prescription Use: X
(21 CFR 801 Subpart D)


AND/OR

Over-The-Counter ____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1


(Division Sign-Off)
Division of Surgical, Orthopedic,
and Restorative Devices

510(k) Number K113723