



GE Healthcare

Special 510(k) Premarket Notification

GE EchoPAC Review station

January 20, 2012

510(k) Summary

In accordance with 21 CFR 807.92 the following summary of information is provided:

Date: January 20, 2012

Submitter: GE Healthcare
9900 Innovation Drive
Wauwatosa, WI, USA 53226

Primary Contact Person: Bryan Behn
Regulatory Affairs Manager
GE Healthcare
Phone: (414) 721-4214
Fax: (414) 918-8275

Secondary Contact Person: Charlotte Kaas Munthe Jørgensen
Regulatory Affairs Specialist
GE Healthcare, GE Vingmed Ultrasound AS
Phone: +47 33 02 12 80
Fax: +47 33 02 13 50

Device: Trade Name: GE EchoPAC
Common/Usual Name: Workstation Software for ultrasound image review, analysis and reporting

Classification Names:

21 CFR 892.2050

Product Code:

LLZ

Predicate Device(s):

K101324 - GE EchoPAC

Device Description: GE EchoPAC provides image processing, annotation, analysis, measurement, report generation, communication, storage and retrieval of ultrasound images that are acquired via GE Vivid family of ultrasound scanners, primarily for cardiology ultrasound applications but also for general imaging. The EchoPAC software is an integral component of each Vivid system, providing the post-acquisition image management and reporting functions of the scanner. EchoPAC will be offered as SW-only to be installed directly on customer PC hardware, or as an accessory to selected 3rd party image management workstations. EchoPAC is DICOM compliant, transferring images and data via LAN between scanners, hard copy devices, file servers and other workstations. The modified or added software features for GE EchoPAC are substantially equivalent to the predicate device and functionality cleared on GE EchoPAC K101324.

Intended Use: The GE EchoPAC workstation is indicated for diagnostic review and analysis of ultrasound images acquired under various modes of operation including B, M, Color M modes, Color, Power, Pulsed &



GE Healthcare

Special 510(k) Premarket Notification
GE EchoPAC Review station
January 20, 2012

K120221
Page 2 of 2

CW Doppler modes, Coded Pulse, Harmonic and Realtime 3D. Clinical applications include: Fetal; Abdominal; Urology (including prostate); Pediatric; Small Organ (breast, testes, thyroid); Neonatal and Adult Cephalic; Cardiac (adult and pediatric); Peripheral Vascular; Transesophageal (TEE); Musculo-skeletal Conventional; Transrectal (TR); Transvaginal (TV); and Intraoperative (abdominal, thoracic, & vascular).

Technology: The EchoPAC employs the same fundamental scientific technology as its predicate device.

Determination of Substantial Equivalence: Summary of Non-Clinical Tests:

The device has been evaluated for conformance to its design specifications and applicable industry standards for software development. It is further verified for system compatibility with the devices with which it communicates, including conformance to DICOM standard.

Summary of Clinical Tests:

The subject of this premarket submission, EchoPAC, did not require clinical studies to support substantial equivalence.

Conclusion: GE Healthcare considers the EchoPAC to be as safe, as effective, and performance is substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room -- WO66-G609
Silver Spring, MD 20993-0002

Mr. Bryan Behn
Regulatory Affairs Manager
GE Healthcare
9900 Innovation Drive, RP-2138
WAUWATOSA WI 53226

MAR 30 2012

Re: K120221
Trade/Device Name: GE EchoPAC
Regulation Number: 21 CFR 892.2050
Regulation Name: Picture archiving and communications system
Regulatory Class: II
Product Code: LLZ
Dated: March 21, 2012
Received: March 22, 2012

Dear Mr. Behn:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of

medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely Yours,

A handwritten signature in black ink, appearing to read "Janine M. Morris", is written over the typed name and title.

Janine M. Morris
Acting Director
Division of Radiological Devices
Office of In Vitro Diagnostic Device
Evaluation and Safety
Center for Devices and Radiological Health

Enclosure



GE Healthcare

Special 510(k) Premarket Notification
GE EchoPAC Review station
January 20, 2012

510(k) Number: K120221

Device Name: GE EchoPAC

Indications for Use:

The GE EchoPAC workstation is indicated for diagnostic review and analysis of ultrasound images acquired under various modes of operation including B, M, Color M modes, Color, Power, Pulsed & CW Doppler modes, Coded Pulse, Harmonic and Realtime 3D. Clinical applications include: Fetal; Abdominal; Urology (including prostate); Pediatric; Small Organ (breast, testes, thyroid); Neonatal and Adult Cephalic; Cardiac (adult and pediatric); Peripheral Vascular; Transesophageal (TEE); Musculoskeletal Conventional; Transrectal (TR); Transvaginal (TV); and Intraoperative (abdominal, thoracic, & vascular).


Prescription Use: YES
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use: NO
(Part 21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)


(Division Sign-Off)
Division of Radiological Devices
Office of In Vitro Diagnostic Device Evaluation and Safety
K120221
510K