



PRODUCT: MILAGRO INTERFERENCE SCREWS  
SUBMISSION DATE: FEBRUARY 27<sup>TH</sup>, 2012  
SUBMISSION TYPE: TRADITIONAL

## ATTACHMENT 1

### 510(k) SUMMARY - DEPUY MITEK MILAGRO® INTERFERENCE SCREWS

#### SUBMITTER'S NAME AND ADDRESS

DePuy Mitek, Inc.  
a Johnson & Johnson company  
325 Paramount Drive  
Raynham, MA 02767

#### CONTACT PERSON

Deep Pal  
Senior Regulatory Affairs Specialist

TELEPHONE 508-977-3998  
FACSIMILE 508-977-6911  
E-MAIL [dpal3@its.jnj.com](mailto:dpal3@its.jnj.com)  
DATE PREPARED 02/17/2012

#### NAME OF MEDICAL DEVICE

##### CLASSIFICATION NAME

Fastener, Fixation, Biodegradable, Soft Tissue

##### COMMON/USUAL NAME

Bone Anchor

##### PROPRIETARY NAME

DePuy Mitek Milagro® Interference Screws

#### SUBSTANTIAL EQUIVALENCE

The proposed DePuy Mitek Milagro® Interference Screws are substantially equivalent to the following devices.


- K103831, K060830, K032717: Milagro Interference Screws
- K051726, K041356, K020043: Arthrex Tenodesis Screws
- K102410: Linvatec's Matryx Screws

#### FDA PRODUCT CODE

MAI, HWC

#### DEVICE CLASSIFICATION

This type of fixation screw was originally classified as a Class II medical device by the Orthopedic Review Panel, regulated as 21 CFR 888.3040 Smooth or Threaded Metallic Bone Fixation Fastener.



PRODUCT: MILAGRO INTERFERENCE SCREWS  
 SUBMISSION DATE: FEBRUARY 27<sup>TH</sup>, 2012  
 SUBMISSION TYPE: TRADITIONAL

## ATTACHMENT 1

*Continues...*

### 510(k) SUMMARY - DEPUY MITEK MILAGRO® INTERFERENCE SCREWS

#### DEVICE DESCRIPTION

The DePuy Mitek Milagro® Interference Screws are absorbable, tapered, cannulated, threaded fasteners for use in interference fixation of soft tissue grafts or bone-tendon grafts. The Interference Screw is made from a composite made of absorbable Poly (lactide-co-glycolide) polymer and Tricalcium Phosphate (TCP).

#### INDICATIONS FOR USE

The Milagro® BR interference screws are designed to attach soft tissues to bone in orthopedic surgical procedures. The screws may be used for interference fixation of soft tissues (such as ligaments or tendons) to bone, when the implant sizes offered are patient appropriate. The implant operates, in conjunction with the appropriate postoperative immobilization, throughout the healing period.

#### TECHNOLOGICAL CHARACTERISTICS

The design specifications of the proposed DePuy Mitek Milagro® Interference Screws are substantially equivalent to the existing DePuy Mitek Milagro® Interference Screws cleared under 510(k) K103831, K060830 and K032717. Technological characteristics including design construct, packaging and indications are similar to the predicate devices and use similar or identical material and packaging as the predicates.

#### NONCLINICAL TESTING

Product Design Verification activities, such as, Insertion Torque, Anchor Pullout (at T=0, 3, 6 and 12 week in-vitro physiological aging), and Torque to Failure were performed on the implant.

#### SAFETY AND PERFORMANCE

Results of performance and safety testing have demonstrated that the proposed device is substantially equivalent to the predicate devices.

Based on the indications for use, technological characteristics, and comparison to predicate devices, the proposed Small Size DePuy Mitek Milagro® Interference Screws have been shown to be substantially equivalent to predicate devices under the Federal Food, Drug and Cosmetic Act.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

DePuy Mitek, Incorporated  
a Johnson & Johnson Company  
% Mr. Deep Pal  
Senior Regulatory Affairs Specialist  
325 Paramount Drive  
Raynham, Massachusetts 02767

APR 24 2012

Re: K120589

Trade/Device Name: Milagro® Interference Screws  
Regulation Number: 21 CFR 888.3040  
Regulation Name: Smooth or threaded metallic bone fixation fastener  
Regulatory Class: Class II  
Product Code: MAI, HWC  
Dated: February 27, 2012  
Received: February 28, 2012

Dear Mr. Pal:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.


You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Mark N. Melkerson', is written over the typed name and title.

Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure



PRODUCT: MILAGRO INTERFERENCE SCREWS  
SUBMISSION DATE: FEBRUARY 27<sup>TH</sup>, 2012  
SUBMISSION TYPE: TRADITIONAL

## ATTACHMENT 2

### INDICATIONS FOR USE

510(k) Number (if known): \_\_\_\_\_


Device Names: Milagro® Interference Screws

**Indications for Use:** The Milagro® BR interference screws are designed to attach soft tissues to bone in orthopedic surgical procedures. The screws may be used for interference fixation of soft tissues (such as ligaments or tendons) to bone, when the implant sizes offered are patient appropriate. The implant operates, in conjunction with the appropriate postoperative immobilization, throughout the healing period.

Prescription Use ☒ AND/OR Over-The-Counter Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)  
(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE).

Page 1 of \_\_\_\_\_

  
(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number K120589