

Page 1 of 2

**6. 510(k) Summary**

This 510(k) summary information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

**APPLICANT:** Pinnacle Spine Group, LLC  
**DATE PREPARED:** May 16, 2012  
**CONTACT PERSON:** Rebecca K Pine  
1601 Elm Street, Suite 300  
Dallas, TX 75201  
Phone: 760.809.5178  
Fax: 760.290.3216  
**TRADE NAME:** InFill™ Graft Delivery System  
**COMMON NAME:** Piston Syringe  
**CLASSIFICATION NAME:** Intervertebral Body Fusion Device  
**DEVICE CLASSIFICATION:** Class 2, per 21 CFR 880.5860  
**PRODUCT CODE:** FMF  
**PREDICATE DEVICES:** InFill™ graft Delivery System (K111632)

**Substantially Equivalent To:**

The modified InFill™ Graft Delivery System is substantially equivalent in intended use, principal of operation and technological characteristics to the InFill™ Graft Delivery System cleared under premarket notification K111632.

**Description of the Device Subject to Premarket Notification:**

The modified InFill™ Graft Delivery System is comprised of a disposable medical piston syringe, a cannulated applicator tip and accessories for mixing bone graft materials and filling of the system. The modified InFill™ Graft Delivery System is provided sterile, for single use only.

**Indication for Use:**

The InFill™ Graft Delivery System is intended to be used for the delivery of hydrated allograft, autograft or synthetic bone graft material to an orthopedic surgical site.

**Technical Characteristics:**

The modified InFill™ Graft Delivery System has the same technological characteristics and is similar in overall design, materials and configuration compared to the current InFill™ Graft Delivery System.

**Performance Data:**

All necessary testing has been performed for the InFill™ Graft Delivery System to assure substantial equivalence to the predicate device and demonstrate the device performs as intended. All testing was performed on test units representative of finished devices.

The device design was qualified through the following tests:

- Simulated Use Testing
- Volume Verification
- Separation Force Testing
- Liquid Leak Testing
- Biocompatibility Testing

The modified InFill™ Graft Delivery System met all specified criteria and did not raise new safety or performance questions.

**Basis for Determination of Substantial Equivalence:**

The Indication/Intended Use and the fundamental scientific technology of the modified device have not been changed and are the same as those described in the unmodified predicate device. The modified InFill™ Graft Delivery System is determined by Pinnacle Spine Group LLC, to be substantially equivalent to the InFill™ Graft Delivery System (K111632).



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

AUG 29 2012

Pinnacle Spine Group, LLC  
% Ms. Rebecca K. Pine  
Consultant  
1601 Elm Street, Suite 300  
Dallas, Texas 75201

Re: K121476  
Trade/Device Name: InFill Graft Delivery System  
Regulation Number: 21 CFR 880.5860  
Regulation Name: Piston syringe  
Regulatory Class: Class II  
Product Code: FMF  
Dated: August 2, 2012  
Received: August 3, 2012

Dear Ms. Rebecca K. Pine:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21


Page 2 – Ms. Pine

CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "for DTH" followed by a stylized signature, and "Det. Cur. Dir." written below it.

Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## 5. Indications for Use Statement

## INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K121476Device Name: **InFill™ Graft Delivery System**

Indications for Use:


The InFill™ Graft Delivery System is intended to be used for the delivery of hydrated allograft, autograft or synthetic bone graft material to an orthopedic surgical site.

AND/OR

Prescription Use X  
(Part 21 CFR 801 Subpart D)Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)  
Division Sign-Off)Division of Surgical, Orthopedic,  
and Restorative DevicesPage 1 of 1510(k) Number K121476