

K122727

Abbreviated 510(k) Premarket Notification
Sterile Nitrile Powder Free Examination Gloves

Section 5, 510(k) Summary
Page 1 of 3

5. 510(k) SUMMARY

DATE: September 20, 2012

OWNER: Northstar Healthcare Holdings
70 Sir John Rogerson's Quay
Dublin 2, Ireland

FEB 22 2013

OFFICIAL CORRESPONDENT: Michael Riordan
Operations Manager
Telephone: 00353-21-4548255
Fax: 00353-21-4548294
Email: michael.riordan@mckesson.ie

DEVICE NAME:

Trade Name: Textured Blue Sterile Powder Free Nitrile Examination Gloves

Common Name: Patient Examination Gloves

Classification: Patient Examination Gloves

Class: Class I

Product Code: LZA

Regulation Number: 880.6250

PREDICATE DEVICE(S):

Predicate 510(k)	Device Name	Indication	Clearance Date	Company
K081466	Sterile Nitrile Powder Free Examination Gloves	The examination glove is a disposable device intended for medical and dental purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.	30 Jul 2008	Smart Glove Corp., SDN, BHD

DEVICE DESCRIPTION: Textured Blue Sterile Powder Free Nitrile Examination Gloves

**STATEMENT OF
INTENDED USE:**

The Textured Blue Sterile Powder Free Nitrile Examination Gloves is a Disposable device intended for medical and dental purposes to be worn on the examiner's hand or finger to prevent contamination between patient and examiner.

**TECHNOLOGICAL
CHARACTERISTICS:**

The Textured Blue Sterile Powder Free Nitrile Examination Gloves is substantially equivalent to the predicate device with regard to physical characteristics, design, product features, and intended use. Both gloves are made with nitrile using similar manufacturing processes.

Feature	Sterile Nitrile Powder Free Examination Gloves K081466 Predicate	Textured Blue Sterile Powder Free Nitrile Examination Gloves Proposed																								
Intended Use	Intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.	Same																								
Indications for Use Statement	The examination glove is a disposable device intended for medical and dental purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.	The examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.																								
Description	Sterile, powder free, examination gloves made of nitrile. The gloves are provided in sizes extra small, small, medium, large and extra large.	Sterile, powder free, examination gloves made of nitrile. The gloves are provided in sizes small, medium, and large.																								
Presentation	Sterile gloves are provided in dispenser boxes.	Same																								
Material	Nitrile	Same																								
Sterilization	Sterile	Same																								
Single Use	Yes	Same																								
Dimensions	Meets ASTM D6319-00a(2005)	<table> <tr> <td>Length</td><td>Small</td><td>230mm min.</td></tr> <tr> <td></td><td>Medium</td><td>230mm min.</td></tr> <tr> <td></td><td>Large</td><td>230mm min.</td></tr> <tr> <td>Width</td><td>Small</td><td>70-90mm</td></tr> <tr> <td></td><td>Medium</td><td>85-105mm</td></tr> <tr> <td></td><td>Large</td><td>101-120mm</td></tr> <tr> <td>Thickness</td><td>Finger</td><td>0.08mm min.</td></tr> <tr> <td></td><td>Palm</td><td>0.06mm min.</td></tr> </table>	Length	Small	230mm min.		Medium	230mm min.		Large	230mm min.	Width	Small	70-90mm		Medium	85-105mm		Large	101-120mm	Thickness	Finger	0.08mm min.		Palm	0.06mm min.
Length	Small	230mm min.																								
	Medium	230mm min.																								
	Large	230mm min.																								
Width	Small	70-90mm																								
	Medium	85-105mm																								
	Large	101-120mm																								
Thickness	Finger	0.08mm min.																								
	Palm	0.06mm min.																								
Physical Properties	Meets ASTM D6319-00a(2005)	<table> <tr> <td></td><td colspan="2">Before aging/after aging</td></tr> <tr> <td>Elongation</td><td>500%</td><td>400%</td></tr> <tr> <td>Tensile Strength</td><td>14MPa</td><td>14MPa</td></tr> </table>		Before aging/after aging		Elongation	500%	400%	Tensile Strength	14MPa	14MPa															
	Before aging/after aging																									
Elongation	500%	400%																								
Tensile Strength	14MPa	14MPa																								
Freedom from Pinholes	Meets ASTM D5151-06	Same																								
Residual Powder	Meets ASTM D6124-06	Same																								

**ASSESSMENT OF
NONCLINICAL DATA:**

Characteristic	Standard	Device Performance
Dimension	ASTM Standard D6319-10	Meets
Physical Properties	ASTM Standard D6319-10	Meets
Freedom from Pinholes	21 CFR 800.20; ASTM D5151-06	Meets
Powder Residual	ASTM Standard D6124-06	Meets Results generated values below 2mg of residual powder
Biocompatibility	Primary Skin Irritation in rabbits (ISO 10993-10:2010)	Gloves are non-irritating
	Dermal Sensitization in the guinea pig (ISO 10993-10:2010)	Gloves do not display any potential for sensitization

CONCLUSIONS:

The Textured Blue Sterile Powder Free Nitrile Examination Gloves meet the requirements of established standards ASTM D6319-10, ASTM D5151-06, ASTM D6124-06 and ISO 10993-10:2010.

Based on the comparison of intended use, design, materials and performance, the Textured Blue Sterile Powder Free Nitrile Examination Gloves are substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

February 22, 2013

Northstar Healthcare Holdings
C/O Mr. Ned Devine
Responsible Third Party Official
Underwriters Laboratories, Incorporated
333 Pfingsten Road
NORTHBROOK IL 60062

Re: K122727

Trade/Device Name: Textured Blue Sterile Powder Free Nitrile Examination Gloves
Regulation Number: 21 CFR 880.6250
Regulation Name: Patient Examination Glove
Regulatory Class: I
Product Code: LZA
Dated: February 7, 2013
Received: February 8, 2013

Dear Mr. Devine:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Anthony D. Watson", written over a stylized graphic that resembles a medical cross or a set of scales.

Anthony D. Watson, B.S., M.S., M.B.A.
Director
Division of Anesthesiology, General Hospital,
Respiratory, Infection Control and
Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATIONS FOR USE

510(k) Number (if known): K122727

Device Name: **Textured Blue Sterile Powder Free Nitrile Examination Gloves**

Indications for Use: The examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

Prescription Use _____
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use X
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE
OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Elizabeth F. Claverie
2013.02.20 18:12:38 -05'00'

(Division Sign-Off)
Division of Anesthesiology, General Hospital
Infection Control, Dental Devices

510(k) Number: K122727

Page 1 of 1