

MAR 28 2013

**510(k) SUMMARY**

A summary of 510(k) safety and effectiveness information in accordance with the requirements of 21 CFR 807.92

<b>SUBMITTER INFORMATION</b>	
<b>Name</b>	Biomet Manufacturing Corp.
<b>Address</b>	56 East Bell Drive Warsaw, IN 46582
<b>Phone number</b>	(574) 267-6639
<b>Fax number</b>	(574) 371-1027
<b>Establishment Registration Number</b>	1825034
<b>Name of contact person</b>	Patricia Sandborn Beres Senior Regulatory Specialist Biomet Manufacturing Corp.
<b>Date prepared</b>	March 21, 2013
<b>NAME OF DEVICE</b>	
<b>Trade name</b>	Biomet Reconstructive Wedges
<b>Common name</b>	Wedge
<b>Classification name</b>	Single/multiple component metallic bone fixation appliances and accessories Smooth or threaded metallic bone fixation fastener
<b>Classification panel</b>	Orthopedics
<b>Regulation</b>	21 CFR 888.3030 21 CFR 888.3040
<b>Product Code(s)</b>	HRS HWC
<b>Legally marketed device(s) to which equivalence is claimed</b>	Wright Medical Technology Inc.'s BIOFOAM® Bone Wedges cleared through the following 510(k)s: K070592, K073535 and K093950.
<b>Reason for 510(k) submission</b>	New device
<b>Device description</b>	Biomet Reconstructive Wedges are metallic blocks designed for insertion into voids for alignment of bones in the foot or ankle. There are two styles of wedges which are oval and rectangular in shape. Wedges are available in a variety of widths, lengths and thicknesses. Rectangular cut-out windows allow the surgeon to pack allograft within the device. The wedges are completely porous, to allow for tissue ingrowth. There is a "skin" of solid titanium on the outer edge of the component to facilitate the insertion instrument and allow for identifying part marking to be added to the devices.

Mailing Address:  
P.O. Box 587  
Warsaw, IN 46581-0587  
Tel: (574) 267-6639  
Fax: (574) 371-1027  
www.biomet.com

Shipping Address:  
56 East Bell Drive  
Warsaw, IN 46582

<b>Indications for use</b>	<p>Biomet Reconstructive Wedges are intended to be used for internal bone fixation for bone structures, fusions or osteotomies in the ankle and foot, such as:</p> <ul style="list-style-type: none"> <li>• Open wedge osteotomies of Hallux Valgus</li> <li>• Evans lengthening osteotomies</li> <li>• Metatarsal/cuneiform arthrodesis</li> </ul> <p>Biomet Reconstructive Wedges are intended for use with ancillary fixation.</p> <p>Biomet Reconstructive Wedges are not intended for use in the spine.</p>
<b>SUMMARY OF THE TECHNOLOGICAL CHARACTERISTICS COMPARED TO THE PREDICATE</b>	
<p>Biomet Reconstructive Wedges uses similar technology to the predicate BIOFOAM® Bone Wedge Devices cleared through 510(k)s K070592, K073535, and K093950. The Biomet Reconstructive Wedges are similar in material, design and sizing to the predicates.</p>	
<b>PERFORMANCE DATA</b>	
<b>Non-Clinical Tests Conducted For Determination Of Substantial Equivalence</b>	
Chemical Composition	Static Shear Strength
Static Compression Testing	Dynamic Compression Testing
Expulsion Testing	Compressive Strength
Porosity/Pore Size	Taber Abrasion
Interconnectivity	Roughness
Tensile Adhesion	MR Compatibility
Microstructure	Animal Data
<b>Clinical Tests Conducted for Determination of Substantial Equivalence and/or of Clinical Information</b>	
No clinical data submitted	
<b>CONCLUSIONS DRAWN FROM NON-CLINICAL AND CLINICAL DATA</b>	
<p>No clinical data was necessary for a determination of substantial equivalence.</p> <p>The results of testing indicated the material performed within the intended use, did not raise any new safety and efficacy issues and were found to be substantially equivalent to the predicate devices.</p>	



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-002

March 28, 2013

Biomet Manufacturing Corporation  
% Ms. Patricia Sandborn Beres  
Senior Regulatory Specialist  
56 East Bell Drive  
Warsaw, Indiana 46582

Re: K122770

Trade/Device Name: Biomet Reconstructive Wedges  
Regulation Number: 21 CFR 888.3030  
Regulation Name: Single/multiple component metallic bone fixation appliances and accessories  
Regulatory Class: Class II  
Product Code: HRS, HWC  
Dated: February 18, 2013  
Received: February 20, 2013

Dear Ms. Beres:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucml15809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**ErinFDKeith**

Mark N. Melkerson  
Director  
Division of Orthopedic Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K122770

Device Name: Biomet Reconstructive Wedges

### Indications For Use:

Biomet Reconstructive Wedges are intended to be used for internal bone fixation for bone structures, fusions or osteotomies in the ankle and foot, such as:

- Open wedge osteotomies of Hallux Valgus
- Evans lengthening osteotomies
- Metatarsal/cuneiform arthrodesis

Biomet Reconstructive Wedges are intended for use with ancillary fixation.

Biomet Reconstructive Wedges are not intended for use in the spine.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use   NO    
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Elizabeth Frank -S

Division of Orthopedic Devices

Page 1 of 1