

K 123531

## 510 (k) Summary

### Submitter's information:

DEC 14 2012

**Name:** LeMaitre Vascular, Inc.  
**Address:** 63 Second Avenue  
Burlington, MA USA 01803  
**Phone:** 781-425-1727  
**Contact Person:** Bryan Cowell, MSc., RAC

**Date of preparation:** 14 Nov 2012  
**Device Name:** UnBalloon Non-Occlusive Modeling Catheter  
**Trade Name:** UnBalloon Non-Occlusive Modeling Catheter

**Common/ Classification Name:** Catheter, Percutaneous / Modeling Catheter

**Classification Panel:** 21CFR §870.1250  
**Class:** II (2)

**Product Code:** DQY

**Establishment Registration:** 1220948

**Establishment:** LeMaitre Vascular, Inc., 63 Second Avenue, Burlington, MA USA 01803

**Owner/Operator:** 1220948

### Proposed Device Description:

The UnBalloon Non-Occlusive Catheter is a silicone surface coated (medical grade) modeling catheter with an expandable Nitinol mesh in a 14F retractable sheath. The Nitinol mesh design allows for expansion without occluding blood flow. The Nitinol mesh and radiopaque markers are highly visible under fluoroscopy and assist in the positioning of the device. The inner lumen allows for a 0.035 or 0.038 inch guidewire for over-the-wire access. Side ports and clear handle/luer allow the device and guidewire lumen to be flushed. The blue handle allows the device to be sheathed/unsheathed while the clear handle/luer controls the expansion of the Nitinol mesh. This submission modifies the current UnBalloon Non-Occlusive Modeling Catheter by expanding the product family line with additional nitinol mesh cage sizes and catheter lengths.

**Proposed Intended Use:**

The UnBalloon Non-Occlusive Modeling Catheter is intended to assist in the modeling of self-expanding endoprotheses in large diameter vessels.

**Predicate Device:**

510(k): K121839  
Device Name: UnBalloon Non-Occlusive Modeling Catheter  
SE Date: 08/08/2012  
Regulation Number: 870.1250  
Device Class Name: Catheter, Percutaneous  
Device Class: 2

510(k): K110891  
Device Name: UnBalloon Non-Occlusive Modeling Catheter  
SE Date: 09/13/2011  
Regulation Number: 870.1250  
Device Class Name: Catheter, Percutaneous  
Device Class: 2

**Substantial Equivalence:****Fundamental Scientific Technological Characteristics:**

The UnBalloon Non-Occlusive Modeling Catheter is a silicone surface coated percutaneous/modeling catheter designed for vascular surgeons. This submission is to include additional product sizes.

**Functional/ Safety testing:**

The verification activities conducted indicate that UnBalloon Non-Occlusive Modeling Catheter device meets the product performance requirements of the device specifications and does not raise any additional safety issues.

**Sterilization:**

The device is validated for ethylene oxide (EO) sterilization according to ANSI/AAMI/ISO 11135-1:2007, "Sterilization of Medical Devices – Validation and Routine Control of Ethylene Oxide Sterilization"

**Biocompatibility:**

All blood contact portions of the device were subjected to Biocompatibility testing according to ISO 10993 guidelines for an externally communicating device with limited contact duration (<24 hours), in circulating blood. There is no new material used for the additional catheter sizes as discussed in this submission.

**Summary of Product Testing:**

The following tests have been completed to evaluate the performance of the catheter with shorter cage:

*Radial outward force*

*Apposition Length*

*Simulated Anatomical Use*



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-002

LeMaitre Vascular, Inc.  
c/o Mr. Bryan Cowell  
Principal Regulatory Affairs Specialist  
63 Second Avenue  
Burlington, MA 01803

DEC 14 2012

Re: K123531

Trade/Device Name: UnBalloon Non-Occlusive Modeling Catheter

Regulation Number: 21 CFR 870.1250

Regulation Name: Percutaneous catheter

Regulatory Class: Class II (two)

Product Code: DQY

Dated: November 15, 2012

Received: November 16, 2012

Dear Mr. Cowell:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use Statement

510(k)  
Number  
(if known)

K123531

Device Name

UnBalloon Non-Occlusive Modeling Catheter

Indications  
for Use

The UnBalloon Non-Occlusive Modeling Catheter is intended to assist in the modeling of self-expanding endoprostheses in large diameter vessels.


PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use   X    
(Per 21 CFR 801.109)

OR

Over-The-Counter Use           

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Cardiovascular Devices  
510(k) Number   K123531