

ATTACHMENT E

MAY 31 2013

510(k) SUMMARY

Trade Name: BRM ARS AND SLANT SCREWS

Sponsor: Advanced Interventional Technology LLC
6703 NW 167 Street, Suite C-9
Miami, Fla 33015

Contact: E. March
GAMA Associates LLC
Ph- 240.506.3212
Email- edogama@comcast.net

Device Generic Name: Bone Fixation Screw

Classification: CFR 888.3040, Class II

Product Code: HWC

Product Description:

The ARS and SLANT Screws are made of titanium alloy (Ti6Al4V) according to ISO 5832-3 and ASTM F 136. They are self-drilling, self-tapping and designed with a double thread for better grip on tissue. The ARS and SLANT Screws are intended for single use only.

The titanium ARS and SLANT SCREWS are cannulated for use with a guide wire and have double-threading and different pitch for the proximal and distal threads. The headless screws can be completely embedded in the bone thus avoiding any protrusion that could cause irritation to the tissue. The double-thread design assures compression and stability of the fractures. Both screws have reverse-cutting flutes in the head threaded portion to facilitate removal of the screw.

The code to clearly identify each product is laser-marked on the unthreaded cylindrical screw part. The screws are color coded (anodized) to be easily recognized ARS (blue) from SLANT (green).

Date: 29 April 2013

Indications for Use:

The ARS and SLANT SCREW medical devices are indicated for the fixation of bone fractures and bone reconstruction. Examples include but are not limited to:

- Fixation of bone fragments in long or small bone fractures
- arthrodesis in hand or foot surgery
- mono or bi-cortical osteotomies in the foot or hand
- distal or proximal metatarsal or metacarpal osteotomies
- fixation of osteotomies for hallux valgus treatment (such as Scarf, Chevròn, etc.)

Predicate Devices

- NewDeal QWIX Screws (K071639, K050346)

Substantial Equivalence Information

Based on available 510k information including engineering analyses, AIT bone fixation screws (BRM ARS and Slant) are deemed substantially equivalent to the predicate devices in terms of indications for use, material, technology and design specifications.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

Advanced Interventional Technology LLC
% GAMA Associates LLC
Mr. Eduardo March
7000 Cashell Manor Court
Derwood, Maryland 20855

Letter dated: May 31, 2013

Re: K130954

Trade/Device Name: BRM ARS and SLANT Screws
Regulation Number: 21 CFR 888.3040
Regulation Name: Smooth or threaded metallic bone fixation fastener
Regulatory Class: Class II
Product Code: HWC
Dated: May 2, 2013
Received: May 3, 2013

Dear Mr. March:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

~~You may, therefore, market the device, subject to the general controls provisions of the Act.~~

~~The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.~~

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Mark N. Melkerson -S

Mark N. Melkerson
Director
Division of Orthopedic Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Page 1 of 1

510(k) Number (if known): K130954

Device Name: BRM ARS and SLANT SCREWS

Indications for Use:

The ARS and SLANT SCREW medical devices are indicated for the fixation of bone fractures and bone reconstruction. Examples include but are not limited to:

- Fixation of bone fragments in long or small bone fractures
- arthrodesis in hand or foot surgery
- mono or bi-cortical osteotomies in the foot or hand
- distal or proximal metatarsal or metacarpal osteotomies
- fixation of osteotomies for hallux valgus treatment (such as Scarf, Chevron, etc.)

Prescription Use X AND / OR
(PART 21 CFR 801.Subpart D)

Over-the -Counter Use _____
(PART 21 CFR 801.Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, OFFICE OF DEVICE EVALUATION

CONFIDENTIAL

29 April 2013

Elizabeth L. Frank -S

Division of Orthopedic Devices