

JUL 03 2013

**510(k) Summary**

---

**NAME OF FIRM:** Emerge Medical  
790 S. Colorado Blvd.  
Suite 550-S  
Denver, CO 80246

**DATE PREPARED:** May 20, 2013

**510(K) CONTACT:** Victoria Trafka  
Vice President of Engineering & Quality  
Tel: (303) 225-7909  
Fax: (800) 698-1440

**PROPOSED TRADE NAME:** Emerge Medical Locking One-third Tubular Plate System  
**DEVICE CLASSIFICATION:** Class II; 21 CFR 888.3030 and 888.3040

**CLASSIFICATION NAME:** Plate, Fixation, Bone; Smooth or threaded metallic bone fixation fastener

**PRODUCT CODE:** HRS and HWC

**DEVICE DESCRIPTION:** The System consists of stainless steel components including locking plates and standard cortex screws, locking cortex screws and standard cancellous screws and washers. The plates are available in a variety of lengths with the number of holes varying depending on plate length. The screws and plates are provided non-sterile.

**INDICATIONS FOR USE:** The Emerge Medical Locking One-third Tubular Plate System is intended for fixation of fractures, osteotomies, and non-unions of the clavicle, scapula, olecranon, humerus, radius, ulna, pelvis, distal tibia, and fibula, including osteopenic bone.

**MATERIALS:** Stainless steel (ASTM F138 and F139)

**PREDICATE DEVICES:** Synthes K011335

**TECHNOLOGIC CHARACTERISTICS:** The fundamental scientific principles and technological characteristics, including the intended use, material, general design, and sizes of the device are equivalent to the predicate device.

**PERFORMANCE DATA:** Mechanical testing performed according to ASTM F543 and ASTM F382 demonstrated that the device performs as well as or better than the predicate device. Clinical data were not needed to demonstrate substantial equivalence.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

July 3, 2013

Emerge Medical  
% Ms. Victoria Trafka  
Vice President of Engineering and Quality  
790 South Colorado Road, Suite 550-S  
Denver, Colorado 80246

Re: K131480

Trade/Device Name: Emmerge Medical Locking One-Third Tubular Plate System  
Regulation Number: 21 CFR 888.3030  
Regulation Name: Single/multiple component metallic bone fixation appliances and accessories  
Regulatory Class: Class II  
Product Code: HRS, HWC  
Dated: May 31, 2013  
Received: June 13, 2013

Dear Ms. Trafka:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Erin Keith**

For

Mark N. Melkerson  
Director  
Division of Orthopedic Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use Statement**

---

510(k) Number (if known): K131480

Device Name: **Emerge Medical Locking One-Third Tubular Plate System**

Indications for Use:

The EmERGE Locking One-third Tubular Plate System is intended for fixation of fractures, osteotomies, and non-unions of the clavicle, scapula, olecranon, humerus, radius, ulna, pelvis, distal tibia, and fibula, including osteopenic bone.

Prescription Use   X   or Over-The-Counter Use             
(21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Elizabeth L. Frank -S**

Division of Orthopedic Devices