



never stop moving

**PRODUCT:** HEALIX ADVANCE KNOTLESS  
PEEK ANCHOR (6.5mm)  
**SUBMISSION DATE:** June 6, 2013  
**SUBMISSION TYPE:** SPECIAL

## 510(k) SUMMARY

JUN 27 2013

**Submitter:**

DePuy Mitek  
a Johnson & Johnson company  
325 Paramount Drive  
Raynham, MA 02767

**Contact Person**

Kristine Christo  
Manager, Regulatory Affairs  
DePuy Mitek  
a Johnson & Johnson company  
325 Paramount Drive  
Raynham, MA 02767, USA

Telephone: 508-828-3359  
Facsimile: 508-977-6911  
e-mail: [kchristo@its.jnj.com](mailto:kchristo@its.jnj.com)

**Date Prepared**

June 6, 2013

**Name of Medical Device**

**Proprietary Name:** HEALIX ADVANCE KNOTLESS PEEK ANCHOR (6.5mm)  
**Classification Name:** Fastener, Fixation, Biodegradable, Soft tissue  
**Common Name:** Bone Anchor

**Substantial Equivalence**

The HEALIX ADVANCE KNOTLESS PEEK ANCHOR (6.5mm) is substantially equivalent to:  
▪ K130539 Mitek Healix Advance Knotless PEEK Anchor (4.75 and 5.5mm)

**Device Classification**

Fastener, Fixation, Nondegradable, Soft Tissue, classified as Class II, product code MBI regulated under 21 CFR 888.3040.

**Device Description**

The proposed Healix Advance Knotless Anchor is a one piece implantable cannulated, threaded anchor designed to secure soft tissue to bone. The anchor is provided loaded on a disposable inserter driver device. The proposed anchors will be offered in a 6.5 mm size. The proposed 6.5 mm Healix Advance Knotless PEEK Anchor is manufactured from PEEK (Polyetheretherketone) material.

**Indications for Use**

The Healix Advance Knotless Anchors are indicated for use in the following procedures for reattachment of soft tissue to bone:

## Shoulder

- Rotator Cuff
- Biceps Tenodesis



never stop moving™

**PRODUCT:** HEALIX ADVANCE KNOTLESS  
PEEK ANCHOR (6.5mm)  
**SUBMISSION DATE:** June 6, 2013  
**SUBMISSION TYPE:** SPECIAL

---

**Comparison of  
Technological  
Characteristics**

The proposed Healix Advance Knotless PEEK Anchors will have the same design as compared to the predicate devices (4.75mm and 5.5mm) but will be larger in size (6.5mm). Both the proposed and predicate Healix Advance Knotless PEEK Anchors are molded from the same PEEK (polyetheretherketone) material. No new technological characteristics were introduced as a result of the proposed changes.

---

**Safety and  
Performance****Non-clinical Testing**

Product Design Verification and Design Validation activities, such as, Insertion Torque, Torque to Failure and Anchor Pullout were performed on the proposed implant device. Results of performance and safety testing have demonstrated that the proposed device is substantially equivalent to the predicate devices.

Based on the indications for use, technological characteristics, and comparison to predicate devices, the proposed 6.5mm Healix Advance Knotless Anchors have been shown to be substantially equivalent to predicate devices under the Federal Food, Drug and Cosmetic Act.

---



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

June 27, 2013

Depuy Mitek Incorporated, a Johnson & Johnson Company  
% Ms. Kristine Christo  
Manager, Regulatory Affairs  
325 Paramount Drive  
Raynham, Massachusetts 02767

Re: K131683

Trade/Device Name: Healix Advance Knotless PEEK Anchor  
Regulation Number: 21 CFR 888.3040  
Regulation Name: Smooth or threaded metallic bone fixation fastener  
Regulatory Class: Class II  
Product Code: MBI  
Dated: June 7, 2013  
Received: June 10, 2013

Dear Ms. Christo:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

For   
Mark N. Melkerson  
Director  
Division of Orthopedic Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K131683

---

Device Names: Healix Advance Knotless PEEK Anchor

**Indications for Use:** The Healix Advance Knotless Anchors are indicated for use in the following procedures for reattachment of soft tissue to bone:

Shoulder

- Rotator Cuff
- Biceps Tenodesis

Prescription Use   √  

AND/OR

Over-The-Counter Use                     

(Part 21 CFR 801 Subpart D)

(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of   1  

Casey L. Hanley, Ph.D.

Division of Orthopedic Devices