

**510(K) SUMMARY**

**A. Submitter Information**

**Manufacturer:** Medos International Sàrl  
Chemin-Blanc 38  
2400 Le Locle, Switzerland

**Submitter:** DePuy Spine, Inc.  
325 Paramount Drive  
Raynham, MA 02767

**Contact Person:** Kirsten Lehmuller  
325 Paramount Drive  
Raynham, MA 02767  
**Telephone number:** 508-828-3291  
**Fax number:** 508-828-3797  
**Email:** klehmull@its.jnj.com

SEP 06 2013

**B. Date Prepared** July 25, 2013

**C. Device Name**

**Trade/Proprietary Name:** SKYLINE® Anterior Cervical Plate System,  
UNIPLATE® Anterior Cervical Plate System,  
and UNIPLATE®2 Anterior Cervical Plate  
System

**Common/Usual Name:** Spinal System

**Classification Name:** Spinal intervertebral body fixation orthosis  
per 21 CFR §888.3060

**D. Predicate Device Name**

**Trade name:** SKYLINE® Anterior Cervical Plate System (K052552)  
UNIPLATE® Anterior Cervical Plate System (K042544)  
UNIPLATE®2 Anterior Cervical Plate System (K082273, K100070)

**E. Device Description**

The UNIPLATE® and UNIPLATE® 2 Anterior Cervical Plate Systems consist of an assortment of titanium alloy plates and screws. The anterior cervical plates are available in various lengths to accommodate one to three segments of fixation. The screws are available in various sizes and tip geometries.

The SKYLINE® Anterior Cervical Plate System consists of an assortment of titanium alloy plates and screws. The plates have two to six screw hole pairs in various lengths. The screws are available in various sizes and screw-tip geometries. Both constrained and variable screws are available to create a constrained, variable or hybrid configuration.

**F. Intended Use**

The SKYLINE Anterior Cervical Plate System is indicated for stabilization of the cervical spine from C2 to C7 employing unicortical screw fixation at the anterior face of the vertebral bodies. Specific clinical indications for anterior plating include: instability caused by trauma; instability associated with correction of cervical lordosis and kyphosis deformity; instability associated with pseudoarthrosis as a result of previously failed cervical spine surgery; instability associated with major reconstructive surgery for primary tumors or metastatic malignant tumors of the cervical spine; instability associated with single or multiple level corpectomy in advanced degenerative disc disease, spinal canal stenosis and cervical myelopathy.

The UNIPLATE Anterior Cervical Plate System and UNIPLATE 2 Anterior Cervical Plate System are intended for anterior cervical intervertebral body fixation. These systems are indicated for patients in which stability is desired following anterior cervical fusion for the indications listed below. The intended levels for treatment range from C2 to T1. Indications include symptomatic cervical spondylosis, trauma, fracture, post-traumatic kyphosis or lordosis, tumor, degenerative disc disease (defined as discogenic pain with degeneration of the disc confirmed by history and radiographic studies), re-operation for failed fusion, or instability following surgery for the above indications.

The UNIPLATE and UNIPLATE2 Anterior Cervical Plate Systems are also indicated for treatment of spinal stenosis.

**F. Summary of Similarities and Differences in Technological Characteristics, Performance and Intended Use**

The proposed modifications to the SKYLINE®, UNIPLATE®, and UNIPLATE® 2 Anterior Cervical Plate Systems are identical to the predicate devices (K042544, K052552, K082273, and K100070) except for the proposed devices will be terminally sterilized via gamma radiation. The design, materials, indications, and technology remain identical to the predicate systems.

**G. Materials**

Manufactured from ASTM F-136 implant grade titanium alloy.

**H. Performance Data**

Performance data is not provided in this submission.

**I. Conclusion**

The Substantial Equivalence Justification demonstrates that the devices are as safe, as effective, and perform as well as the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

DePuy Spine, Incorporated  
% Ms. Kirsten Lehmuller  
Regulatory Affairs Associate  
325 Paramount Drive  
Raynham, Massachusetts 02767

September 6, 2013

Re: K132324

Trade/Device Name: SKYLINE® Anterior Cervical Plate System, UNIPLATE® Anterior Cervical Plate System, and UNIPLATE®2 Anterior Cervical Plate System

Regulation Number: 21 CFR 888.3060

Regulation Name: Spinal intervertebral body fixation orthosis

Regulatory Class: Class II

Product Code: KWQ

Dated: August 6, 2013

Received: August 7, 2013

Dear Ms. Lehmuller:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Mark N. Melkerson -S**

Mark N. Melkerson  
Director  
Division of Orthopedic Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K132324

Device Name: SKYLINE® Anterior Cervical Plate System, UNIPLATE® Anterior Cervical Plate System, and UNIPLATE®2 Anterior Cervical Plate System

### Indications For Use:

The SKYLINE Anterior Cervical Plate System is indicated for stabilization of the cervical spine from C2 to C7 employing unicortical screw fixation at the anterior face of the vertebral bodies. Specific clinical indications for anterior plating include: instability caused by trauma; instability associated with correction of cervical lordosis and kyphosis deformity; instability associated with pseudoarthrosis as a result of previously failed cervical spine surgery; instability associated with major reconstructive surgery for primary tumors or metastatic malignant tumors of the cervical spine; instability associated with single or multiple level corpectomy in advanced degenerative disc disease, spinal canal stenosis and cervical myelopathy.

The UNIPLATE Anterior Cervical Plate System and UNIPLATE 2 Anterior Cervical Plate System are intended for anterior cervical intervertebral body fixation. These systems are indicated for patients in which stability is desired following anterior cervical fusion for the indications listed below. The intended levels for treatment range from C2 to T1. Indications include symptomatic cervical spondylosis, trauma, fracture, post-traumatic kyphosis or lordosis, tumor, degenerative disc disease (defined as discogenic pain with degeneration of the disc confirmed by history and radiographic studies), re-operation for failed fusion, or instability following surgery for the above indications.

The UNIPLATE and UNIPLATE 2 Anterior Cervical Plate Systems are also indicated for treatment of spinal stenosis.

Prescription Use ☒ X ☐

AND/OR

Over-The-Counter Use

☐ ☐

(Part 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Caroline Pham -S**

(Division Sign-Off)

Division of Orthopedic Devices

510(k) Number: K132324