

KINGTECH ENTERPRISES LIMITED

ROOM 2016, 20/F., BLOCK B, REGENT CENTER, 70 TA CHUEN PING ST., KWAI CHUNG, N.T. HONG KONG

TEL : 852-24810188 FAX : 852-24253939

510(K) SUMMARY

for KINGTECH Forehead Thermometer, TF1265, TF1267

This summary of 510(K) safety and effectiveness information is being submitted in accordance with the requirements of 21 CFR 807.92.

The assigned 510(K) number is: K 133113

Submission Date: September 22, 2013

Submitter: Kingtech Enterprises Limited
Room 2016, 20/F., Block B, Regent Center 70 Ta Chunen Ping
Street, Kwai Chung, New Territories, HK
Tel: +852-24810188 Fax: +852-2423939

Manufacturer: Kingtech (Dong Guan) Enterprises Limited
Farm Village, Da Ling Shan Town, DongGuan, China
Tel: +86-76985636260 Fax: +86-76985636350

**Establishment
Registration No.:** 3008808166

Official Contact: Dr. Jen, Ke-Min
Tel: +886-3-5208829 Fax: +886-3-5209783
Email: ceirs.jen@msa.hinet.net

**Common /
Usual Name:** Forehead Thermometer

Trade Name: KINGTECH Forehead Thermometer, Models: TF1265, TF1267

**Classification
Code:** FLL, Class II, 21 CFR 880.2910

Intended Use: TF1265, TF1267 Forehead Thermometer is an electronic thermometer using an infrared sensor to detect human body temperature from the surface of human skin without contact. It is for use by people of all ages in the homecare environment

**Predicated
Devices:** K122221, Nexus IR30 Thermometer, Model: TD-1265
TaiDoc Technology Corporation

Device Description: The Forehead Thermometer TF1265, TF1267 is characterized by measuring human body temperature from the surface of human skin. It utilizes infrared technology to measure infrared energy emitted from the skin surface when making a temperature measurement.

Test Principle: The thermometer Measures temperature by reading infrared radiation emitting from the skin and converts it into a temperature value.

KINGTECH ENTERPRISES LIMITED

ROOM 2016, 20/F., BLOCK B, REGENT CENTER, 70 TA CHUEN PING ST., KWAI CHUNG, N.T. HONG KONG

TEL : 852-24810188 FAX : 852-24253939

Performance Tests:

Safety Test:

- IEC 60601-1 – Medical electrical equipment Part 1. General requirements for safety, 2005.

Electromagnetic Compatibility Test:

- EN/IEC 60601-1-2 – Medical electrical equipment, Part 2. Electromagnetic compatibility – Requirements and tests, 2007.

Clinical Tests:

In accordance with:

EN12470, ASTM E 1965-98, and ASTM E1112-00

Comparison and Conclusion:

This 510k submission only need change the application name of the predicate K122221 from “Nexus IR30 Thermometer, TD-1265” to “KINGTECH Forehead Thermometer, TF1265” there is the entire identical specifications and only need to separate into different 510k. And all of the test reports and documentation for this 510k submission were prepared by the applicant of the predicate device: TaiDoc Technology Corporation. We, Kingtech Enterprises Limited, are authorized by BioCare Co., Ltd. who is a branch office of TaiDoc Technology Corporation. We also present the authorized “Statement” as the following page.

In addition, another model TF1267 of the subject device, it is only minor difference to the exterior dimension for the TF1265.

Thus the new device is substantially equivalent to the predicate devices in this aspect.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

December 26, 2013

Kingtech Enterprises Limited
C/O Dr. Ke-Min Jen
Official Correspondent
Room 2016, 20/F., Block B, Regent Center
70 Ta Chunen Ping Street
Kwai Chung, New Territories
HONG KONG

Re: K133113
Trade/Device Name: KINGTECH Forehead Thermometer, Model: TF1265, TF1267
Regulation Number: 21 CFR 880.2910
Regulation Name: Clinical electronic thermometer
Regulatory Class: II
Product Code: FLL
Dated: September 22, 2013
Received: September 30, 2013

Dear Dr. Jen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Kwame O. for
Ulmer-S

Erin I. Keith, M.S.
Director
Division of Anesthesiology, General Hospital,
Respiratory, Infection Control and
Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K133113

Device Name

KINGTECH Forehead Thermometer, Models: TF1265, TF1267

Indications for Use (Describe)

TF1265 and TF1267 Forehead Thermometer is an electronic thermometer using an infrared sensor to detect human body temperature from the surface of human skin without contact. It is for use by people of all ages in the homecare environment.

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)



Digitally signed by Richard C.
Chapman

Date: 2013.12.26 11:50:34 -05'00'

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."