

APPENDIX A: 510(K) SUMMARY

Sponsor/Submitter: Acclarent, Inc.
1525-B O'Brien Drive
Menlo Park, California 94025

Contact Person: James Patrick Garvey II
Senior Manager, Regulatory Affairs
Phone: 650 687-4807
Fax: (650) 687-4449

Date of Submission: January 17, 2014

Device Trade Name: *Relieva Scout*TM Sinus Dilation System

Common Name: Sinus Dilation System

Device Classification: Class I

Regulation Number: 21 CFR 874.4420

Classification Name: Ear, Nose, and Throat Manual Surgical Instrument

Product Code: LRC

Predicate Device: Relieva Seeker Balloon Sinuplasty System (K120280)

Device Description: The *Relieva Scout* Sinus Dilation System is an integrated device with a low-profile rail, balloon catheter, sinus illumination system with an illuminated ball tip, and an ergonomic handle. The sinus balloon may be inflated to dilate the frontal recess, frontal sinus ostia, and spaces within the frontal sinus cavity.

Indications for Use: For patients aged 18 and older, the *Relieva Scout*TM Sinus Dilation System is intended to provide a means to access the frontal sinus space and to dilate the frontal recess, frontal sinus ostia and spaces within the frontal sinus cavity for diagnostic and therapeutic procedures. In addition, the device is intended to illuminate within and transilluminate across nasal and sinus structures.

**Technological
Characteristics:**

The Scout Balloon Sinuplasty System combines features of a frontal ostium seeker with the tissue expansion effect of balloon dilation. The distal end of the device is permanently curved to optimize frontal ostium access. Light from an extendable integrated illumination system can be seen via transillumination. Additionally, the Scout Balloon Sinuplasty System allows for the manual wire adjustment to provide balloon navigational assistance.

Performance Data:

Bench testing met all acceptance criteria for attributes such as dimensional attributes, cycle fatigue, balloon burst, and bond separation. Testing also showed that the Scout Balloon Sinuplasty System is biocompatible.

The sterilization process was validated per AAMI/ANSI/ISO 11135-1: 2007 and demonstrated a sterility assurance level of 10^{-6} . The method used for sterilization validation will be the overkill (half-cycle approach) in a fixed chamber. Testing of ethylene oxide residuals met ISO 10993-7:2008 requirements. The subject device is not tested nor labeled as “non-pyrogenic”.

Packaging shelf life was established at one year via accelerated aging per ASTM F1980-07.

The performance data demonstrate that the device performs as intended.

**Summary of Substantial
Equivalence:**

The *Relieva Scout* Sinus Dilation System is substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

February 20, 2014

Acclarent, Inc.
% Mr. James Patrick Garvey II
Sr. Manager, Regulatory Affairs, Acclarent Inc.
1525-B O'Brien Drive
Menlo Park, California 94025

Re: K140160

Trade/Device Name: Relieva Scout Sinus Dilation System
Regulation Number: 21 CFR 874.4420
Regulation Name: Ear, nose, and throat manual surgical instrument
Regulatory Class: Class I
Product Code: LRC
Dated: January 22, 2014
Received: January 22, 2014

Dear Mr. Garvey:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Eric A. Mann -S

for Malvina B. Eydelman, M.D.
Director
Division of Ophthalmic and Ear, Nose
and Throat Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

APPENDIX B: INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K140160

Trade Name: *Relieva Scout™* Balloon Dilation System

Common Name: Sinus Dilation System

Indications For Use: For patients aged 18 and older, the Relieva Scout Sinus Dilation System is intended to provide a means to access the frontal sinus space and to dilate the frontal recess, frontal sinus ostia and spaces within the frontal sinus cavity for diagnostic and therapeutic procedures. In addition, the device is intended to illuminate within and transilluminate across nasal and sinus structures.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Sageev George -S

2014.02.20 15:23:02 -05'00'

Page of