



MAY 14 2014 Headquarters
Wright Medical Technology, Inc.

K140408 page 1 of 2 1023 Cherry Road
Memphis, TN 38117

901 867 9971
wmt.com

510(K) SUMMARY

In accordance with the Food and Drug Administration Rule to implement provisions of the Safe Medical Devices Act of 1990 and in conformance with 21 CFR 807.92, this information serves as a Summary of Safety and Effectiveness for the use of the WMT SALVATION™ Osteopenic Screw.

(a)(1). Submitted By:

Wright Medical Technology, Inc.
1023 Cherry Road
Memphis, TN 38117

Date:

May 1, 2014

Contact Person:

Val Myles
Regulatory Affairs Specialist
Office - (901) 290-5162
Fax - (901) 867-4190

(a)(2). Proprietary Name:

SALVATION™ Osteopenic Screw

Common Name:

Bone Screw

Classification Name and Reference:

21 CFR 888.3040 – Class II

Device Product Code, Device Panel:

HWC: Screw, Fixation, Bone

(a)(3). Predicate Devices:

K110670: VLP FOOT Talus and
Percutaneous Calcaneus Bone Plates,
VLP Bone Screws; PERI-LOC Ankle
Fusion Bone Plates and Instruments

K090675: VLP FOOT Plating, Screw
System, and Accessories

(a)(4). Device Description

The SALVATION™ Osteopenic Screws are non-locking, fully threaded screws offered in various diameters and lengths. All described implants are manufactured from titanium alloy and have a solid core. The implants are single use only devices.

(a)(5). Intended Use

The Salvation™ Osteopenic Screw is indicated for the treatment of fracture fixation, osteotomies, and reconstruction/arthrodeses of small bones, as well as patients with osteopenic bone.

(a)(6). Technological Characteristics Comparison

The Salvation™ Osteopenic Screw is technologically substantially equivalent to the predicate device.

(b)(1). Substantial Equivalence – Non-Clinical Evidence

Torsion testing and testing rationales related to bending and pull-out strength were provided to support the substantial equivalence of the subject device and show that no new worst-case devices are introduced in this system.

The safety and effectiveness of the SALVATION™ Osteopenic Screw is adequately supported by the mechanical testing, testing rationales, substantial equivalence information, materials information and comparison of design characteristics provided within this premarket notification. Through the analysis of technical characteristics the new devices are substantially equivalent to the predicate devices.

(b)(2). Substantial Equivalence – Clinical Evidence

N/A

(b)(3). Substantial Equivalence – Conclusions

The design characteristics of the subject system do not raise any new types of questions of safety or effectiveness. From the evidence submitted in this 510(k), the subject devices can be expected to perform at least as well as the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

May 14, 2014

Wright Medical Technology, Incorporated
Ms. Val Myles
Regulatory Affairs Specialist
1023 Cherry Road
Memphis, Tennessee 38117

Re: K140408

Trade/Device Name: SALVATION™ Osteopenic Screw
Regulation Number: 21 CFR 888.3040
Regulation Name: Smooth or threaded metallic bone fixation fastener
Regulatory Class: Class II
Product Code: HWC
Dated: May 2, 2014
Received: May 5, 2014

Dear Ms. Myles:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Ronald P. Jean -S for

Mark N. Melkerson
Director
Division of Orthopedic Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES
Food and Drug Administration

Form Approved: OMB No. 0910-0120

Expiration Date: January 31, 2017

See PRA Statement on last page.

Indications for Use

510(k) Number (if known)

K140408

Device Name

SALVATION™ Osteopenic Screw

Indications for Use (Describe)

The Salvation™ Osteopenic Screw is indicated for the treatment of fracture fixation, osteotomies, and reconstruction/arthrodeses of small bones, as well as patients with osteopenic bone.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

Elizabeth L. Frank -S

Division of Orthopedic Devices