

aap Implantate AG
aap LOQTEQ® Olecranon Plate

Summary of Safety and Effectiveness

MAY 08 2014

Sponsor: aap Implantate AG
Lorenzweg 5
D-12099 Berlin Germany

Company Contact: Dr. Christian Zietsch
Phone: +49-30-750-19 -193
Fax: +49-30-750-19 - 111

Date February 26 / 2013

Trade Name: aap LOQTEQ® Olecranon Plate

Common Name: Olecranon Plate

Classification Name and Reference: 21 CFR 888.3030 Single/multiple component metallic bone fixation appliances and accessories – Class II

Product Code HRS

Predicate device: Synthes (USA) 3.5mm LCP Olecranon Plates, as part of the 3.5 / 4.5 mm LCP System under the premarket notification K082807 (Jan 13, 2009) and Synthes (USA) Small Fragment Dynamic Compression Locking (DCL) System under the premarket notification K000684 (Apr 28, 2000)

Device Description: The aap LOQTEQ® Olecranon Plate consists of bone plates and bone screws, to be implanted by a surgeon in order to achieve an internal fixation of bone fragments typically after fractures. If the plates are used in conjunction with locking screws, a so called internal fixator will be realized (internal fixation).

The aap LOQTEQ® Olecranon Plate System consists of:

- LOQTEQ® Olecranon Plate (left and right)
to be used with
- LOQTEQ® Cortical Screw 2.7, small head T8, self-tapping
- Cortical Screw 2.5, small head T8, self-tapping
- LOQTEQ® Cortical Screw 3.5, T15, self-tapping
- Cortical Screw 3.5, self-tapping
- Set of Instruments aap LOQTEQ® Olecranon Plate

Material: Implants are made of Ti6Al4V (ASTM F136 or ISO 5832-3)

Indications: The aap LOQTEQ® Olecranon Plates are indicated for fixation of fractures, osteotomies and nonunions of the olecranon, particularly in osteopenic bone.

aap Implantate AG
aap LOQTEQ® Olecranon Plate

Substantial Equivalence The Substantial Equivalence of the new device and the predicate device is based on similar intended use, design, functionality, components and materials in use.

Documentation including mechanical testing to show the substantial equivalence has been provided with this submission.

**Performance Data
(Non-Clinical and /
or Clinical):**

Non-clinical tests have been performed and show the substantial equivalence of the device.

Summary of Non-clinical tests:

Type of test:

Customized 4-point bending tests of bone plates, based on ASTM F382-99 (static and dynamic)

Assessment of test results:

Substantial equivalence with respect to the mechanical performance of the aap system could be stated due to the test results gained. The subject device is substantial equivalent, and whose performance meets the requirements of its pre-defined acceptance criteria and intended uses.

Documentation regarding the mechanical testing to show the substantial equivalence has been provided with this submission.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

May 8, 2014

aap Implantate AG
Dr. Christian Zietsch
Manager Regulatory Affairs
Lorenzweg 5
12099 Berlin
Germany

Re: K140607

Trade/Device Name: *aap* LOQTEQ® Olecranon Plate
Regulation Number: 21 CFR 888.3030
Regulation Name: Single/multiple component metallic bone fixation appliances and accessories
Regulatory Class: Class II
Product Code: HRS
Dated: February 28, 2014
Received: March 10, 2014

Dear Dr. Zietsch:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21

CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Lori A. Wiggins

for
Mark N. Melkerson
Director
Division of Orthopedic Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use Statement

510(k) Number (if known): K140607

Device Name: aap LOQTEQ® Olecranon Plate

Indications for Use:

The aap LOQTEQ® Olecranon Plates are indicated for fixation of fractures, osteotomies and nonunions of the olecranon, particularly in osteopenic bone.

Prescription Use ☒ AND/OR Over-The-Counter Use _____
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Elizabeth M. Frank -S

Division of Orthopedic Devices